**Instruction Manual**  
**Manuale di istruzioni**  
**Manuel d'instructions**  
**Manual de instrucciones**  
**Bedienungsanleitung**  
**操作指南**

## **ECODryBlock**

F101A0700

### **General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 基本信息**

Before using the unit, please read the following instruction manual carefully.  
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.  
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.  
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de instrucciones.  
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch.  
使用前请仔细阅读以下使用说明书。

Caution, hot surface!  
Attenzione, superficie calda!  
Attention, surface chaude!  
Prudencia, superficie caliente!  
Vorsicht, heiße Oberfläche!  
注意：加热面高温

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2012/19/UE.  
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2012/19/UE.  
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2012/19/UE.  
No tirar el equipo en los residuos urbanos, como exige la Directiva 2012/19/UE.  
Dieses Gerät unterliegt der Richtlinie 2012/19/UE und darf nicht mit dem normalen Hausmüll entsorgt werden.  
按照 EEC 指令 2002/96/CE，不要将该设备作为城市垃圾处理。

The product can be used with flammable liquids.  
Il prodotto può essere utilizzato con liquidi infiammabili.  
Le produit peut être utilisé avec des liquides inflammables.  
El producto puede utilizarse con líquidos inflamables.  
Das Produkt kann mit brennbaren Flüssigkeiten verwendet werden.  
本产品可与易燃液体混合使用。

**This unit must be used for indoor laboratory applications only.**  
The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a way not specified by the manufacturer or with unspecified accessories, the product's safety may be compromised.

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**IEC/EN 61010-1**  
**IEC/EN 61010-2-051**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

**CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improve their quality.  
Nell'impegno di migliorare costantemente la qualità dei prodotti, VELP si riserva la facoltà di variarne le caratteristiche.  
Dans le but d'améliorer constamment la qualité de ses produits, VELP se réserve le droit d'apporter des modifications aux caractéristiques de ceux-ci.  
VELP se reserva el derecho de modificar las características de productos con el fin de mejorar constantemente su calidad.  
VELP behält sich zum Zwecke der ständigen Verbesserung der Produktqualität das Recht auf Änderung der Geräteeigenschaften vor.  
VELP 保留修改其产品特性的权利，以不断提高其质量。

## Safety Regulations / Norme di Sicurezza / Consignes de Sécurité / Advertencias de Seguridad / Sicherheitshinweise / 安全法规

The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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Use only the power cable provided with the instrument.

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The socket located on the rear panel of the instrument is intended only for powering a single additional unit (code F101A0700). Connecting more than one additional unit of any other type of device may compromise safety.

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The heating blocks can reach temperatures up to 165°C. When removing them, ensure that the displayed temperature is close to room temperature or use the handle for block removal supplied as accessory (A00000351).

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the placement in a hood and personal protective equipment (masks, gloves, goggles, etc.).

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The vessel must be made of a suitable material to withstand the foreseen temperature.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least) and between instruments of 10 cm (at least).

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Do not use with explosive or dangerous materials for which the equipment is not designed. The instrument must not be used in explosive atmospheres or to heat combustible liquids that have a low combustion temperature.

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It is responsibility of the user to appropriately decontaminate the instrument in the case that dangerous substances fall on or in it, according to the safety datasheet of substances used and to the current laboratories' safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts about the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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12. DECLARATION OF CONFORMITY **UK CA** ..... 65

The ECODryBlock is designed to heat samples placed in suitable containers positioned inside precision-machined aluminum blocks.

This instrument is used for various laboratory applications, including sensitive analytical procedure, particularly where precise temperature control of samples is required.

The aluminum blocks are removable and interchangeable, allowing the operator to select the right block depending on the application. This flexibility enables the ECODryBlock to support a wide range of container shapes and sizes, thanks to the wide range of aluminum blocks available.

Each block features holes specifically and precisely shaped to accommodate the container being used.

Thanks to a range of blocks specifically designed for COD applications, the ECODryBlock can be used for Chemical Oxygen Demand (COD) analysis.

The instrument allows simple and easy programming of time and temperature.

Through the digital display, it is possible to monitor the test at any time, with the analysis time constantly shown.

Using default parameters, to ensure maximum safety, the unit automatically stops at the end of the analysis time.

The instrument features simple and easy programming of time and temperature. The digital TFT touch display continuously shows the remaining time or count up, allowing the operator to monitor the process at any time. The ECODryBlock operates safely at temperatures up to 165°C, ensuring consistent and accurate temperature control.

A key feature of the ECODryBlock is its modularity. Using a single ControlPad accessory (code A00000489), it is possible to control up to two Units from a single user interface. Each Unit operates with completely independent settings and parameters. This capability offers great flexibility in laboratory workflows, allowing users to optimize space and reduce the need for multiple different instruments.

**Note:** the Unit (code F101A0700) must be used in conjunction with the ControlPad accessory (code A00000489) to set operating parameters and monitor the instrument's status. Without the ControlPad, the Unit cannot be used. A single ControlPad can control up to two Units (code F101A0700).

Additionally, the ECODryBlock can work in combination with the Pt100 or Pt1000 external probe for precise and accurate control of the block or sample temperature.

- 1 ECODryBlock Unit
- 2 Heating Plate
- 3 Connecting element for ControlPad
- 4 ControlPad

- A Power Input
- B Power socket for additional Unit
- C Main switch
- D Probe socket
- E ControlPad power socket

## 2. Installation

### Unpacking

- Check the integrity of the unit after unpacking.

### The box includes (code F101A0700)

- ECODryBlock Unit (1);
  - Instruction manual;
  - EU (Schuko) power cable;
  - US (NEMA 5-15P) power cable.
- Note:** blocks are not included and must be purchased separately

### The box includes (code A00000489)

- ControlPad (4) for ECODryBlock;
- Connecting element for ControlPad (3) with three screws and Allen key.

### Installation of a single Unit and ControlPad

- Place the Unit on a non-flammable surface;
- Turn the Unit upside down and fasten the Connecting element for ControlPad (3) on the front part of the bottom of the device using the three hex socket head screws and the supplied Allen key;
- Pass the ControlPad cable through the Connecting element for ControlPad, then plug the circular connector into the appropriate socket (E) on the left side of the Unit. Secure the cable to the bottom of the instrument by placing it into the retaining clips;
- Secure the ControlPad to the Connecting element for ControlPad by simply pressing it into place. An audible "click" confirm that the ControlPad is securely attached.

### Installation of an additional unit with accessory "Kit ECODryBlock junction" (A00000531) (optional)

- Turn the Unit upside down and unplug the ControlPad cable from the connector of the unit (E). Remove the cable from the retaining clips;
- Connect the female socket of the power cable provided with the Kit ECODryBlock junction to the upper connector of the Unit (A) to which the ControlPad is attached;
- Position the metal junction element (included in the kit) on the bottom of the Unit to which the ControlPad is attached and secure it with three hex socket head screws;
- Place the second Unit behind the frontal Unit to which the ControlPad is attached and route the power cable provided with the Kit ECODryBlock junction between the metal junction element and the bottom of the second unit;
- Secure the metal junction element to the second Unit with four hex socket head screws;
- Connect the male plug of the power cable provided with the Kit ECODryBlock junction to the lower socket (B) of the second Unit;
- Route the ControlPad cable through the metal junction element, then plug the circular connector into the appropriate socket on the left side of the first Unit (E). Secure the cable again by placing it into the retaining clips;
- It is also possible to install the second Unit without the Kit ECODryBlock junction. In this case, simply connect the Unit to a separate power outlet using the standard power cable provided;
- The Unit to which the ControlPad cable is connected is designated as Unit 1.

### Installation – final steps

- Rotate the instrument back to its normal position;
- Assure that the heating plate is perfectly clean, free from dust and free from scratches. To clean the plate, use a cloth dampened with distilled water;
- Assure that the bottom face of the blocks is perfectly clean, free from dust and free from scratches. To clean the blocks, use a cloth dampened with distilled water;
- Place the blocks inside the Unit;
- **Note:** the Unit must always be used with at least one single block installed. If using only one block, place it on the left side of the heating plate. Operation without blocks is not permitted;
- Make sure that the technical specification of the unit and of the power supply line are the same;  
**Note:** the instrument can operate at both 115V and 230V voltage levels. However, it is possible to change supply voltage (e.g., from 230V to 115V or from 115V to 230V) only when the heating plate has cooled down to ambient temperature.
- Make sure that the electrical network is grounded;
- Make sure that the main switch (C) is OFF;
- Connect the instrument to an easily accessible socket (compliant with the current safety norms), using only the provided power cable.  
**Note:** Keep the power cable far away from the hot plate and from the blocks.
- Place the "Transparent lid" (A00000493) (optional) over the Unit to have a better thermoregulation of samples.

### 3. Operation

#### Power-up

- Turn the unit on using the main switch (C);
- The display shows the software version;
- If no Unit has been paired yet, the ControlPad displays the **Setup** page to allow wireless pairing of one or two Units;  
**Note:** *configure only one ControlPad at a time. Do not pair Units with more than one ControlPad simultaneously.*
- If the Unit has already been paired to the ControlPad, the **Home** page will be displayed, showing the last temperature and timer settings.

**Note:** *if the set starting mode is Stop ([see chap. 4](#)), the heating and timer functions are not active. Otherwise, the instrument starts working at the last set temperature setpoint with the remaining timer duration.*

#### Unit Pairing and Unpairing with ControlPad

- F1 – Available Units Windows: displays the serial number and status of all units within the range of the ControlPad. The possible status are:
  - UNPAIRED: Unit is available and not yet paired with the ControlPad  
→ The Unit can be paired;
  - OFFLINE: Unit is paired but not currently connected  
→ The Unit is paired, but not available for operation (e.g., powered off);
  - ONLINE: Unit is paired and connected to the ControlPad.  
→ The Unit is ready for use;
- To pair and connect a Unit to the ControlPad (status of the Unit: UNPAIRED):
  - Touch the window of the unit to be paired. The window becomes highlighted;
  - Press the “CONNECT” button (P2);
  - The pairing process will start. Once the Unit is connected, its status changes to “ONLINE”.
- To disconnect and unpair a Unit from the ControlPad (status of the Unit: ONLINE or OFFLINE):
  - Touch the window of the Unit to be deleted. The window becomes highlighted;
  - Press the “DELETE” button (P3);
  - The deletion process will start. If the Unit is powered on (ONLINE), once deleted, its status will change to “UNPAIRED”. If the unit is not powered (OFFLINE), once deleted, it will disappear from the list.
- A maximum of two units can be paired to one ControlPad;
- Touch the “X” button (P1) to exit the **Setup page** and access the **Home page**;
- The instrument is now ready for use.

#### Operation

- Refer to [Chapter 4 - User Menu](#) for a detailed description of the main functions of the instrument;
- The ControlPad automatically reduces the brightness if there is no interaction with the display for a few minutes in order to reduce energy consumption.

#### Operation with external probe

- Turn off the instrument using the main switch (C);
- Screw the threaded support rod into its seat on the back of the instrument (optional);


- Fasten the clamp onto the support rod (optional);
- Place the external temperature probe into the clamp (optional); then place it into the sample or the dedicated hole of the aluminum blocks
- Plug the probe into the dedicated socket on the back of the Unit (D). The instrument automatically recognizes the type of temperature probe (Pt100 or Pt1000);
- Turn the Unit on using the main switch (C). When the external probe is connected correctly, the Probe icon (I1) is on, showing the thermoregulation mode selected ([see chap. 4](#));

**Note:** with external probe, the instrument always has primary control of the heating plate temperature. The temperature control feature of the heating plate can also be used as a safety thermostat. In this case, the maximum temperature of the heating plate will not exceed the temperature setting on the instrument ([see chap. 4](#) – Temp. Limit parameter), meaning that a longer heating time is required in order to reach the temperature setting, thus reducing temperature oscillation at a setpoint value.


## 4. User menu

### Home page

One Unit connected to the ControlPad:

Two Units connected to the ControlPad:

- F2 – Temperature Set Window: displays the current temperature setpoint.  
The temperature setpoint can be adjusted using the “+” (P1) and “–” (P2) buttons. Keeping pressed either button will increase the adjustment speed;
- F1 – Real Temperature Window: displays the current temperature of the heating plate or, if connected, of the external probe;
- I1 – Probe Icon: lights up when the external probe is connected. When the instrument is turned on and the external probe is connected, the icon shows “100” or “1000” depending on the type of probe (Pt100 or Pt1000), then shows “Fine” or “Fast” depending on the selected mode ([see chap. 4](#) – User Menu Page – Control Type);
- I2 – Heating Icon: lights up, when START was pressed, to shows that the analysis is running;
- I3 – Hot plate Icon: lights up when the heating plate temperature exceeds 50°C;
- F4 – Timer Set Window: displays the current set timer value.  
The timer hours and minutes can be adjusted using the “+” (P3) and “–” (P6) buttons. Keeping pressed either button will increase the adjustment speed.
- F3 – Timer Window: displays the elapsed time (count-up) or remaining time (count-down) of the timer, depending if the count-down timer is enabled.
- P4 – Count-down Enable Button: activates the count-down timer (button turns blue), which starts from the set value (F4) when the analysis begins or when the temperature setpoint is reached ([see chap. 4](#) – User Menu Page – Time

*Start*). If the count-down timer is not enabled, the “Timer Window” (F3) displays a count-up, regardless of the set count-down value.

The count-down timer can be enabled or disabled at any time during the analysis by simply pressing the “Count-down Enable Button” (P4) (the Timer Window (F3) switches between displaying the count-down and the count-up).

At the end of the count-down, the “Timer Window” (F3) displays “End”, the buzzer of the ControlPad emits an acoustic signal and the timer set value is retained for the next analysis. Heating is stopped only if the “Time End” parameter in the user menu is set to “OFF” ([see chap. 4 – User Menu Page – Time End](#))

- P8 – Lock button: to disable the operation of all buttons. Press and hold the button for 3 seconds. The button turns blue, indicating that the lock is active. To unlock, press and hold the button again for 3 seconds.
- I4 – Start Mode Icon: displays “S” (Stop) or “R” (Run) depending on the selected start mode ([see chap. 4 – User Menu Page – Start Mode](#))
- P7 – Start/Stop Button: to start and stop the analysis. Press the button to start the analysis and activate heating. Press the button while the analysis is in progress to stop it and deactivate heating.
- P5 – Methods Button: to open the **Methods page**. When the button is blue, it indicates that the unit is operating according to a method. To deselect the method, press and hold the button for three seconds.
- P9 – User Menu Button: to open the **User Menu page**.

### User Menu Page

The following parameters are available:

| Parameter name   | Default value | Range                                                      | Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|------------------|---------------|------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Temp. Limit      | ON            | OFF / 30 °C ÷ 165 °C / ON<br>(5 °C steps)                  | <u>Heating plate temperature limit</u> <ul style="list-style-type: none"> <li>➤ OFF: Disables the heating function.</li> <li>➤ Value set between 30 °C and 165 °C: The maximum temperature setpoint is limited to the set value. The heating plate will not exceed this temperature, even when an external probe is connected. When using an external probe, setting the <i>Temp. Limit</i> parameter slightly above the desired setpoint can improve temperature control.</li> <li>➤ ON: When an external probe is connected, the heating plate temperature is limited to the maximum safety value (approximately 185 °C) to ensure that setpoints close to 165 °C can be reached. Suggested for C.O.D. analyses.</li> </ul> |
| Start Mode       | Stop          | Stop – Run                                                 | <u>Start mode</u> <ul style="list-style-type: none"> <li>➤ Stop: when the instrument is powered, the heating function is not active;</li> <li>➤ Run: when the instrument is powered, it starts working at the last working conditions.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Control Type     | Fine          | Fine – Fast                                                | <u>Thermoregulation mode with external probe</u> <ul style="list-style-type: none"> <li>➤ Fine: optimized temperature control for minimized overshoot and oscillations, slow rise in temperature;</li> <li>➤ Fast: fast rise in temperature. Suggested for C.O.D. analysis.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Time Start       | Setpoint      | Snap – Setpoint                                            | <u>Set countdown timer start</u> <ul style="list-style-type: none"> <li>➤ Snap: timer countdown starts as soon as the analysis is started;</li> <li>➤ Setpoint: timer countdown starts when temperature set point is reached.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Time End         | OFF           | ON – OFF                                                   | <u>Continue heating after time end</u> <ul style="list-style-type: none"> <li>➤ If ON, heating continues after timer ends.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Probe AL6        | 07:00         | Enabled / Disabled<br>("Clock" Icon)<br>-<br>00:30 ÷ 30:00 | <u>External probe safety alarm</u><br>It allows to modify the time to show the AL 6 error message ( <a href="#">see chap. 6</a> ) or to disable it.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Probe AL7        | ON            | ON – OFF                                                   | <u>External probe safety alarm</u><br>It allows enabling/disabling AL 7 error message ( <a href="#">see chap. 6</a> ).                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Heat. plate Cal. | 0.0 °C        | -10 ÷ 10 °C<br>0,1 °C Step                                 | <u>Heating plate calibration</u><br>It allows to align the internal sensor reading to a reference thermometer.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Probe Cal.       | 0.0 °C        | -10 ÷ 10 °C<br>0,1 °C Step                                 | <u>Probe calibration</u><br>It allows to align the external temperature probe reading to a reference thermometer.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Key Press sound  | ON            | ON - OFF                                                   | <u>Key Press sound Enable/Disable</u><br>It allows enabling/disabling the buzzer sound when keys are pressed.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |

|             |    |         |                                                                                                                 |
|-------------|----|---------|-----------------------------------------------------------------------------------------------------------------|
| Temp. Scale | °C | °C – °F | <u>Temperature scale</u><br>It allows to switch the temperature scale between Celsius (°C) and Fahrenheit (°F). |
|-------------|----|---------|-----------------------------------------------------------------------------------------------------------------|

Select the parameter to be modified by touching it, then press the “+” and “–” buttons to adjust the value.

The following buttons are available:

- “RESET” Button: to restore all system settings (listed in the table above) to their default values. To perform a reset, press and hold the button for three seconds, until it turns blue;
- “SETUP” Button: to open the **Setup page** (see [Chapter 3. Operation](#)).

To exit the “User Menu Page”, press the “X” button.

### Methods Page

The instrument features a list of standard and custom methods. Each method consists of one to ten programmable ramps, each with individually adjustable temperature and time values.

From the Method Page, it is possible to access the **Standard Methods Page** or the **Custom Methods Page**.

### Custom/Standard Methods Page


This page allows the user to select a method or access the **Method settings Page**. On the **Method settings Page**, users can view or modify the method name, the number of ramps and temperature and time of each ramp.

Method details can be modified only for custom methods. Standard methods are view-only.

To select and import a method, touch the desired method and press the “CONFIRM” button. The display then returns to the **Home Page**, showing the parameters of the selected method. The “Methods” button (P5) on the **Home Page** turns blue indicating that the instrument is operating according to a method.

To access the **Method settings Page**, touch the method and then press the “DETAILS” button.

### Method settings Page

- F2 – Method Name Window: displays the method name.  
Only for custom methods, it is possible to modify the method name:
  - Select the window by touching it. The first letter is highlighted;
  - Use the “+” (P4) and “–” (P5) buttons to select the letter;
  - Touch again the window to switch to the following letter;
  - Repeat the previous two point for all the letters of the name.
- F3 – Number of Ramps Window: displays the number of ramps of the selected method.  
Only for custom methods, it is possible to modify the number of ramps:
  - Select the window by touching it;
  - Use the “+” (P4) and “–” (P5) buttons to select the number of ramps.
- F1 – Ramp Details Window: displays the temperature and timer values associated to each ramp.  
Only for custom methods, it is possible to modify the temperature and timer of the ramps:
  - Select the temperature / hours (hh) or minutes (mm) window by touching it;
  - Use the “+” (P4) and “–” (P5) buttons to select the temperature/timer value.
- P1 – Exit Button: to exit the page and go back to the **Custom/Standard Methods Page**, without saving the changes applied to the method

- P2 – Save Button: to exit the page and go back to the **Custom/Standard Methods Page**, saving the changes applied to the method
- P3 – Confirm Button: to exit the page and go back to the **Home Page**, saving the changes applied to the method and importing the method for an analysis. The “Methods” button (P5) on the **Home Page** turns blue indicating that the instrument is operating according to a method.

## 5. Guidelines for optimal temperature control

### Sample and blocks preparation

- To prevent sample loss due to evaporation and ensure precise thermoregulation of the samples, use sealed containers during the process;
- Ensure the liquid volume inside tubes or vials is below the top edge of the block;
- Select the correct block for the container being used. A block that closely matches the dimensions and shape of the container maximizes heat transfer efficiency and ensures temperature stability within the samples;
- Avoid big air gaps between the vessels and block wells.

### Accurate temperature reading

- When using an external temperature probe, ensure it is fully immersed to the bottom of the vessel for accurate temperature measurement;
- Ensure proper contact between the heating plate and blocks to achieve optimal thermal conductivity. Both surfaces should be smooth, flat (not deformed), and free of scratches. Improper contact can lead to unstable temperature control and slower heating;
- Use two single blocks or one double block to ensure optimal temperature control. If only one single block is used, position it on the left side of the heating plate;
- If the block is supplied with screws (only for “COD” double blocks), it can be secured to the heating plate to prevent movement, minimize the risk of damage, and enhance thermal performance. Before fastening the block to the plate, ensure both contact surfaces are perfectly clean and not damaged. Tighten the screws with low torque to not damage the plate’s threads. In these blocks, the probe should be inserted into the central hole only;
- During operation, consider the physical and chemical properties of the sample, as they can influence temperature stability and measurement accuracy.

**Note:** *only place compatible blocks on the heating plate*

### Temperature control with internal sensor

- When the external probe is not plugged in the unit, the temperature control relies on the internal temperature sensor, placed on the heating plate;
- When using the Internal sensor for temperature control, the set temperature is the one of the heating plate, regardless of the actual temperature of the sample or the blocks;
- The real temperature of the block, and therefore of the samples, may differ from the set temperature. To compensate for this, the device allows the user to apply an adjustable temperature offset to the internal probe reading. To set the offset value, [see chap. 4 – User Menu Page – Heat. plate Cal.](#)

### Temperature control with external sensor

- For more precise temperature regulation, use the external probe (accessories: Pt1000 code A00000458; Pt100 code A00000268) connected to the dedicated socket (D). In this configuration, temperature control is based on the external probe, and the set temperature corresponds to the temperature measured by the probe;
- The use of the external probe allows to thermoregulate directly the sample or the block, ensuring optimal accuracy;
- The external probe can be positioned in the designated hole of the blocks or directly in the sample. Proper probe placement is essential to prevent overheating or inaccurate readings;
- To further improve temperature control, it is possible to adjust the *Temp. Limit* parameter slightly above the setpoint in order to increase stability ([see chap. 4 – User Menu Page – Temp. Limit](#));
- Use sealed vessels to maintain temperature stability during measurement with the external probe.

### External probe adjustment

- Certain sample types or vessel characteristics may cause heat dissipation, causing the sample temperature to differ from the setpoint. To compensate for this, the device allows the user to apply an adjustable temperature offset to the external probe reading;
- To set the offset value, [see chap. 4 – User Menu Page – Probe Cal.](#)

## 6. Error and warning messages

When the display shows an error message, the instrument functions automatically stop.

|            |                                                         |                                                                 |
|------------|---------------------------------------------------------|-----------------------------------------------------------------|
| <b>AL1</b> | Heating plate overtemperature (T > 205 °C)              |                                                                 |
| <b>AL2</b> | Excessive heating time                                  |                                                                 |
| <b>AL4</b> | External probe overtemperature                          | Only if the external probe is connected                         |
| <b>AL5</b> | Hardware malfunction                                    |                                                                 |
| <b>AL6</b> | Slow increase in temperature read by the external probe | Only if the external probe is connected and Probe AL6 is active |
| <b>AL7</b> | Fast decrease in temperature read by the external probe | Only if the external probe is connected and Probe AL7 is ON     |

To remove the error message, restart the instrument from the main switch.

If the alarm persists on the display, please contact VELP Scientifica's technical assistance service.

## 7. Maintenance

No routine or extraordinary maintenance is required except periodic cleaning.

### Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

### Repair

Repairs must be carried out by authorized Velp personnel only.

The transport of the instrument by freight forwarders, couriers or others must be carried out using the original shockproof packaging. Follow any instructions on the original packaging (e.g., palletizing).

It is the responsibility of the user to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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|                           |                                                                                                                                                        |
|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| EN 61010-1                | Safety requirements for electrical equipment for measurement, control, and laboratory use<br>Part 1: General requirements                              |
| EN 61010-2-010            | Particular requirements for laboratory equipment for the heating of material                                                                           |
| ETSI EN 301 489-1 V2.2.3  | Electromagnetic compatibility (EMC) standard for radio equipment and services;<br>Part 1: Common technical requirements                                |
| ETSI EN 301 489-17 V3.3.1 | Electromagnetic Compatibility (EMC) standard for radio equipment and services;<br>Part 17: Specific conditions for Broadband Data Transmission Systems |
| EN 61326-1                | Electrical equipment for measurement, control and laboratory use - EMC requirements -<br>Part 1: General requirements                                  |

and satisfies the essential requirements of the following directives:

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|                        |                                                        |
|------------------------|--------------------------------------------------------|
| 2006/42/EC             | Machinery directive                                    |
| 2014/30/EU             | Electromagnetic compatibility directive                |
| 2015/863/EU (RoHS III) | Restriction of the use of certain hazardous substances |
| 2012/19/EU (WEEE)      | Waste of electric and electronic equipment             |

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung /  
加修改.

## 5. Declaration of conformity <sup>UK</sup> <sub>CA</sub>

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2017/1206                      Radio Equipment Regulations 2017

according to the relevant designated standards:

EN 61010-1                      Safety requirements for electrical equipment for measurement, control, and laboratory use  
Part 1: General requirements

EN 61010-2-010                      Particular requirements for laboratory equipment for the heating of material

ETSI EN 301 489-1 V2.2.3                      Electromagnetic compatibility (EMC) standard for radio equipment and services;  
Part 1: Common technical requirements

ETSI EN 301 489-17 V3.3.1                      ElectroMagnetic Compatibility (EMC) standard for radio equipment and services;  
Part 17: Specific conditions for Broadband Data Transmission Systems

EN 61326-1                      Electrical equipment for measurement, control and laboratory use - EMC requirements -  
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597                      Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032                      Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic  
Equipment Regulations 2012

S.I. 2013/3113                      Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

## Thank you for having chosen Velp!

Established in 1983, Velp is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. Velp works according to **ISO 9001**, **ISO 14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

### Analytical instruments

Elemental Analyzers  
Digestion Units  
Distillation Units  
Solvent Extractors  
Fiber Analyzers  
Dietary Fiber Analyzers  
Oxidation Stability Reactor  
Consumables

### Laboratory Equipment

Magnetic Stirrers  
Heating Magnetic Stirrers  
Heating Plates  
Overhead stirrers  
Vortexers and Shakers  
Dispersers  
Dry Block Heaters and COD  
BOD and Respirometers  
Cooled Incubators  
Flocculators  
Turbidimeter  
Open Circulating Baths  
Pumps

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**Velp Scientifica Srl**