

FILTER APPARATUS

POROSITIES

The porosity for the pore diameter of the filter determined by the method specified in ASTM E-128 "MAXIMUM PORE DIAMETER AND PERMEABILITY OF RIGID POROUS FILTERS FOR LABORATORY USE."

<u>Porosity</u>	<u>Designation</u>	<u>Pore Size Microns</u>	<u>Principal Uses</u>
Fine	F	4-5.5	Filtration of fine precipitates. Extraction
Medium	M	10-15	Extraction. Filtration of crystalline precipitates.
Coarse	C	40-60	Gas Dispersion. Gas Washing.
Extra Coarse	XC	170-220	Extraction. Mercury Filtration. Used as a support for other filter material. Gas dispersion. Gas Washing.

PRESSURE

Fritted glassware is designed for low pressure vacuum filtration or gas flow. The **MAXIMUM** differential pressure should not exceed 15lbs. per square inch.

THERMAL LIMITS

Fritted ware should not be exposed to excessive temperature change or direct exposure to flame due to the fact that it has less resistance to thermal shock than standard borosilicate glassware.

Dry ware may be put in a furnace at the setting of 150°C. However, most drying to constant weight is done at 110°C. Fritted ware can be heated safely to 500°C provided that the cycle of heating and cooling is gradual.

CLEANING

New filters should be cleaned by suction with hot hydrochloric acid followed by a water rinse.

Used filters can usually be cleaned by rinsing with water, passed thru the underside of the filter disc at pressure not exceeding 15lbs. per square inch. There are some precipitates that tend to clog the pores of a fritted filter. Following are some suggestions for dealing with these precipitates that have to be removed by chemical means.

MATERIAL

Albumen

Aluminous and siliceous residues

Copper or Iron Oxides

Glucose

Fatty Materials

Mercuric Sulfide

Mercury

Organic Matter

Silver Chloride

CLEANING AGENT

Hot ammonia or hydrochloric acid

2% hydrofluoric acid followed by concentrated sulfuric acid. Rinse immediately with water until no trace of acid can be detected.

Hot hydrochloric acid plus potassium chlorate

Hot mixed acid ($H_2SO_4 + HNO_3$)

Carbon Tetrachloride

Hot aqua regia

Hot Nitric Acid

Hot concentrated cleaning solution, or hot concentrated sulfuric acid with a few drops of sodium nitrate.

Ammonium or sodium hyposulfide