by VELP Scientifica

Instruction Manual

CG-1999-V-10

CG-1999-V-10E

Heating Magnetic Stirrer

General Information

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 Caution, hot surface!
Attention, surface chaude!

 Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2002/96/CE.

This unit must be used for laboratory applications only. The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a not specified way by the manufacturer or with non specified accessories, product's safety may be compromised.

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This unit has been designed and manufactured in compliance with the following standards:
L'instrument a été conçu et fabriqué conformément aux normes suivantes:

Safety requirements for electrical equipment for measurement, control and for laboratory use	IEC/EN 61010-1
Règles de sécurité pour appareils électriques de mesurage, de régulation et de laboratoire	IEC/EN 61010-2-051
Electrical equipment for laboratory use	UL 61010-1
General requirement - Canadian electrical code	CAN/CSA-C22.2 No.61010-1

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it. It is also responsibility of the user to use safety substances for cleaning or decontaminating, which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor. Est responsabilité de l'utilisateur la décontamination en cas de déversement de matières dangereuses sur ou à l'intérieur de l'équipement. Est responsabilité de l'utilisateur à utiliser des substances qui ne produisent pas de danger pour le nettoyage ou de décontamination, qui ne réagissent pas avec les parties internes de l'appareil ou avec la matière qu'il contient. En cas de doute sur la compatibilité d'une solution de nettoyage, contactez le fabricant ou le distributeur local.

Contents

1.	INTRODUCTION.....	3
2.	ASSEMBLY AND INSTALLATION.....	3
	2.1 ELECTRICAL CONNECTIONS.....	3
	2.2 START-UP	3
3.	OPERATING CONTROLS	3
	3.1 TURNING ON AND REGULATING	3
	3.2 USING THE THERMOREGULATOR VTF.....	3
	3.3 USING THE EXTERNAL PROBE.....	4
	3.4 SETTING MODE.....	4
	3.5 ERROR MESSAGES	4
4.	MAINTENANCE	4
	4.1 CLEANING	4
5.	TECHNICAL DATA	4
6.	WIRING DIAGRAM	5
7.	ACCESSORIES	5
8.	WARRANTY	5
9.	DECLARATION OF CONFORMITY 	6

The AREC.X heating magnetic stirrer with ceramic plate, is used in the laboratory for general use and for all those applications that require a precise regulation of the stirring speed and of the heating plate temperature.

The ceramic heating plate (1) means that higher temperatures can be obtained compared to the traditional aluminium heating plate and it takes less time to heat the sample. Ceramic is an inert and very hard material and is resistant to almost any type of chemical or mechanical aggression; its surface characteristics remain unaltered over time.

The AREC.X is used to mix liquids using a magnetically driven stir bar placed inside the vessel. The size and shape of the magnetic stir bar determines the stirring efficiency at any given speed. It is therefore advisable to choose the most suitable magnetic stir bar in relation to the quantity and characteristics of the liquid and the type of vessel used.

The stirring bar which satisfies most applications is code A00001060 (10x40 mm).

The instrument is turned on using the ON/OFF button. The stirring speed and the temperature (4) of the heating plate are adjusted using the relative knobs on the front panel, the values are indicated on the digital display (5).

NOTE: If the heating function is used the vessel must be made of a suitable material to withstand the foreseen temperature.

2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- Heating Magnetic Stirrer
- Power supply cord
- Instruction manual

2.1 Electrical connections

After having unpacked the instrument, place the unit on the laboratory bench. Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. The equipment shall be connected to a power socket with protective earth connection, using only the power supply cord provided with the instrument. Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

NOTE: the mains cable must remain far away from the hot plate.

2.2 Start-up

Rotate the speed (right) and temperature (left) knobs completely to the left. Place the flask containing the sample and a suitable magnetic stirring bar on the stirring plate. Then, set the speed and temperature by turning the dedicated knobs.

3. Operating controls

3.1 Turning on and regulating

Connect the unit to mains. In the next 5 seconds, displays shows the software version and the restarting mode (A or b, chapter "Setting mode"). Turn the unit on using the On-Off button. The green led indicates that the instrument is On.

Start the stirring function by turning the stirrer knob "Stirrer rpm" on the front panel. Speeds of from 50 to 1500 rpm can be selected using the analogical scale around the knob. Bright bars rotating on the display indicate that the instrument is running. The speed of rotation of these bars is proportional to the stirring speed.

NOTE: in case of black-out, once the power is back the device will restart in the restarting mode set (OPt tyP), selected in chapter "Setting mode".

To start the heating function turn the heating knob "Heating °C" on the front panel. Temperatures of from 5 to 550 °C can be selected at 5 °C intervals. The temperature setting is shown on the display. When heating is turned Off (Heating knob OFF) and the heating plate temperature exceeds 50 °C, the message "Hot" flashes on the display until the temperature of the heating plate falls to below 50 °C (or until the instrument is turned off using the ON/OFF).

NOTE: this warning is not active if the instrument is not connected to the power supply.

The knobs on the front panel are easily accessible and are well away from the heating plate in order to ensure maximum operator safety as well as safeguarding the electronic components inside the unit.

3.2 Using the thermoregulator VTF

Screw the threaded support rod into its seat on the back of the instrument and fasten the VTF thermoregulator onto the support rod. Place the temperature probe in the flask making sure that it is completely immersed in the sample.

Connect the two instruments (AREC.X and VTF) by plugging the VTF into the dedicated socket on the back of the AREC.X. Select the operating temperature required on the VTF thermoregulator.

Turn the temperature control knob on the front panel of the AREC.X to the desired temperature.

NOTE: when the two instruments are connected, the display of the AREC.X does not show the temperature but a dashed line (---). The operating temperature is shown on the display of the VTF thermoregulator.

When using the VTF thermoregulator always select the max temperature on the AREC.X.

The temperature control function of the heating plate can also be used as a safety thermostat.

In this case the maximum temperature of the heating plate will not exceed the temperature setting on the AREC.X meaning that a longer heating time is required in order to reach the VTF thermoregulator temperature setting.

3.3 Using the external probe

Screw the threaded support rod (optional accessory) into its seat on the back of the instrument and fasten the clamp (optional accessory) onto the support rod.

Fit the temperature probe onto the clamp making sure that it is completely immersed in the sample.

Plug the temperature probe into the dedicated socket on the back of the AREC.X with the instrument off. A led on the left-hand side of the display indicates that the probe is connected, the display indicates the probe temperature reading. Select the operating temperature using the heating knob on the AREC.X. A led at the center of the display indicates that the set-point temperature is shown; after 5 seconds the led goes out and the display indicates the probe temperature.

3.4 Setting mode

Turn the left knob twice from the minimum to the maximum during the visualization of the software version in order to enter into the setting mode.

Once there, by turning the left knob, it is possible to pass from one of the following parameter to the next ones:

Parameter shown		Default value	Range	Description
CAL	P10	0,0	-10,0 ÷ +10,0°C	External probe calibration
FS	tCJ	550	0 ÷ 550°C	Maximum temperature settable
OPt	tyP	b	A – b	(*) Restarting mode
OPt	P10	FInE	FInE - FASt	(**) Type of control with external probe

Turn the right knob from the minimum to the maximum and then back to the minimum in order to pass to the parameter you want to modify. The display will show the value saved for 3 seconds, and then the new value according to the position of the right knob.

The new value is automatically saved after 4 seconds, if during this time the right knob is not moved (no change of value). "End" is shown and then the system goes back to the main menu.

Once the new values are set, switch off and restart the instrument.

(*) When **A** is selected: in case of black-out, once the power is back the device will remain off.

When **b** is selected: in case of black-out, once the power is back the device will automatically restart with the latest settings.

(**) Select **FInE** to minimize overshoot and fluctuations of temperature.

Select **FASt** to save time.

3.5 Error messages

If an error message appears on the display, please contact VELP Scientifica's technical service department.

When the display shows an error message, the stirring and heating functions stop automatically.

NOTE: To remove the error message disconnect the instrument from the power supply.

AL1 Thermocouple is out of range

AL2 Excessive heating time

AL3 The stirring system doesn't run correctly

AL4 Plate temperature higher than 252 °C or the probe is disconnected with instrument ON

Should any of the above occur, please contact your nearest VELP Scientifica service centre.

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual.

In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors.

The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.


5. Technical data

Power supply	230V / 50-60Hz or 115V / 60Hz
Dimensions WxHxD	203x93.5x344 mm (8x3.68x13.5 in)
Weight	3.5 Kg (7.7 lbs)
Overall power	800 W
Dimension of the heating plate	180 x 180 mm
Programmable temperature range	5 – 550 °C
Type of temperature control	Digital
Overtemperature protection	Yes

Stirring capacity	15 liters of H ₂ O
Programmable speed range	50 – 1500 rpm
Type of motor control	Electronic
Connection	VTF - External probe Pt100
Electric rating for external probe	5 Vdc – 1W (max)
Accuracy of temperature with VTF	+/-0,5°C*
Accuracy of temperature with external probe	+/-1°C*
Noise level	<< 80 dBa
Temperature range	+5...+40 °C
Storage temperature range	-10...+60 °C
Max humidity	80%
Level of electrical protection CEI EN60529	IP 42
Overvoltage category	II
Pollution degree CEI EN61010-1	2
Max altitude	4000 m

* in the following conditions: 800ml water in 1 liter glass beaker (diameter 105mm), stirring bar 10x40mm, 600rpm, 50°C.

6. Wiring diagram

7. Accessories

Please get in contact with Chemglass Life Sciences for more details about accessories.
www.cglifesciences.com Phone: 1-800-843-1794

8. Warranty

The unit is guaranteed against production defects for 25 months from our invoice date.

In accordance with this guarantee Chemglass Life Sciences undertakes to repair any units resulting as faulty due to the quality of the materials used or poor workmanship.

Units rendered faulty due to inexpert handling/use or carelessness will not be replaced or repaired under warranty.

Exclusions:

The guarantee will be considered null and void for faults resulting from:

- inexperience and carelessness of the operator
- repairs, maintenance or replacement of parts carried out by personnel or Companies not authorized by the manufacturer
- use of the instrument that does not comply to the instructions/recommendations given in the present operating manual
- use of non-original spare parts.

9. Declaration of conformity

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

EN 61010-1 (2001)
EN 61326-1 (2006)
2011/65/EU (RoHS)
2002/96/CE (RAEE)

and satisfies the essential requirements of the following directives:

- Machinery directive 2006/42/EC
- Low voltage directive 2006/95/EC
- Electromagnetic compatibility directive 2004/108/EC
- plus modifications