

• BelloCell-500 Technical Report VI



Cultivation of suspend rCHO cells for protein secretion

Description BelloCell-500 provides a powerful tool to achieve high cell density and high productivity of target bioproducts in a cell culture because it has a unique feature of offering high oxygen transfer and low shear stress culture environment. Users can easily collect highly concentrated cells, virus or secreted products from one 500 ml BelloCell-500 bottle. In this study, the applications of BelloCell-500 for growth of suspend CHO cells and production of human IgG is illustrated. 7.9×10^7 CHO cells were seeded in one BelloCell-500 unit. A final 4.5×10^9 cell population in one BelloCell -500 unit was obtained. For the protein production, a total of 157.5 mg IgG protein was harvested in 5 L conditioned culture media. In comparison with 1.1 mg IgG protein in a 200 ml spinner flask. It is 5.7 folds productivity increase on an equal volume basis. This technical sheet provides a general protocol for users to start up their culture. However, the optimum condition of each cell culture for each case may require the users to determine.

Material

Device	Cell Line/Product	Medium	Seed
BelloCell-500	CHO/IgG	Ex-CELL 301 (JRH)	7.9×10^7 cells/bottle

Protocol *Detail protocol is in General Instruction Manual

Inoculum preparation Prepare one 250 ml spinner flasks and inoculate 1.8×10^5 suspend cells/ml in 100 ml Ex-CELL 301 culture medium. Culture at 50 rpm, 37°C for 3 days. After cell density reaches above 8×10^5 cells/ml and viability remain above 95%, it is ready for the preparation of inoculation. Collect 5×10^7 to 1.0×10^8 suspend cells from the spinner flask by centrifugation and collect in one 50 ml centrifuge tube with 50 ml fresh media.

Preparation before cell seeding Place BelloStage-3000 controller in a 37°C incubator. Set up the inoculation parameters (See below). Warm up Ex-CELL 301 medium in 37°C water bath. Take out one BelloCell-500 bottle aseptically and place it in a biosafety cabinet. Open the cap and add 450 ml fresh culture medium in the bottle.

Inoculation Open the cap and dispense 50 ml media containing 5×10^7 to 1.0×10^8 suspend cells on top of the matrix of BelloCell-500. Bring the bottle and lock up on the BelloStage-3000 controller in incubator at 37°C. Press “START” button to start the controller. Avoid swirling or shaking the bottle before start compression.



Culture The inoculation parameters were set as below:

Rising rate	Top holding time	Down rate	Bottom holding time
2.0 mm/s	20 sec	2.0 mm/s	0 sec

After 2 to 3 hours, reset the parameters for culture condition. Usually, above 90% cells will be immobilized in the matrices within 2 hours. The culture parameters were set as below:


Rising rate	Top holding time	Down rate	Bottom holding time
1.5 mm/s	0 sec	1.5 mm/s	1 min

For optimal culture, monitor the pH, residual glucose concentration and other metabolic in order to predict the time for medium replenishment.

When pH drop below 7.0 during culture, add 5 ml 1M HEPES and/or 10 ml 7.5% NaHCO₃ to re-adjust the pH to above 7.5. When glucose concentration below 1.5 g/L, add 5 ml 18 w/w% glucose and 125 mM glutamine concentrate. Replenish culture medium once a day from the day 3rd of culture. Add glucose, glutamine, HEPES and NaHCO₃ supplements once a day from the day 5th of culture. We suggest that users can replenish culture medium in the afternoon, and add supplement the next day in the morning. For simpler operation, user can replenish medium twice a day in the morning and in the afternoon.


Results

1. Immobilization Efficiency






2. Glucose uptake rate (mg/hr) and glutamine uptake rate (umol/hr)




3. Glucose and Lactate concentration






4. IgG Production




5. IgG and GUR



GUR shows a good correlation with IgG production.




6. IgG concentration between BelloCell-500 and Spinner flask



Max. IgG concentration in BelloCell-500 is about 6 folds compared with spinner flask


7. IgG production comparison with spinner flask



Total harvest IgG in BelloCell-500 is 71 folds compared with a 500 ml spinner flask
 Average IgG concentration in BelloCell-500 is around 7 folds compared with a 500 ml spinner flask




8. Total cell population and cell density based on same volume



Total cell population in BelloCell-500 is 7.5 folds compared with a 500 ml spinner flask; cell density in BelloCell-500 is 8.2 folds compared with a 500 ml spinner flask.

9. Medium consumption based on same productivity



For the same production, it will consume 28 L culture medium in a 500 ml spinner to get the same productivity as one BelloCell-500 bottle.

- BelloCell-500 Technical Report VI •
- Cultivation of CHO cells for protein secretion



Summary

Seed	Inoculum volume	Medium volume	Medium
7.9×10 ⁷ cells/total	50ml	500 ml	Ex-CELL 301
Total culture age	Total medium consumed	Total medium replenish frequency	Total cell counted
15 days	5000 ml	9 times	4.5×10 ⁹ cells/total
Max. GUR	Total Protein produced	Max. Protein concentration	Multiplication of cells
121.71 mg/hr	157.5 mg	61.2 ug/ml	57 folds

Please contact Cescobio Bioengineering Technical Support for any questions or comments

[http:// www.cescobio.com.tw](http://www.cescobio.com.tw)

e-mail: info@cescobio.com.tw