

# FINDENSER

conserves water and minimizes risk of lab flooding

- No risk of flooding from running water
- Eliminate water purchase and disposal costs
- For solvent volumes from 5mL up to 1L
- Helps meet sustainable water reduction targets

Water is a precious resource. **Findenser™ requires no running water to operate.** It makes little economic or environmental sense to waste thousands of liters just to cool a single condenser.

**Patented Design**  
sealed-in thermofluid and air-cooling fins combine for excellent thermal conductivity


**Ground Glass Connections**  
fit standard tapered glass joints

**Anti-Roll Design**  
contoured edges prevent rolling and fin damage


**Easy Clamping**  
with standard laboratory clamps

**Reaction Blocks and Inserts**  
See Part No. CG-1992-series

**Choice of Flasks**  
10mL to 2L flasks  
1L maximum solvent volume



**Findenser consists of an internal borosilicate glass condenser wrapped in a finned, aluminum, heat sink jacket.**



The area between the glass condenser and the anodized aluminum jacket is filled with a small amount of water, permanently sealed in. The design of the glass condenser maximizes the internal surface area compared to conventional glass condensers; improving the overall heat transfer capacity. Findenser replaces the need for water cooled condensers in more than 95% of common chemistry applications.

## Findenser Compared with an Air Condenser

For synthesis with 'low' boiling point solvents, Findenser showed a significant improvement in solvent retention. With acetone or DCM the reaction boiled dry when using an air condenser, yet Findenser retained 95% of the solvent under the same conditions.

For synthesis with 'medium' boiling point solvents, Findenser delivered improved solvent retention particularly with larger volumes and high temperatures.

**Note:** Always ensure vessel and Findenser are clamped. Do not expose to temperatures below 0°C. Do not autoclave. Do not exceed 60°C when oven drying. Do not exceed 50°C when using in a dishwasher.


Part No.	Description
CG-1218-F114	Findenser Mini, 275mm, B14 Outer and Inner Joints
CG-1218-F124	Findenser Mini, 275mm, B24 Outer and Inner Joints
CG-1218-F224	Findenser, 400mm, B24 Outer and Inner Joints

\*Not Subject To Discount.



1-800-843-1794  
www.cglifesciences.com  
tech@cglifesciences.com

Average Water Consumption Compared



Findenser TM Radleys



See the other side for our economic data logger