Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung

OHS 200 ADVANCE Overhead stirrers

F20100480, F20110480

OHS 100 ADVANCE Overhead stirrers

F20100481, F20110481

OHS 60 ADVANCE Overhead stirrers

F20100482, F20110482

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de funcionamiento.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2002/96/CE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2002/96/CE.
No tirar el aparato en los desechos urbanos, como exige la Directiva 2002/96/CE.
Dieses Gerät unterliegt der Richtlinie 2002/96/EG und darf nicht mit dem normalen Hausmüll entsorgt werden.

This unit must be used for laboratory applications indoor only. The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a not specified way by the manufacturer or with not specified accessories, product's safety may be compromised.

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**IEC/EN 61010-1 IEC/E
61010-2-051**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The values indicated on the rating plate of the instrument must correspond to those of the power supply.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Safe working conditions are ensured only when the accessories described in the dedicated chapter are used.

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The working speed set on the instrument must be such as to avoid wobbling and/or splashes.

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Do not use with explosive or dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie or to stir harmful liquids prior using protective measures according to the safety standards of the processed products and/or in force in the laboratories including personal protective equipment and the presence of an extraction hood which ensures at least 10-fold air change in accordance with the standards EN 14175 and DIN 12924.

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It is dangerous to run the unit with the stirring blade turning in free air. Always place the stirring shaft in the receptacle before turning the unit on.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it.

It is also responsibility of the user to use safety substances for cleaning or decontaminating, which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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The solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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Switch off the the stirring before opening chuck locking ring. Switch off the instrument before removing chuck.

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The instrument contains a battery. Its substitution must be carried out by authorized Velp personnel only.

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OHS overhead stirrers with electronic speed control, brushless motor, and advanced safety features are able to satisfy the most difficult laboratory applications in terms of viscosity and volume. The new chuck ensures higher safety for the operator and allows to use a passing rod of up to 8,5mm in diameter.

Figure 1. OHS 200 Advance

Figura 1. OHS 200 Advance – Frontal view

- A Display LCD
- B Display keys
- C Speed control knob
- D Lock key
- E Main switch
- F Smart Chuck


2. Assembly and installation

- Unpacking
 - Check the integrity of the unit after unpacking.

- The box includes
 - OHS advance overhead stirrer
 - Power supply cable
 - Instruction manual
 - Allen key
 - Allen head screw
 - Shaft support overhead stirrer
 - Pt100 probe Ø 3mm with cable

- First installation
 - Place the unit on non-flammable surface
 - Fix the shaft support to the instrument with the allen head screw
 - Fasten the unit to the support rod (A00000366 or A00000369) using the double clamp (A00001301)
 - Secure the receptacle using the ribbon clamp (A00001302)
 - Slide the stirring rod into the chuck and tight it
 - Make sure that the rating value of the instrument corresponds to the one of the power supply
 - Ensure that the socket provided, with grounding, is compliant with the current safety norms and easy to reach. Use only the cable provided with the instrument.
 - Insert the power cable into the socket.

3. Display symbols

4. Working

Commissioning	<ul style="list-style-type: none"> ➤ Switch on the instrument using the main switch (E) ➤ Display (A) shows Welcome page and the main screen <div style="display: flex; justify-content: space-around;"> <div data-bbox="344 232 660 651"> </div> <div data-bbox="794 232 1110 651"> </div> </div> <div style="display: flex; justify-content: space-around; margin-top: 10px;"> <div data-bbox="320 667 683 696"> <p>Figure 4. Main OHS 60 Advance</p> </div> <div data-bbox="783 667 1161 696"> <p>Figure 5. Main OHS 200 Advance</p> </div> </div>
Stirring	<ul style="list-style-type: none"> ➤ Adjust speed set point by turning the speed control knob (C). As soon as the knob is moved, set rpm (6) becomes blue. ➤ Click the speed control knob (C) to start stirring. ➤ Speed increases until set point is achieved. ➤ A microprocessor ensures constant speed even when the viscosity changes (counter-reaction). ➤ Switch off the stirring by clicking the knob (C).
Timer	<ul style="list-style-type: none"> ➤ Click Timer button (12) to select the timer. ➤ Set timer (8) time becomes blue. Adjust timer by turning the speed control knob (C). ➤ Click the knob (C) to confirm. ➤ If the instrument is already working, timer countdown starts immediately, otherwise timer (7) is fixed as set timer (8) until stirring begins. ➤ If timer is not set, set timer (8) shows hh:mm:ss and timer (7) works as a counter.
Gear (only for OHS 200)	<div style="display: flex;"> <div data-bbox="320 1182 687 1664" style="flex: 1;"> </div> <div data-bbox="730 1182 1541 1473" style="flex: 2;"> <ul style="list-style-type: none"> ➤ Set the operating speed range by clicking Gear button (11) and rotating the speed control knob (C). ➤ Click the knob (C) to confirm the speed range. Once confirmed, Current Gear symbol becomes black. ➤ It's possible to select between two gears: <ul style="list-style-type: none"> ➤ Gear I: Low speed (6 – 400 rpm) High torque (till 200 Ncm). ➤ Gear II: High speed (30 – 2000 rpm) Low torque (till 40 Ncm). </div> </div> <div style="margin-top: 10px;"> <p>Figure 6. Set Gear</p> </div>
Torque	<ul style="list-style-type: none"> ➤ The intensity of the torque is indicated on the main screen (9).
Lock	<ul style="list-style-type: none"> ➤ Holding the Lock key (D) for 3 seconds, the instrument will lock its settings during operations. ➤ Unlock the control panel by holding the Lock key (D) for 3 seconds. ➤ If other buttons are clicked while the instrument is locked, the two LEDs aside lock button blink for many seconds.

5. External Connections

USB	<ul style="list-style-type: none"> ➤ Each model has USB connection in the back for PC controlling, data logging and software upgrading by PC. ➤ Software version, dedicated program to be installed in the PC, and software installation guide have to be requested by e-mail to service@velp.it
Pt100	<ul style="list-style-type: none"> ➤ All models have Pt100 connection in the back for fluid temperature measurement (measuring range from -200°C to +550°C)

Controller register assignment

Register address	Register name	Write	Byte Number	Function/explanation
1	Instrument	NO	2	Instrument model
2	Serial number	NO	6	Instrument serial number
5	Product code	NO	12	Instrument product code
11	Main board sw	NO	8	Main board sw version
15	Display board sw	NO	8	Display board sw version
24	State	NO	2	Instrument state
25	Alarm	NO	2	The instrument notifies an alarm
26	Gear	NO	2	Gear set (only for OHS200)
27	Speed	NO	2	rpm measured
28	Torque	NO	2	Ncm measured
29	Timer	NO	4	Residual timer or counter
31	Timer set point	NO	4	Timer set value
33	Speed set point	NO	2	Speed set value
48	Pt100 connection	NO	2	Pt100 connected to the instrument
49	Temperature	NO	2	Temperature indication
50	Speed limit	NO	2	Speed maximum value
51	Torque limit	NO	2	Torque maximum value
52	Ramp	NO	2	Setting of acceleration
53	Method	NO	2	Reading of method number
54	Method n steps	NO	2	Total number of steps
55	Method step	NO	2	Current step visualized
56	Method n loops	NO	2	Total number of loops
57	Method loop	NO	2	Current loop visualized
102	Gear	SI	2	Setting of gear (only for OHS200)
103	Speed	SI	2	Setting of the speed
104	Timer	SI	4	Setting of the timer
106	Motor stop	SI	2	Motor stops stirring
107	Motor start	SI	2	Motor starts stirring

Sampling time: 1s or more

Examples (CRC 16 bit)

Reading single register: Speed (rpm)

Request		Reply		
Field	(Hex)	Field	(Hex)	Description
Address	0x64	Address	0x64	
Control command	0x03	Control command	0x03	
High start address	0x00	Number of bytes	0x02	
Low start address	0x1B	High register value	0x01	340 (RPM)
Number of High registers	0x00	Low register value	0x54	
Number of Low registers	0x01	High CRC	0xF4	
High CRC	0xFD	Low CRC	0x23	
Low CRC	0xF8			

Reading multiple register: Serial number

Request		Reply		
Field	(Hex)	Field	(Hex)	Description
Address	0x64	Address	0x64	
Control command	0x03	Control command	0x03	
High start address	0x00	Number of bytes	0x06	
Low start address	0x02	High register value	0x31	12345 0x00 0x31 0x32 0x33 0x34 0x35
Number of High registers	0x00	Low register value	0x00	
Number of Low registers	0x03	High register value	0x33	
High CRC	0xAD	Low register value	0x32	
Low CRC	0xFE	High register value	0x35	
		Low register value	0x34	
		High CRC	0x0A	
		Low CRC	0x0A	

Writing single register: STOP rotation

Request				
Field	(Hex)			Description
Address	0x64			
Control command	0x06			
High start address	0x00			
Low start address	0x6A			
High register value	0x00			Stop rotation
Low register value	0x01			
High CRC	0x61			
Low CRC	0xE3			

Virtual serial port

Baudrate	9600
Bits	8
Stop Bit	1
Parity	None
Maximum number of registers for single request	24

6. VELP Ermes Configuration

VELP Ermes is a revolutionary cloud platform that transforms and improves your laboratory experience by creating an ecosystem of instruments, people and data. The VELP Ermes platform is able to reduce distances and accelerate scientific processes in total safety. In order to access on ERMES, you need to enable your VELP account by selecting "Configure your VELP ERMES account" at <http://www.velp.com/en/login>.

To be able to communicate, the instrument needs to be in the operating range of laboratory Wi-Fi (2.4 GHz) and be configured as follows:

- Switch on the OHS Advance and select AP in the menu Set-up Wi-Fi (see chapter 7.5).
- Using the PC/Tablet/Mobile phone, select the OHS_SERIAL NUMBER available on the Wi-Fi list, in order to connect directly to the instrument.
- Open a browser on the PC/Tablet/Mobile phone and insert the address 192.168.240.1 to reach the configuration page. Insert "admin" "admin" when requested as user name and password.
- Set the parameters required to connect to your Wi-Fi (network name, password, security, mac address, etc.) according to your internal procedure and save. If necessary, contact your IT administrator.
- Select menu Ermes from the Service menu (see chapter 7.6) and proceed with the product registration from the VELP Ermes platform. For more information see FAQ on VELP website.

NOTE: To access to VELP Ermes is necessary to have a VELP account.

7. Menu

Clicking Menu button (13) the following figure appears


Figure 7. Menu OHS 200 Advance



Figure 8. Menu OHS 100 - 60 Advance

It's possible to move among submenus rotating the speed control knob (C).

Enter in a submenu with a click of the speed control knob (C) once it's highlighted in blue. Press Home to go back to the main screen.

7.1 Menu Structure

Service →

- Reset Torque
- Real Torque
- Pt100 Alignment
- Check Locking Ring
- Reset Parameters
- Update Software
- Counter
- Ermes
- Motor Calibration

Timer

Gear *

* Graph for models OHS 100 and OHS 60

** Not visualized for models OHS 100 and OHS 60

7.2 Methods

In this submenu is possible to set 10 different methods. For OHS 200, 10 different methods for each speed range. Once set, a method becomes white in the method list and a blue bar appears on the left. Press START to begin the method.

Figure 9. Methods


Figure 2. Methods set


Figure 11. Stop Method

When a method is working, on the main screen the current method (15) is shown.

It's not possible to modify speed, gear, or the timer.

It's possible to navigate in menu.

To stop a method before its end, enter in the method list and click STOP.

7.2.1 Method set-up

In each method it's possible to set parameters in Figure 12. Once method parameters are set, they are highlighted in white with a blue bar on the left as shown in Figure 13.

Figure 3. Method parameters

12 Figure 4. Method parameters set

- 5 steps can be set for each method.
- All steps are programmable with speed, timer, ramp, and intermittent mode.
- A method is considered set when at least speed and timer are set for one step.


Figure 5

Steps

- Rpm: set point speed
- Time: countdown performed for each step visualized also on the main screen
- Ramp: if is selected, OHS stirs for all the time at the speed set
If is selected, OHS reaches the speed set in time selected
- Interm. Mode: if YES is selected, but no intermittent mode is set for the method, OHS works in continuous mode.

Start Time

- It allows to set an hour at which the method starts.
- Enter in the menu, rotate speed control knob (C) till the desired hour. Click the knob to confirm.
- Even if a Start Time is set, the method begins if START button in method list is clicked.

Intermittent Mode

- It allows to set stirring period and stop period alternatively.
- Enter in the menu, rotate speed control knob (C) till the desired working time. Click to confirm.
- Rotate speed control knob (C) till the desired pause time. Click to confirm.
- Click again the knob to modify set values.

Loop

- It allows to repeat the whole method N times.
- Enter in the menu, rotate speed control knob (C) till the desired value. Click to confirm.

Click RESET to erase all method values.

7.3 Graph

In this submenu is possible to visualize the graph of the current analysis. For models OHS 100 and OHS 60 is possible to access to the graph directly through the main screen.

Click RESET to erase the graph. Only the last 60 minutes are shown. Click ZOOM to see the last 3 minutes.

Figure 6. Graph with method

Figure 7. Graph without method


When a method is set graph appears as in Figure 15. When a method is not set, graph appears as in Figure 16.



7.4 Safety

In this submenu it's possible to set all limits linked to the safety depending on the working conditions

Speed Limit	<ul style="list-style-type: none"> ➤ It allows to set the speed full scale. ➤ The maximum speed value is set by default (for OHS 200 it depends on the gear selected at the moment). ➤ Enter in the menu, rotate speed control knob (C) by 100 rpm steps till the desired value. Click to confirm.
Torque Limit	<ul style="list-style-type: none"> ➤ It allows to set the torque limit. ➤ The maximum value is set by default (for OHS 200 it depends on the gear selected at the moment). ➤ Enter in the menu, rotate speed control knob (C) by 10 Ncm steps till the desired value. Click to confirm.
Temperature Limit	<ul style="list-style-type: none"> ➤ It allows to set a speed reduction if Pt100 detects a certain temperature increase (or decrease) in a defined time.
	<p>Figure 8. Temperature Limit</p> <ul style="list-style-type: none"> ➤ Enter in the menu, rotate speed control knob (C) till the delta temperature desired value. Click to confirm. ➤ Rotate the knob (C) to set the interval time in which the delta temperature has to be considered. Click to confirm. ➤ Rotate the knob (C) to set the percentage of speed reduction desired if the selected delta temperature is detected in the set interval time. Click to confirm. ➤ Speed reduction can be set by 10% steps.
Acceleration	<ul style="list-style-type: none"> ➤ It allows to choose among 3 different acceleration types when the instrument starts to stir or when a higher set point value is set. ➤ Enter in the menu, select Slow, Medium, or Fast depending on the customer application. ➤ Click speed control knob (C) to confirm.
Vibration sensor	<ul style="list-style-type: none"> ➤ It allows to set a vibration sensitivity level. ➤ Enter in the menu, rotate speed control knob (C) to select OFF, Low, Medium or High. Click to confirm. ➤ When vibration sensor level is set different than OFF,  appears on the blue upper bar (3).

7.5 Set-up

Start Mode	<ul style="list-style-type: none"> ➤ The decision of the instrument's restart mode in case of blackout or power loss. ➤ Enter in the menu, rotate speed control knob (C) to select Stop or Run. Click to confirm. <ul style="list-style-type: none"> ➤ Stop: when the instrument is switched on, it's requested a click of the control knob to start stirring. ➤ Run: when the instrument is switched on, it restarts to work with the last set point set. ➤ If Stop is selected, A appears on the blue upper bar (2). ➤ If Run is selected, B appears on the blue upper bar (2).
Wi-Fi Set-up	<ul style="list-style-type: none"> ➤ It allows to switch on wi-fi module for IoT transmission. ➤ Enter in the menu, rotate speed control knob (C) to select ON or OFF. Click to confirm. ➤ When wi-fi is ON,  appears on the blue upper bar (1).
Wi-Fi Information	<ul style="list-style-type: none"> ➤ It describes all wi-fi information (MAC address, Power of signal, wi-fi name).
Display	<ul style="list-style-type: none"> ➤ It allows for brightness to be set or if torque and temperature should be displayed on the main screen. ➤ Enter in the menu and scroll with speed control knob (C). Click to enter in submenus. ➤ Brightness: enter in this submenu, rotate speed control knob (C) to select the desired brightness value. Click to confirm. ➤ Brightness can be set by 10% steps. ➤ Torque: enter in this submenu, rotate speed control knob (C) to choose between ON or OFF if the display of torque on the main screen is desired or not. Click to confirm. ➤ Temperature: enter in this submenu, rotate speed control knob (C) to choose between ON or OFF if the display of temperature on the main screen is desired or not. Click to confirm ➤ If torque and temperature are displayed on the main screen, the display menu appears as in Figure 18. ➤ When temperature is displayed but Pt100 probe is not inserted in the instrument, the main screen appears as in Figure 19. <div style="display: flex; justify-content: space-around; align-items: flex-start; margin-top: 20px;"> <div style="text-align: center;">  <p>Figure 9. Torque and temperature shown</p> </div> <div style="text-align: center;">  <p>Figure 10. Main without temperature probe</p> </div> </div>
Time & Date	<ul style="list-style-type: none"> ➤ It allows to set hour and date. ➤ Enter in the menu and select Set Time. Rotate speed control knob (C) till the right time. Click to confirm. ➤ Return to the previous page clicking BACK button. ➤ Select Set Date. Rotate speed control knob (C) till the right day. Click to confirm ➤ Repeat the operation for month and year.
Temperature Unit	<ul style="list-style-type: none"> ➤ It allows to choose temperature unit that has to be visualized between °C and °F.
Language	<ul style="list-style-type: none"> ➤ It allows to select the interface language. ➤ Enter in the menu, rotate speed control knob (C) to select language. Click to confirm.
7.6 Service	
Reset Torque	<ul style="list-style-type: none"> ➤ It allows to reset the current torque. ➤ Click with the speed control knob (C) to reset torque.
Real Torque	<ul style="list-style-type: none"> ➤ It allows to return to real torque value. ➤ Click with the speed control knob (C) to return to real torque value. ➤ Negative torques can't be visualized anymore.
Pt100 Alignment	<ul style="list-style-type: none"> ➤ It allows for the alignment of the Pt100 probe to a reference thermometer. ➤ Enter in the menu, rotate speed control knob (C) to select the desired alignment value (from -10.0°C to 10.0°C). Click to confirm.

<p>Check Locking Ring</p>	<ul style="list-style-type: none"> ➤ Alignment sensitivity 0.1°C. ➤ It allows to control the right functioning of Smart Chuck system. ➤ Enter in the menu. Figure 20 is visualized. <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  <p>Figure 12. Check Locking Ring</p> </div> <div style="text-align: center;">  <p>Figure 11. Check Locking Ring performed</p> </div> </div> <ul style="list-style-type: none"> ➤ Follow instructions described in the screen. ➤ Once the operation is performed, Check becomes white and a blue bar appears on the left.
<p>Reset Parameters</p>	<ul style="list-style-type: none"> ➤ It allows to return to default value for all functions. ➤ Click with the speed control knob (C) to reset parameters.
<p>Update Software</p>	<ul style="list-style-type: none"> ➤ It allows to update the device with a new software version. ➤ Software version, dedicated program to be installed in the PC, and software installation guide have to be requested by e-mail to service@velp.it.
<p>Counter</p>	<ul style="list-style-type: none"> ➤ It allows to see the number of instrument working hours.
<p>Ermes</p>	<ul style="list-style-type: none"> ➤ It allows to connect the instrument to Ermes Cloud.
<p>Motor Calibration</p>	<ul style="list-style-type: none"> ➤ It allows to reset the zero torque value.

8. SmartChuck™

This new mechanism (F) allows to change the stirring shaft with one hand.

- When instrument is off, hold the stirring shaft with one hand.
- From work position (Figure), turn Locking Ring 90° to the left (Figure).
- The main screen displays Figure 22.
- Pull Locking Ring down (Figure 25).
- Open the chuck.
- Remove the shaft.
- Introduce a new stirring shaft and place it in the operating position.
- Ensure the stirring shaft is centered as best as possible. Then close the chuck.
- Push Locking Ring up (Figure 26) and turn it 90° to the right.
- The yellow upper bar on the display becomes blue and the OHS is ready to be used.


Figure 13. Locking Ring open

Figure 23. Working position

Figure 24. Intermediate position

Figure 15. Open position

Figure 14. Intermediate position

9. Error messages

When the display shows an error message, the stirring function stops automatically.

Error code	Cause
AL1	Motor doesn't start stirring
AL2	High internal motor temperature
AL3	Motor overload
AL4	High driver temperature
AL5	Safety relay fault
AL10	Vibrations too high
AL11	Temperature too high (only with Pt100 inserted)
AL12	Temperature too low (only with Pt100 inserted)

To remove the error message, disconnect the instrument from the power supply.

If alarm persists on the display, please contact VELP Scientifica's technical service department. service@velp.it

10. Maintenance

Maintenance	➤ No routine or extraordinary maintenance is necessary;
	➤ Repairs must be carried out by authorized Velp personnel only;
	➤ Instrument must be transported in its original packaging any indications present on the original packaging must be followed (e.g. palletized);
	➤ It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.
Cleaning	➤ Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

11. Technical data

Models		F20100480, F20100481, F20100482	F20110480, F20110481, F20110482	
General features	Power supply	230 V / 50-60 Hz (+/-10%)	115V / 60 Hz (+/-10%)	
	Dimensions (WxHxD)	OHS 200	90x315x235 mm (3.54x12.40x9.25 in)	90x315x235 mm (3.54x12.40x9.25 in)
		OHS 100 – OHS 60	90x285x235 mm (3.54x11.22x9.25 in)	90x285x235 mm (3.54x11.22x9.25 in)
	Weight	OHS 200	4,6 kg (10,14 lb)	4,6 kg (10,14 lb)
		OHS 100	4,1 kg (9,04 lb)	4,1 kg (9,04 lb)
		OHS 60	4,1 kg (9,04 lb)	4,1 kg (9,04 lb)
	Power input	190 W	190 W	
	Construction material (structure)	Aluminum	Aluminum	
	Working in continuous	Admitted	Admitted	
	Settable restart modality	Stop or work	Stop or work	
	Noisiness	<< 60 dBa	<< 60 dBa	
	Environmental temperature admitted	+5...+40 °C	+5...+40 °C	
	Storage temperature admitted	-10...+60 °C	-10...+60 °C	
	Max humidity	80%	80%	
	Level of electrical protection CEI EN60529	IP 54	IP 54	
	Overvoltage category	II	II	
Pollution degree CEI EN61010-1	2	2		
Max altitude	2000 m	2000 m		
Stir	Stirring capacity	OHS 200	100 l H ₂ O	100 l H ₂ O
		OHS 100	100 l H ₂ O	100 l H ₂ O
		OHS 60	40 l H ₂ O	40 l H ₂ O
	Programmable speed range	OHS 200	6-400rpm (I) – 30-2000rpm (II)	6-400rpm (I) – 30-2000rpm (II)
		OHS 100	30-1300rpm	30-1300rpm
		OHS 60	30-2000rpm	30-2000rpm
	Motor type	BLDC	BLDC	
	Speed selection	1 rpm step	1 rpm step	
Stirring alarm	Motor fault	Motor fault		
Motor rating output	150 W	150 W		
Torque	Max torque admitted	OHS200	200 Ncm (I) – 40 Ncm (II)	200 Ncm (I) – 40 Ncm (II)
		OHS100	100 Ncm	100 Ncm
		OHS60	60 Ncm	60 Ncm
Counters	Motor counter	Working hours	Working hours	

12. Accessories / Spare parts

A00000002	Probe extension cable, length 1 m	A00001306	Stirring shaft with fixed blade
A00000363	Glass Temperature Probe OHS Advance	A00001307	Stirring shaft with propeller
A00000366	Support rod and base OHS	A00001308	Stirring shaft with paddle, 6 holes
A00000369	Telescopic support rod and base OHS	A00001309	Stirring shaft with turbine
A00000391	OHS/ControllerSoft	A00001310	Stirring shaft with turbo propeller
A00000372	Support rod for probe	A00001311	Stirring shaft with anchor
A00001300	Support rod and base	A00001312	Stirring shaft protection OHS
A00001301	Double clamp		
A00001302	Ribbon clamp	10007175	Knob 35D blue
A00001304	Stirring shaft with floating blade	40002632	Pt100 probe Ø 3 OHS Advance with cable
A00001305	Stirring shaft with folding blade		

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-1 (2010)	Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use - Part 1: General Requirements
EN61010-2-051 (2015)	Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-051: Particular requirements for laboratory equipment for mixing and stirring
EN 62479:2010	Assessment of electronic and electrical equipment related to human exposure restrictions for electromagnetic fields (10M Hz - 300 GHz).
EN 301 489-1 V2.2.0	Electromagnetic compatibility and Radio spectrum Matters (ERM); ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Common technical requirements.
EN 301 489-17 V3.2.0	Electromagnetic compatibility and Radio spectrum Matters (ERM); ElectroMagnetic Compatibility (EMC) standard for radio equipment; Specific conditions for Broadband Data Transmission Systems.
EN 300 328 V2.1.1	Electromagnetic compatibility and Radio spectrum Matters (ERM); Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz ISM band and using wide band modulation techniques; Harmonized EN covering the essential requirements of article 3.2 of the R&TTE Directive
2011/65/EU (RoHS)	Restriction of the use of certain hazardous substances
2012/19/EU (RAEE)	On waste electrical and electronic equipment
FCC CFR 47 Part 15 Subpart	Electronic Code of Federal Regulations – Industrial, Scientific and Medical equipment applications and Authorizations.
UL 61010-1	Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use – Part 1: General Requirements.
CAN/CSA-C22.2	Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use – Part 1: General Requirements.

and satisfies the essential requirements of the following directives:

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Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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VELP Scientifica Srl