# Guardian™ 3000 User Guide

English

Español

Français

Deutsch

Italiano

Português

Svenska

Nederlands

Dansk

Polski

Čeština

Magyar

## 1. SAFETY INFORMATION

### Definition of Signal Warnings and Symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

WARNING	For a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.
CAUTION	For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or minor or medium injuries if not avoided.
ATTENTION	For important information about the product. May lead to equipment damage if not avoided.
NOTE	For useful information about the product.

### Warning Symbols

General hazard

Explosion hazard

Caution, hot surface

Protective conductor terminal

Alternating current

Electrical shock hazard

### Safety Precautions

**WARNING!** The protection provided by the unit may be impaired if used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer.

- Always operate unit on a level surface for best performance and maximum safety.
- DO NOT lift unit by the top plate.
- To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplugging from the wall outlet.
- Disconnect unit from the power supply prior to maintenance and servicing.
- Spills should be removed promptly, after the unit has cooled down.
- Alkalis spills, hydrofluoric acid or phosphoric acid spills may damage the unit and lead to thermal failure.
- DO NOT immerse the unit for cleaning.
- DO NOT operate the unit at high temperatures without a vessel/sample on the top plate.
- DO NOT operate the unit if it shows signs of electrical or mechanical damage.
- Protective earthing of the equipment is achieved via connection of the provided power cord to a compatible grounded power outlet.

**WARNING:** unit is not explosion proof. Use caution when unit is on or when heating volatile materials.

**WARNING!** DO NOT use the unit in explosive atmospheres or with materials that could cause a hazardous environment from processing. Keep in mind the material flash point relative to the target temperature that has been set.

**CAUTION!** The top plate can reach 500°C, DO NOT touch the heated surface. Use caution at all times. Keep the unit away from explosive vapors and clear of papers, drapery and other flammable materials. Keep the power cord away from the heater plate.

**CAUTION! Beware of the following risks when heating.**

- Flammable materials
- Low boiling point combustible substances
- Glass breakage as a result of mechanical shaking power
- Incorrect container size
- Too much medium
- Unsafe condition of container

Earth Ground - Protective Conductor Terminal

Alternating Current

### Intended Use


This instrument is intended for use in laboratories, pharmacies, schools, businesses and light industry. It must only be used for processing materials as described in these operating instructions. Any other type of use and operation beyond the limits of technical specifications, without written consent from OHAUS, is considered as not intended. This instrument complies with current industry standards and the recognized safety regulations; however, it can constitute a hazard in use. If the instrument is not used according to these operating instructions, the intended protection provided by the instrument may be impaired.

## 2. BASIC FUNCTIONS


Quick Press

Long Press




Press &amp; Hold

Rotate

Note:  refers to beep in the following table.

 <b>I/O – Standby Switch</b>	<b>Turn unit on</b> (Press switch)	
	<b>Turn unit off</b> (Press switch)	
 <b>Temp Knob Functions</b>	<b>Change temp settings</b> (Rotate temp knob)	
	<b>Turn heater on</b> (Long press temp knob until unit beeps & green light appears)	
	<b>Change temperature settings, while heating</b> (Rotate temp knob, the setting will blink, then quick press knob)	 (If setting is idle for 6 sec, setting will return to prior value)
	<b>Turn heater off</b> (Long press temp knob until unit beeps & green light fades)	
	<b>Enter settings menu</b> <i>[heater &amp; stirrer must be off]</i> (Press & hold temp knob until "MENU" appears)	
 <b>Speed Knob Functions</b>	<b>Change speed settings</b> (Rotate speed knob)	
	<b>Turn stirrer on</b> (Long press speed knob until unit beeps & green light appears)	
	<b>Change speed settings, while stirring</b> (Rotate speed knob, the setting will blink, then quick press knob)	 (If setting is idle for 6 sec, setting will return to prior value)
	<b>Turn stirrer off</b> (Long press speed knob until unit beeps & green light fades)	

### 3. CONTROL PANEL

1. **Display Screen**
2. **Left Knob:** Controls temperature and settings menu.
3. **Standby Switch**
4. **Heater Indicator:** Illuminates when heater is running.
5. **Hot Top Caution Indicator:** Illuminates once the heater temperature reaches 40 °C.
6. **Speed Indicator:** Illuminates when stirrer is running.
7. **Right Knob:** Controls speed.

**Note:** The top plates are of different dimension and material. Please refer to the actual product for detail.

### 4. MENU

To access the settings menu, the heater and stirrer must be off. Then press and hold the Temp/Menu knob until “MENU” appears on the screen. (The heater will not turn on unless the knob is released before “MENU” appears.)

To navigate the settings menu, rotate the knob then quick press the knob to choose the selection. Once in the selection, rotate the knob to change the setting and quick press the knob to confirm.

To return to the top level of the settings menu, rotate the knob to “BACK” then quick press the knob. To return to the main operating screen from the top level of the settings menu, rotate the knob to “ESC” then quick press the knob.

To exit the menu at any time, flip the standby switch off and back on. The unit will return to the main operating screen and will be ready for normal use.

#### Menu Features:


##### 1. “CAL” – Single Point Calibration

Single Point Calibration (SPC) improves the accuracy of the heater at user-selected temperature points. SPC can be used during plate-controlled and probe-controlled heating.

##### 2. “SYS” – System Settings

- i. Enabling / Disabling the Beeper
- ii. Changing the Power Recovery Setting
- iii. Resetting to Factory Default Settings

## 5. DISPLAY

1. **Heater Temperature:** Switches to external probe temperature when the probe is plugged in and external prob icon is illuminated.
2. **Heat Setting Indicator:** Switches heater temperature to heat setting when illuminated.
3. **Single Point Calibration Icon**
4. **External Probe Icon**
5. **Speed Setting Indicator:** Illuminates until stirrer reaches the speed setting.
6. **Stir Speed**

## 6. MAINTENANCE

### 6.1 Cleaning

**WARNING:** Electric Shock Hazard. Disconnect the equipment from the power supply before cleaning.

Make sure that no liquid enters the interior of the instrument.

**Attention:** Do not use solvents, harsh chemicals, ammonia or abrasive cleaning agents.

The housing may be cleaned with a cloth dampened with a mild detergent if necessary.

### 6.2 Replacing Power Fuse

The fuse is located at the back of the instrument. Choose the correct fuse to replace according to the fuse specification listed in the specification table.

**WARNING:** Electric Shock Hazard. Disconnect the equipment from the power supply before replacing the fuse.

The use of a fuse of a different type or with a different value, or bridging or shunting the fuse is not allowed and can possibly cause a hazard to your safety and lead to instrument damage!

### 6.3 Troubleshooting

The following table lists common problems and possible causes and remedies. If the problem persists, contact OHAUS or your authorized dealer.

Table 6-1 Troubleshooting

Error*	Cause of Error	How to Fix
Unit fails to power on	Missing or blown fuse	Add or replace fuse as necessary.
E1	Plate RTD open	Not fixable by user, please contact Ohaus.
E2	Plate RTD short	Not fixable by user, please contact Ohaus.
E3	No stirring motion / cannot reach speed	Not fixable by user, please contact Ohaus.
E4	Probe RTD open (Removing the probe while the unit is heating)	Switch unit to standby, then return to normal operating mode.
E5	Probe RTD short (Malfunctioning probe)	Switch unit to standby, remove the probe from the unit, then return to normal operating mode.
E6	A/D lock error	Not fixable by user, please contact Ohaus.
E7	User Probe Error (Plugging the probe into the unit while it is heating)	Switch unit to standby, then return to normal operating mode.
E8	Plate over temperature	Not fixable by user, please contact Ohaus.
E9	Plate under temperature	Not fixable by user, please contact Ohaus.
E10	Triac fault	Not fixable by user, please contact Ohaus.
E11	Heater damage	Not fixable by user, please contact Ohaus.
E12	Heater overheating	Switch unit to standby, then return to normal operating mode.

**\*Note:** Error code instances will stop equipment operation by default.

### 6.4 Service Information

If the troubleshooting section does not resolve your problem, contact an authorized OHAUS Service Agent. For Service assistance in the United States, call toll-free 1-800-526-0659 between 8:00 AM and 5:00 PM Eastern Standard Time. An OHAUS Product Service Specialist will be available to assist you. Outside the USA, please visit our website [www.ohaus.com](http://www.ohaus.com) to locate the OHAUS office nearest you.

## 7. TECHNICAL DATA

Equipment Ratings:

Indoor use only

Altitude: 0 to 2000m

Operating temperature: 5°C to 40°C

Storage temperature: -20°C to 65°C

Operating humidity: 20 to 80% relative humidity, non-condensing

Storage humidity: 20 to 80% relative humidity, non-condensing

Voltage fluctuations: Mains supply voltage fluctuations up to  $\pm 10\%$  of the nominal voltage.

Overvoltage category II

(Installation category):

Pollution degree: 2

Table 7-1 Specification

Model	e-G31HSRDS	e-G31HS04C	e-G31HS07C
Overall Dimensions (L x W x H)	264 x 163 x 101 mm (10.4 x 6.4 x 4.0")	264 x 163 x 110 mm (10.4 x 6.4 x 4.3")	318 x 223 x 112 mm (12.5 x 8.8 x 4.4")
Top Plate Dimensions	Ø 13.5 cm (5.3")	10.2 x 10.2 cm (4 x 4")	17.8 x 17.8 cm (7 x 7")
Top Plate Material	Ceramic Coated Stainless Steel	Ceramic	
Electrical (50 / 60 Hz)	100-120V~ 5A / 220-240V~ 4A		100-120V~ 10A / 220-240V~ 7A
Fuses	10 A, 5x20 mm, 250 VAC / 6.3 A, 5x20 mm, 250 VAC		15 A, 5x20 mm, 250 VAC / 10 A, 5x20 mm, 250 VAC
Temperature Range	Ambient +5°C to 380°C	Ambient +5°C to 500°C	
Temperature Stability of Top Plate*	$\pm 3\%$ (> 100°C), $\pm 2^\circ\text{C}$ ( $\leq 100^\circ\text{C}$ )		
Temperature Stability with Temperature Probe**	$\pm 2\%$ (> 100°C), $\pm 2^\circ\text{C}$ ( $\leq 100^\circ\text{C}$ )		
Stir Capacity	15 L		
Speed Range	80 to 1600 rpm		
Speed Stability	$\pm 2\%$		
Weight Capacity	Up to 15 kg (33 lbs)		
Shipping Dimensions (W x D x H)	330 x 250 x 140 mm (13.0 x 9.8 x 5.5")		390 x 290 x 160 mm (15.4 x 11.4 x 6.3")
Shipping Weight	3.1 kg	3.1 kg	4.6 kg

Note:

- + The parameter in the specification table are applicable to 2" diameter center of the top plate.
- ++ Conditions permitting. Variations in temperature measurement process, vessel, ambient and sample will impact actual performance. To improve the accuracy of the system, please use the Single Point Calibration feature.

## COMPLIANCE

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
	This product complies with the applicable harmonized standards of EU Directives 2011/65/EU (RoHS), 2014/30/EU (EMC) and 2014/35/EU (LVD). The EU Declaration of Conformity is available online at <a href="http://www.ohaus.com/ce">www.ohaus.com/ce</a> .
	This product complies with the applicable statutory standards of the Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012, UK Electromagnetic Compatibility Regulations 2016 and Electrical Equipment (Safety) Regulations 2016. The UK Declaration of Conformity is available online at <a href="http://www.ohaus.com/uk-declarations">www.ohaus.com/uk-declarations</a> .
	This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to <a href="http://www.ohaus.com/weee">www.ohaus.com/weee</a> .
	EN 61326-1
	CAN/CSA-C22.2 No. 61010-1, CAN/CSA-C22.2 No. 61010-2-010, CAN/CSA-C22.2 No. 61010-2-051
	UL 61010-1, UL 61010-2-010, UL 61010-2-051

### ISED Canada Compliance Statement:

CAN ICES-003(A) / NMB-003(A)

### ISO 9001 Registration

The management system governing the production of this product is ISO 9001 certified.

Scan the QR code to download the full instruction manual or visit <http://www.ohaus.com/support> to view the full user manual.

Guardian™ 3000  
Instruction Manual

**FCC Supplier Declaration of Conformity**

Unintentional Radiator per 47CFR Part B

Trade Name: OHAUS CORPORATION

Model or Family identification: e-G31

**Party issuing Supplier's Declaration of Conformity:**

Ohaus Instruments (Changzhou) Co., Ltd.

2F, 22 Block, 538 West Hehai Road, Xinbei District, Changzhou

Jiangsu 213022

China

Phone: +86 519 85287270

**Responsible Party – U.S. Contact Information:**

Ohaus Corporation

7 Campus Drive, Suite 310

Parsippany, NJ 07054

United States

Phone: +1 973 377 9000

Web: [www.ohaus.com](http://www.ohaus.com)

**FCC Compliance Statement:**

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

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P/N 30642337B © 2022 Ohaus Corporation, all rights reserved / todos los derechos reservados / tous droits réservés / alle Rechte vorbehalten / tutti i diritti riservati / todos os direitos reservados / Alla rättigheter förbehållna / alle rechten voorbehouden / alle rettigheder forbeholdes / wszelkie prawa zastrzeżone / všechna práva vyhrazena / Minden jog fenntartva