**Instruction Manual**  
**Manuale di istruzioni**  
**Manuel d'instructions**  
**Manual de instrucciones**  
**Bedienungsanleitung**  
**操作指南**

## AREX 5 Digital Heating Magnetic Stirrer

F20500580, F20510580

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 基本信息

Before using the unit, please read the following instruction manual carefully.  
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.  
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.  
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de instrucciones.  
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch.  
使用前请仔细阅读以下使用说明书。

Caution, hot surface!  
Attenzione, superficie calda!  
Attention, surface chaude!  
Prudencia, superficie caliente!  
Vorsicht, heiße Oberfläche!  
注意：加热面高温

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2012/19/UE.  
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2012/19/UE.  
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2012/19/UE.  
No tirar el equipo en los residuos urbanos, como exige la Directiva 2012/19/UE.  
Dieses Gerät unterliegt der Richtlinie 2012/19/UE und darf nicht mit dem normalen Hausmüll entsorgt werden.  
按照 EEC 指令 2002/96/CE，不要将该设备作为城市垃圾处理。

The product can be used with flammable liquids.  
Il prodotto può essere utilizzato con liquidi infiammabili.  
Le produit peut être utilisé avec des liquides inflammables.  
El producto puede utilizarse con líquidos inflamables.  
Das Produkt kann mit brennbaren Flüssigkeiten verwendet werden.  
本产品可与易燃液体混合使用。

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**This unit has been designed and manufactured in compliance with the following standards:**

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**IEC/EN 61010-1  
IEC/EN 61010-2-051  
IEC/EN 61010-2-010**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

**CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improve their quality.

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The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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Use only the power cable provided with the instrument.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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The vessel must be made of a suitable material to withstand the foreseen temperature.

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Magnetic field doesn't interfere with the function of cardiac pacemakers or data media if they are farther than 20 cm to the instrument.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Do not use with explosive or dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie or to stir combustible liquids that have a low combustion temperature. The product is intended for use with very small quantities of flammable liquids or flammable liquids that have a fire point higher than 625°C and a flash point higher than 600°C.

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It is responsibility of the user to appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories' safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts about the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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
## 1. Introduction

The AREX 5 Digital heating magnetic stirrer is used to heat and mix liquids inside a suitable container placed on the aluminium plate. The stirring is possible thanks to a magnetic drive stir bar placed into the vessel.

The instrument is used in the laboratory for general use and for all those applications that require precise regulation of the stirring speed and of the heating plate temperature.

- 1 Aluminum alloy plate with ceramic coating
- 2 Main switch
- 3 Speed control knob
- 4 Temperature control knob
- 5 Display

- A Probe led
- B Temperature setpoint led

The Aluminium alloy heating plate with ceramic coating (1) ensures:

- Optimum heat distribution and a high specific power thanks to the circular configuration
- Temperature homogeneity
- High resistance to thermal stress and thermal shock
- High resistance to chemical agents and surface abrasion

**Note:** Using the heating plate at high temperatures may cause discolouring. This does not alter the thermal, mechanical and chemical resistance of the plate in any way.

The powerful brushless motor is able to stir from 50 up to 1500 rpm. SpeedServo technology maintains constant speed as viscosity changes.

**Note:** it's important to choose the most suitable magnetic stirring bar related to the type and the quantity of the liquid to be stirred as well as to the type of the beaker. Size and shape of the magnetic stir bar determines the stirring efficiency at any given speed. The stirring bar which satisfies most applications is code A00000356 (8 x 40 mm).

The AREX 5 Digital can work in combination with the Pt100 or Pt1000 external probe or with thermoregulator VTF for precise and accurate control of the sample temperature.

## 2. Installation

- Unpacking
  - Check the integrity of the unit after unpacking
- The box includes
  - AREX 5 Digital heating magnetic stirrer
  - Instruction manual
  - Power cable
- Installation
  - Place the unit on a non-flammable surface
  - Make sure that the technical specification of the unit and of the power supply line are the same
  - Make sure that the electrical network is grounded
  - Make sure that the main switch is on position "0" (OFF)
  - Connect the instrument to an easily accessible socket (compliant with the current safety norms), using only the provided power cable

**Note:** keep the power cable far away from the hot plate.

- Place the container with the liquid and the stir bar on the stirring plate

## 3. Operation

### Power-up

- Turn the unit on using the main switch
- The display shows the software version and the last setpoint values

If the set starting mode is 5 $\pm$ 0 $\bar{P}$  (see chap. 4), the heating and the stirring functions are not active and both displays show 0 $\bar{F}$ F. Otherwise, the unit starts working at the last set temperature and speed setpoints.

### Stirring

- Rotate the right knob to set the stirring speed
- Push the knob to start stirring
- The real speed increases until set point achievement
- A microprocessor ensures constant speed even if the viscosity changes (SpeedServo)
- Push the knob to stop stirring

### Heating

- Rotate the left knob to set the heating plate temperature
- Push the knob to start heating
- The display shows the current heating plate temperature for 3 seconds and the setpoint for 1 second
- When the display shows the temperature setpoint, the Set led is on
- Push the knob to stop heating

**Note:** if the heating function is not active and the heating plate temperature is higher than 50 °C, the display shows the blinking message "Hot". This warning is active even if the main switch is off but not if the instrument is not powered.

### Operation with thermoregulator VTF

- Turn off the instrument using the main switch
- Screw the threaded support rod into its seat on the back of the instrument and fasten the VTF thermoregulator onto the support rod. Place the temperature probe in the flask making sure that it is completely immersed in the sample
- Plug the VTF into the dedicated socket on the back of the instrument
- Turn the unit on using the main switch. When the VTF is correctly connected, the Probe led is on. The working temperature is shown on the display of the VTF
- Select the temperature set point on the VTF
- Push the left knob to start heating. The left display of the AREX 5 Digital shows ---

**Note:** The AREX 5 Digital always has primary control of the heating plate temperature. The temperature control feature of the heating plate can also be used as a safety thermostat. In this case, the maximum temperature of the heating plate will not exceed the temperature setting on the magnetic stirrer (see option ΕΕΠΡ L ΠΕ chap. 4), meaning that a longer heating time is required in order to reach the VTF thermoregulator temperature setting, thus reducing temperature oscillation at a setpoint value.

#### Operation with external probe

- Turn off the instrument using the main switch
- Screw the threaded support rod into its seat on the back of the instrument (optional)
- Fasten the clamp onto the support rod (optional)
- Place the external temperature probe into the clamp (optional); then place it into the sample contained inside the flask
- Plug the probe into the dedicated socket on the back of the instrument. The AREX 5 Digital automatically recognize the type of temperature probe (Pt100 or Pt1000)
- Turn the unit on using the main switch. When the external probe is connected correctly, the Probe led is on
- Turn the left knob to set the working temperature; the Set led indicates that the display is showing the set point temperature
- Push the left knob to start heating. The left display shows the temperature read by the probe for 3 seconds and the setpoint for 1 second

## 4. Setting mode

Press both knobs for 3 seconds when heating and stirring are OFF to access the setting mode.

Turn the left knob to browse through the available menus, listed below.

Press the right knob to see the set parameter. Then, turn the knob to modify the value of the parameter (when it is possible).

Press the knob to confirm.

To exit the setting mode, do not press/turn the knobs for 10 seconds or press both knobs at the same time.

Parameter shown		Default value	Range	Description
Display 1	Display 2			
TEMP	LIMT	310	OFF - 50 - 310 10 °C steps	<u>Heating plate temperature limit</u> It limits the maximum value of the temperature setpoint for the heating plate. ➤ If OFF is selected, the heating function is "disabled".
Start	MODE	STOP	STOP - RUN	<u>Starting mode</u> ➤ STOP: when the instrument is powered, the heating and the stirring functions are not active and both displays show OFF. ➤ RUN: when the instrument is powered, it starts working at the last temperature and speed setpoints.
CTRL	TYPE	FINE	FINE - FAST	<u>Thermoregulation mode with external probe</u> ➤ FINE: optimized temperature control for minimized overshoot and oscillations, slow rise in temperature. ➤ FAST: fast rise in temperature, increased overshoot and oscillations.
PrOb	ALAR	On	On - OFF	<u>External probe safety alarm</u> ➤ It allows enabling/disabling AL 6 and AL 7 error messages (see chap. 5).
PrOb	CAL	0.0 °C	-10 10 °C 0,1 °C steps	<u>Probe calibration</u> ➤ It allows to align the external temperature probe reading to a reference thermometer.
rot	dir	CC	C - CC	<u>Stirring direction:</u> ➤ C: stirring direction is clockwise ➤ CC: stirring direction is counter clockwise
PrOb	CAL	0.0 °C	-10 10 °C 0,1 °C steps	<u>Heating plate probe calibration</u> ➤ It allows to align the heating plate probe reading to a reference thermometer.
COun	MODE	---		<u>Motor and heating element operating time</u> It shows motor and heating element operating times. Operating times are shown in hours up to 9999 hours (around 416 days). Then they are shown in days (a digits' decimal point is turned on).
COun	HEAT			
rSEt		n0	YES - n0	<u>Reset</u> ➤ It allows to restore all system settings to their default values.

## 5. Error and warning messages

When the display shows an error message, the instrument functions automatically stop.

<b>AL1</b>	Heating plate overtemperature (T > 330 °C)	
<b>AL2</b>	Excessive heating time	
<b>AL3</b>	The stirring system doesn't run correctly	
<b>AL4</b>	External probe overtemperature	Only if the external probe is connected
<b>AL5</b>	Heating plate safety probe overtemperature (T > 370 °C)	
<b>AL6</b>	Slow increase in temperature read by the external probe	Only if the external probe is connected and $P_{r_{\text{ob}}}$ is $\square_{r}$
<b>AL7</b>	Fast decrease in temperature read by the external probe	

To remove the error message, restart the instrument from the main switch. For AL5 is necessary to disconnect the instrument from the power supply.

If the alarm persists on the display, please contact Velp Scientifica's technical assistance service.

## 6. Maintenance

No routine or extraordinary maintenance is required except periodic cleaning.

### Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

### Repair

Repairs must be carried out by authorized Velp personnel only.

The transport of the instrument by freight forwarders, couriers or others must be carried out using the original shockproof packaging. Follow any instructions on the original packaging (e.g., palletizing).

It is the responsibility of the user to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

## 7. Technical data

	Model	F20500580	F20510580
	<b>General features</b>	Voltage	230 V – 50/60 Hz
	Power input	630 W	
	Current consumption	2,8 A	5,5 A
	Dimensions (WxHxD)	160x85x270 mm (6,3x3,3x10,6 in)	
	Weight	1,7 kg (3,7 lb)	
	Construction material	Aluminium body – Technopolymer enclosure	
	Working in continuous	Admitted	
	Maximum load on the plate	25 kg (55 lb)	
	Noisiness	<< 80 dBa	
	Environmental temperature admitted	+5...+40 °C	
	Storage temperature admitted	-10...+60 °C	
	Level of electrical protection CEI EN60529	IP 42	
	Max humidity	80%	
	Overvoltage category	II	
	Pollution degree CEI EN61010-1	2	
	Max altitude	4000 m	
<b>Heating</b>	Heating plate power output	600 W	
	Heating plate dimensions	Ø 135 mm	
	Temperature range	0 ÷ 310 °C (1 °C step)	
	Temperature resolution	1 °C	
	Heating plate construction material	CerAlTop™	
	Safety circuit	Separated with dedicated probe	
	Hot plate alarm	Over 50°C	
	Overtemperature alarm	Over 330°C	

<b>External probe</b>	Type	Pt100 – Ø 3mm or Pt1000 – Ø 3mm (optional)
	Temperature range	0 ÷ 310 °C (1 °C step)
	Temperature resolution	1 °C
	Accuracy	± 1 °C *
	Electrical data	3.3 VDC – 1 W (max)
<b>VTF</b>	Temperature range	-10 ÷ 300 °C (1 °C step)
	Temperature resolution	0,2 °C
	Accuracy	+/-0,5°C*
	Electrical data	12Vdc – 1.2W (max)
<b>Stir</b>	Stirring capacity	20 l H <sub>2</sub> O
	Speed range	50 ÷ 1500 rpm (5 rpm step)
	Speed resolution	1 rpm
	Motor type	BLDC
	Motor rating output	10 W
<b>Counters</b>	Motor	Operating hours
	Heating plate	

\* in the following conditions: 800 ml of water in 1 liter glass beaker (diameter 105 mm), stirring bar 8 x 40 mm, 600 rpm, 50°C.

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements
EN 61010-2-010	Particular requirements for laboratory equipment for the heating of material
EN 61010-2-051	Particular requirements for laboratory equipment for mixing and stirring
EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements

and satisfies the essential requirements of the following directives:

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2006/42/EC	Machinery directive
2014/30/EU	Electromagnetic compatibility directive
2015/863/EU (RoHS III)	Restriction of the use of certain hazardous substances
2012/19/EU (WEEE)	Waste of electric and electronic equipment

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung /  
加修改.

## 11. Declaration of conformity <sup>UK</sup>CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use  
Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -  
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic  
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

## Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO 14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

### Analytical instruments

Elemental Analyzers  
Digestion Units  
Distillation Units  
Solvent Extractors  
Fiber Analyzers  
Dietary Fiber Analyzers  
Oxidation Stability Reactor  
Consumables

### Laboratory Equipment

Magnetic Stirrers  
Heating Magnetic Stirrers  
Heating Plates  
Overhead stirrers  
Vortex mixers  
Homogenizers  
COD Thermoreactors  
BOD and Respirometers  
Cooled Incubators  
Flocculators  
Overhead Shakers  
Turbidimeter  
Radiation Detector  
Open Circulating Baths  
Pumps

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VELP Scientifica Srl