Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
指导手册

OV 725 Digital Dispenser

F20900475, F20910475

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 一般信息

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de funcionamiento.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch.
在使用本装置之前，请仔细阅读以下使用说明书。

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2002/96/CE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2002/96/CE.
No tirar el aparato en los desechos urbanos, como exige la Directiva 2002/96/CE.
Dieses Gerät unterliegt der Richtlinie 2002/96/EG und darf nicht mit dem normalen Hausmüll entsorgt werden.
根据 EEC 指令 2002/96/CE，请不要将本设备作为城市垃圾处理。

This unit must be used only for indoor laboratory applications. The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a not specified way by the manufacturer or with specified accessories, product's safety may be compromised.

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This unit has been designed and manufactured in compliance with the following standards:

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**IEC/EN 61010-1
IEC/EN 61010-2-051**

Electrical equipment for laboratory use
General requirement - Canadian electrical code

**UL 61010-1
CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improve their quality.
Nell'impegno di migliorare costantemente la qualità dei prodotti, VELP si riserva la facoltà di variarne le caratteristiche.
Dans le but d'améliorer constamment la qualité de ses produits, VELP se réserve le droit d'apporter des modifications aux caractéristiques de ceux-ci.
VELP se reserva el derecho de modificar las características de sus productos con el objetivo de mejorar constantemente su calidad.
VELP behält sich zum Zwecke der ständigen Verbesserung der Produktqualität das Recht auf Änderung der Geräteeigenschaften vor.
VELP 保留修改其产品特性的权利，目的是不断提高其质量。

Safety Regulations / Norme di Sicurezza / Consignes de Sécurité / Advertencias de Seguridad / Sicherheitshinweise / 安全条例

The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.
La spina è il mezzo di disconnessione dell'apparecchio. Pertanto, non posizionare l'apparecchio in modo che sia difficile azionare il mezzo di disconnessione.
Le bouchon est le moyen de déconnexion de l'appareil. Par conséquent, placer l'appareil où il peut être rapidement débranché.
El tapón es el medio de desconexión del dispositivo. No coloque el dispositivo en una forma que es difícil de desconectar.
Der Stecker trennt das Gerät. Daher Stellen Sie das Instrument, wo es schnell getrennt werden kann.
该插头可以断开仪器的连接。因此，要把仪器放在可以快速断开的地方。

The values indicated on the rating plate of the instrument must correspond to those of the power supply.
I valori di tensione indicato sulla targhetta del modello e quello di rete devono coincidere.
Les valeurs indiquées sur la plaque signalétique de l'appareil doivent correspondre à ceux de l'alimentation.
Los valores de tensión indicados en la placa y que de la red debe ser los mismos.
Die angegebene Spannung Wert auf dem Typenschild und das Netzwerk muss gleich sein.
仪表的铭牌上标明的数值必须与电源的数值一致。

Position the Double Rod Mini H-Stand Set (A00000484) on a flat, non-slip, dry and non-flammable surface at 30 cm (at least) from the walls.
Posizionare il Mini H-Stand con doppia asta completo (A00000484) su una superficie piana, antiscivolo, asciutta e non infiammabile ad una distanza dalle pareti di almeno 30 cm.
Positionner le support mini H (A00000484) sur une surface plane, antidérapant, sec et ininflammable avec une distance de la paroi de 30 cm (au moins).
Coloque el soporte Mini H (A00000484) sobre una superficie plana, antideslizante, seca y no inflamable con una distancia de la pared de 30 cm (por lo menos).
Stellen Sie der H-Stativ (A00000484) auf einer ebenen, rutschfesten, trockenen und nicht brennbaren Oberfläche auf Abstand zur Wand von 30 cm (mindestens).
将双杆迷你 H 型支架套件 (A00000484) 放置在距离墙壁 (至少) 30 厘米的平坦、防滑、干燥且不易燃的表面上。

Fasten the unit to the Double Rod Mini H-Stand Set (A00000484). Secure the receptacle using the vessel holder with safety belt.
Fissare saldamente lo strumento al Mini H-Stand con doppia asta completo (A00000484) e il recipiente di agitazione con l'apposito supporto contenitore con cinghia.
Fixez l'unité à le statif (A00000484). Fixer le réceptacle à l'attache souple avec le support de récipient avec sangle.
Asegurar firmemente la herramienta a la barra de soporte (A00000484) y el recipiente con la Abrazadera y una correa para sujetar recipientes.
Befestigen Sie das Gerät an der H-Stativ (A00000484). Sichern Sie die Aufnahme mit der Spannhalter.
使用双夹子 将仪器固定在支撑杆 (A00000484) 上。用带状夹子 固定插座。

Safe working conditions are ensured only when the accessories described in the dedicated chapter are used. Read their instructions manual carefully.

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Do not use the device without a dispersing tool.

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Before installing the dispersing tools, disconnect the plug from the power supply.

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It is dangerous to run the unit with the dispersing tool turning in free air. Always place the dispersing tool in the receptacle before turning the unit on.

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Never use the dispersing tool dry because tools used without coolant will destroy the gasket.

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The working speed set on the instrument must be such as to avoid wobbling and/or splashes.

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If the noise and vibration of the instrument increase significantly (resonance), change the speed setpoint quickly. This is a physical phenomenon that depends on the configuration in which the disperser operates (dispersing tool, sample type, volume and speed). If this condition persists, check the wear of the dispersing tool.

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Do not use with explosive or dangerous materials for which the equipment is not designed. The disperser must not be used in explosive atmospheres, in bain-marie or to stir harmful liquids prior using protective measures according to the safety standards of the processed products and/or in force in the laboratories including personal protective equipment and the presence of an extraction hood which ensures at least 10-fold air change in accordance with the standards EN 14175 and DIN 12924.

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It is responsibility of the user to appropriately decontaminate the instrument in case of dangerous substances fall on or in it.

It is also responsibility of the user to use safety substances for cleaning or decontaminating, which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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The solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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The OV 725 Digital Dispenser, equipped with an electronic speed control system, a brushless motor and an advanced safety system, it is the ideal solution for dispersing, homogenizing, mixing and crushing samples of biological materials (cells, animal and plant tissues), pharmaceutical, cosmetic and food products.

The term dispersion means the subdivision and distribution of a solid or fluid in a liquid phase that is not completely miscible with the other phase.

The dispersing tool is mainly composed of the rotor and the stator. The high number of revolutions of the rotor allows the tool to axially suck the sample into the dispersion head and then push it radially through the slots between the rotor and the stator. Normally, the processing time is a few minutes to obtain the desired fineness, longer processing times only slightly improve the fineness of the sample but significantly increase its temperature.

| | Description | Function |
|-------|----------------------------------|--|
| A | OV 725 Digital | D20 Drive unit |
| B | Dispersing tool fixing knob | To fix the dispersing tool into dispenser |
| C | Mini H-Stand | To fix the dispenser |
| D | Vessel holder with safety belt | To fix the receptacle with the sample |
| E | Main switch | To turn the dispenser on |
| F | Control knob | Rotate the knob to adjust the speed and the timer when the function is active. Press the knob to start/stop the dispersing process |
| G | 4-digit LED display | To show the speed value (\div 1000 rpm) |
| H | LED light bar | Indicates the remaining time of maintenance for dispersing tool |
| I / L | Timer LEDs | I or L blink when the countdown is active |
| M | Timer key | To show on display the remaining time of dispersing process selected by the control knob (F) |
| N | Dispersing tool maintenance LEDs | Blink when the remaining time of maintenance for dispersing tool is less than 1h |
| O | Dispersing tool maintenance key | Shows on display the remaining time of maintenance for dispersing tool |

2. Assembly and installation

- Unpacking
 - Check the integrity of the unit after unpacking.
- The box includes
 - OV 725 Digital Dispenser (A)
 - Power cable
 - Tool disassembly kit
 - Instruction manual
 - Allen key and allen head screw
 - Support shaft
- Installation
 - Place the Mini H-Stand (C) (accessory) on a surface as described on safety regulations
 - Fix the support shaft to the instrument with the allen head screw
 - Firmly fix the instrument to the Mini H-Stand (C)
 - Secure the receptacle with the sample using the vessel holder with safety belt (D) supplied to the Mini H-Stand
 - Slide the dispersing tool (accessory) into the head of the dispenser up to the stop and tight it with the knob (B)
 - Place the dispersing tool into the sample as indicated in Chapter 5
 - Make sure that the rating value of the instrument corresponds to the one of the power supply
 - Ensure that the socket, with grounding, is compliant with the current safety norms and easy to reach
 - Enter the plug into the socket

IMPORTANT NOTE

- Check that the support shaft is firmly seated. The vibrations may loosen the allen head screw, for safety reasons, be sure to tight it periodically.
- Check that the dispenser is securely mounted in position before use and at regular periods. The position of dispenser can only be changed if the instrument is turned off and disconnected from the power supply.

3. Operation

- Power-up
 - Turn the instrument on through the main switch (E)
 - The display (G) shows the software version, last setpoint value and OFF (if "OFF" is set on STOP, see chapter 4)
- Dispensing
 - Adjust speed set point by turning the speed control knob (F)
 - Click the control knob to start dispensing
 - Speed increases until setpoint achievement
 - A microprocessor ensures constant speed even when the viscosity changes (SpeedServo)
 - Stop the dispensing process by pressing the control knob (F)
- Timer
 - Click the Timer key (M). "MM:SS" appears on display and the LED (L) lights up
 - Rotate speed control knob (F) to set the desired value. Click to confirm
 - The LED lights continuously until the instrument starts working
 - When the instrument is working, the LED blink and the countdown starts
 - It is possible to check the countdown at any time during dispensing by clicking the Timer key (M)
 - When the countdown finishes, "END" is visualized on the display and the dispenser stops.
- Timer modality
 - Click the Timer key (M). "MM:SS" appears on display and the LED (L) lights up
 - Click again Timer key (M). "HH:MM" appears on display and the LED (I) lights up
 - Rotate control knob (F) to set the timer. Click to confirm
 - It's possible to select two scales of timer:
 - MM:SS from 00:00 to 99:59
 - HH:MM from 00:00 to 99:59

- LED light bar
 - LED light bar (H) shows an idea of the remaining time before the dispersing tool needs to be maintained
 - If the set maintenance value is 20:00 (HH:MM), each LED lit means 5 hours before carrying out maintenance
- Maintenance timer for dispersing tool
 - Click the dispersing tool maintenance key (O). The display shows the remaining time before the dispersing tool requires maintenance
 - After 3 seconds the display shows the dispersing status again
 - With the motor in OFF status, keeping pressing the key (O) for 3 seconds, it is possible to reset the dispersing tool maintenance timer. When the display shows “r5Et”, click the control knob (F), the display shows “dOnE” for 3 seconds. In this way, the LED light bar lights up completely and the countdown restarts from the set time value (“[RL” value, see Chapter 4). This operation must be carried out once maintenance of the dispersing tool has been completed: changing the gasket and/or replacing any worn part
 - When the remaining time is less than 1 h, the LEDs (N) flash
 - When the remaining time arrives at 00:00, the LED light bar flashes. The disperser continues to work.

Note: the frequency of maintenance of the dispersing tool can be set as described in Chapter 4 according to the user's application requirements.

4. Setting mode

When the motor is stopped (OFF status), hold the control knob (F) pressed for 3 seconds to enter into the setting mode. By turning the control knob (F) it is possible to move from one of the following parameters to the next one. Click the knob to modify:


| Parameter shown | Default value | Range | Description |
|-----------------|---------------|---------------------|---|
| mode | STOP | STOP - run | Set starting mode: <ul style="list-style-type: none"> ➤ STOP: when the instrument is switched on through the main switch, OFF is on display ➤ run: when the instrument is switched on, it restarts to work with the last speed set and the remaining time |
| limit | 25000 | 2000 - 25000 | Set the maximum speed value (rpm) |
| count | | | Running time in hours. If “.” is shown before the number, the running time is in days. |
| reset | no | YES - no | To restore default setting parameters |
| RL | 20:00 | 02:00 - 50:00 - OFF | To modify the frequency of maintenance for the dispersing tool (hh:mm). Less than 02:00 it's possible to select OFF, in this way, the timer is never active. |

To leave settings mode, wait 10 seconds without touching the knob or click the maintenance key (O).

5. Dispersing tools

Insertion of the dispersing tool into the disperser

Unscrew the knob (B) until its head is no more visible through the hole. Slide the dispersing tool from the bottom into the disperser, up to the stop. The dispersing tool enters with an audible click if you insert it with enough force. It's fitted correctly when the notch on the dispersing tool is aligned with the bottom of the disperser.

Now fix the dispersing tool by tightening the knob (B).

If the dispersing tool is not inserted into disperser and the control knob is pressed to start the dispersion, the display shows “-00-” “E00L” and the unit does not start.

Note: before turning on the unit, make sure that the dispersing tool is inserted in the liquid and that the receptacle is well fixed with the vessel holder with safety belt (D).

ABBREVIATIONS

All D20 series tools are allowed (it means with a 20mm diameter head).

Material: stainless steel AISI 316L.

Head of disperser:

D20 For drive unit with a head of 20mm as diameter

Diameter and type of stator:

S20C Coarse stator diam. 20 mm

S20F Fine stator diam. 20 mm

S25C Coarse stator diam. 25 mm

S25F Fine stator diam. 25 mm

S40C Coarse stator diam. 40 mm (only for mixtures up to 20.000 rpm)

Type of gasket:

P PTFE

Diameter and type of rotor:

R20S Standard rotor diam. 20 mm

R20C Coarse rotor diam. 20 mm

R20E Emulsifying rotor diam. 20 mm

R20M Mixing rotor diam. 20 mm

R25C Coarse rotor diam. 25 mm

R25E Emulsifying rotor diam. 25 mm

R25M Mixing rotor diam. 25 mm

Example of dispersing tool configuration

D20-S20C-P-R20S Velp code: A00000478

AVAILABLE DISPERSING TOOLS

See all dispersing tools (accessories) in Chapter 9 or go to www.velp.com.

Use only the original Velp Scientifica dispersing tools after reading their instruction manual with more information about applications and technical data.

TO OBTAIN BETTER RESULTS

Diameter of the receptacle with the sample

It is better to disperse using a container slightly larger than the dispersing tool in order to obtain a faster and more homogeneous solution.

This will prevent part of the sample from remaining on the walls of the container and therefore not being homogenized. Keep in mind that the insertion of the tool into the receptacle causes the liquid level to rise.

Position of the dispersing tool

- Place the dispersing tool not in the middle of the receptacle (Figure 1)
- Place the dispersing tool about 1/3 of the sample height from the bottom of the receptacle (Figure 2).

Figure 1

Figure 2

If you use a dispersing tool that is too large for the sample amount, it could suck air.

Dispersing speed

Always start at a low speed, then gradually increase it. This allows to keep control of the process.

Normally the dispersion occurs in 30 – 60 seconds, if this does not happen, check:

- The dispersing tool position in the container
- If the sample is too viscous
- Incorrect dispersing tool for the sample (size, type)

TO NOT REDUCE THE LIFETIME OF THE DISPERSING TOOL AND THE DISPERSER

Never run the dispersing tool without liquid

The dispersing tool requires a liquid sample for the lubrication of the gasket during machining.

A PTFE gasket is placed between the rotor and stator to prevent the rotor from bending or moving inside the tool. The gasket must be immersed in the liquid sample, otherwise it will burn causing overheating and subsequent damage to the sample and the dispersing tool.

Without the liquid, the PTFE gasket can burn causing damage to the tool and the disperser.

Cleaning and maintenance of dispersing tool

It's necessary to disassemble and wash the dispersing tool after each use, checking the good condition of the gasket and the metal parts and replacing them if necessary. The PTFE gasket is naturally subjected to wear because of the friction with the rotating elements of the tool; therefore it's necessary to periodically replace it in order to guarantee the correct operation of the tool and the dispersing instrument.

Replace the PTFE gasket if the following conditions occur:

- Visible wear of the component (e.g. deformations, thinning),
- Unusual noise and vibration during the operation of the dispersing instrument,
- When the space between the gasket and the internal and external metal parts is excessive.

To disassemble the dispersing tool, use the appropriate tool disassembly kit supplied with the drive unit and read the instruction manual supplied with each tool.

6. Error and warning messages

When the display shows an error message, the motor stops automatically. To restore the functionality of the instrument, proceed as follows:

- Turn the unit off using the main switch
- Use the suggestions described in the column "Solution"
- Turn the unit on

| Code | Cause | Solution |
|------|--|--|
| RL1 | The motor doesn't start or not reach the set speed | <ul style="list-style-type: none">• Check that there are no locks in the dispersing tool or structure |
| RL2 | High internal motor temperature | <ul style="list-style-type: none">• Check that there are no obstructions on the air vents for internal ventilation.• Check that there are no locks in the dispersing tool or structure• Reduce speed or load |
| RL3 | Motor overload | <ul style="list-style-type: none">• Check that there are no locks in the dispersing tool or structure• Reduce speed or load |
| RL4 | High internal temperature | <ul style="list-style-type: none">• Check that there are no obstructions on the air vents for internal ventilation.• Check that there are no locks in the dispersing tool or structure• Reduce speed or load |
| RL5 | Internal hardware error | <ul style="list-style-type: none">• Contact service department |
| RL6 | The instrument does not detect the dispersing tool | <ul style="list-style-type: none">• Check that the dispersing tool is correctly inserted in the disperser and firmly fixed with the knob (B) |
| RL8 | Speed variation too high | <ul style="list-style-type: none">• Check that there are no locks in the dispersing tool or structure |

If an error persists on the display, please contact VELP Scientifica's technical service department.

7. Maintenance and cleaning

No routine or extraordinary maintenance is necessary for the unit. The maintenance is foreseen only on the dispersing tools as described in Chapter 5.

Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent. During cleaning operations, prevent water and humidity residues from entering the unit.

Repair

Repairs must be carried out by authorized Velp personnel only.

The transport of the instrument by freight forwarders, couriers or others must be carried out using the original shockproof packaging. Follow any instructions on the original packaging (e.g., palletizing).

It is the responsibility of the user to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

8. Technical data

| | Models | F20900475 | F20910475 |
|--------------------------|---|---|---------------------------------|
| General | Power supply | 230 V / 50-60 Hz (+/-10%) | 115V / 60 Hz (+/-10%) |
| | Dimensions (WXHxD): | 100x280x235 mm (3.9x11x9.3 in) | |
| | Weight | 3,7 kg (8,1 lb) | |
| | Power consumption | 500 W | |
| | Construction material | Aluminum - Technopolymer | |
| | Noise (without element) | 70 dBa | |
| | Environmental temperature admitted | +5...+40 °C | |
| | Storage temperature admitted | -10...+60 °C | |
| | Max humidity | 80% | |
| | Overvoltage category | II | |
| | Pollution degree CEI EN61010-1 | 2 | |
| | Max altitude | 4.000 m | |
| | Dispersing | Volume range | 0,2 ml – 2,5 l H ₂ O |
| Programmable speed range | | 2.000 – 25.000 rpm (100 rpm step) (at nominal voltage 230V 50/60Hz and 115V 60 Hz) | |
| Motor type | | BLDC | |
| Speed selection | | 100 rpm step | |
| Motor rating output | | 400 W | |
| Maximum viscosity | | 5.000 mPa s | |
| Safety | Alarm if dispersing tool is not present | Yes | |
| | Alarm inside overtemperature | Yes | |
| | Alarm if motor can't start or not reach the speed set | Yes | |
| Timer | Working time | mm:ss | 00:00 ...99:59 |
| | | hh:mm | 00:00...99:59 |
| Counters | Dispersing tool maintenance | hh:mm of operation | |
| | Motor counter | Hours or days of operation | |

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of the following directives:

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2006/42/EC Machinery directive

2014/30/EU Electromagnetic compatibility directive

2015/863/EU (RoHS III) Restriction of the use of certain hazardous substances

2012/19/EU (WEEE) Waste of electric and electronic equipment

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung /
加修改.

12. Declaration of conformity ^{UK}_{CA}

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity.

VELP works according to **ISO 9001**, **ISO 14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortexers and Shakers
Dispersers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Turbidimeter
Open Circulating Baths
Pumps

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VELP Scientifica Srl