**Instruction Manual**  
**Manuale di istruzioni**  
**Manuel d'instructions**  
**Manual de instrucciones**  
**Bedienungsanleitung**  
**操作指南**

## AREC 4 Digital Heating Magnetic Stirrer

F20500670, F20510671

## AREC 10 Digital Heating Magnetic Stirrer

F20500691, F20510691, F20580691

## REC 10 Digital Heating Plate

F20700693, F20710693, F20780693

### General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 基本信息

Before using the unit, please read the following instruction manual carefully.

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Attention, surface chaude!

Prudencia, superficie caliente!

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2012/19/UE.

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The product can be used with flammable liquids.

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**This unit must be used for indoor laboratory applications only.**

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a way not specified by the manufacturer or with non specified accessories, product's safety may be compromised.

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**This unit has been designed and manufactured in compliance with the following standards:**

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**IEC/EN 61010-1  
IEC/EN 61010-2-051  
IEC/EN 61010-2-010**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

**CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improve their quality.

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The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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Use only the power cable provided with the instrument.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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The vessel must be made of a suitable material to withstand the foreseen temperature.

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Magnetic field doesn't interfere with the function of cardiac pacemakers or data media if they are farther than 20 cm to the instrument.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least) and between instruments of 10 cm (at least).

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Do not use with explosive or dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie or to stir combustible liquids that have a low combustion temperature. The product is intended for use with very small quantities of flammable liquids or flammable liquids that have a fire point higher than 625°C and a flash point higher than 600°C.

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It is responsibility of the user to appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories' safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts about the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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The AREC 4 / AREC 10 Digital heating magnetic stirrer and the REC 10 Digital heating plate is used to heat and mix (not available on REC 10 Digital) liquids inside a suitable container placed on the ceramic plate. The stirring (not available on REC 10 Digital) is possible thanks to a magnetic drive stir bar placed into the vessel.

The instrument is used in the laboratory for general use and for all those applications that require precise regulation of the stirring speed (not available on REC 10 Digital) and of the heating plate temperature.

- 1 Ceramic plate
- 2 Main switch
- 3 Speed control knob (not available on REC 10 Digital)
- 4 Temperature control knob
- 5 Display

The ceramic heating plate allows to reach higher temperature than a traditional aluminium heating plate; furthermore, the ceramic is an inert and very hard material, so it offers great resistance to almost any type of chemical or mechanical aggression keeping the surface characteristics unaltered over time.

The heating plate works up to 550 °C.

The powerful brushless motor (not available on REC 10) is able to stir from 100 up to 1500 rpm. SpeedServo™ technology maintains constant speed as viscosity changes.

**Note (not for REC 10 Digital):** it's important to choose the most suitable magnetic stirring bar related to the type and the quantity of the liquid to be stirred as well as to the type of the beaker. Size and shape of the magnetic stir bar determines the stirring efficiency at any given speed. The stirring bar which satisfies most applications is code A00000356 (8 x 40 mm).

The AREC 4 / AREC 10 Digital / REC 10 Digital can work in combination with the Pt100 or Pt1000 external probe or with thermoregulator VTF for precise and accurate control of the sample temperature.

## 2. Installation

- Unpacking
  - Check the integrity of the unit after unpacking
- The box includes
  - AREC 4 / AREC 10 Digital heating magnetic stirrer / REC 10 Digital heating plate
  - Instruction manual
  - Power cable
  - Magnetic stir bar Ø 8 x 40 mm (A00000356) (not for REC 10 Digital)
- Installation
  - Place the unit on a non-flammable surface
  - Make sure that the technical specification of the unit and of the power supply line are the same
  - Make sure that the electrical network is grounded
  - Make sure that the main switch is on position "0" (OFF)
  - Connect the instrument to an easily accessible socket (compliant with the current safety norms), using only the provided power cable

**Note:** keep the power cable far away from the hot plate.

  - Place the container with the liquid and the stir bar (not for REC 10 Digital) on the stirring plate

## 3. Operation

### Power-up

- Turn the unit on using the main switch
- The display shows the software version and the last setpoint values

If the set starting mode is stop (see chap. 4), the heating and the stirring (not available on REC 10 Digital) functions are not active and both displays show OFF. Otherwise, the unit starts working at the last set temperature and speed setpoints.

### Stirring (not available on REC 10 Digital)

- Rotate the right knob to set the stirring speed
- Push the knob to start stirring
- The real speed increases until set point achievement
- A microprocessor ensures constant speed even if the viscosity changes (SpeedServo™)
- Push the knob to stop stirring

### Heating

- Rotate the left knob to set the heating plate temperature
- Push the knob to start heating
- The display shows the current heating plate temperature for 3 seconds and the setpoint for 1 second
- When the display shows the temperature setpoint, the Set led is on
- Push the knob to stop heating

**Note:** if the heating function is not active and the heating plate temperature is higher than 50 °C, the display shows the blinking message "Hot". This warning is active even if the main switch is off but not if the instrument is not powered.

### Operation with thermoregulator VTF

- Turn off the instrument using the main switch
- Screw the threaded support rod into its seat on the back of the instrument and fasten the VTF thermoregulator onto the support rod. Place the temperature probe in the flask making sure that it is completely immersed in the sample
- Plug the VTF into the dedicated socket on the back of the instrument
- Turn the unit on using the main switch. When the VTF is correctly connected, the Probe led is on. The working temperature is shown on the display of the VTF
- Select the temperature set point on the VTF
- Push the left knob to start heating. The left display of the AREC Digital shows ---

**Note:** The AREC Digital always has primary control of the heating plate temperature. The temperature control feature of the heating plate can also be used as a safety thermostat. In this case, the maximum temperature of the heating plate will not exceed the temperature setting on the magnetic stirrer (see option tENP LINT chap. 4), meaning that a longer heating time is required in order to reach the VTF thermoregulator temperature setting, thus reducing temperature oscillation at a setpoint value.

#### Operation with external probe

- Turn off the instrument using the main switch
- Screw the threaded support rod into its seat on the back of the instrument (optional)
- Fasten the clamp onto the support rod (optional)
- Place the external temperature probe into the clamp (optional); then place it into the sample contained inside the flask
- Plug the probe into the dedicated socket on the back of the instrument. The AREC4 / AREC 10 / REC 10 Digital automatically recognizes the type of temperature probe (Pt100 or Pt1000)
- Turn the unit on using the main switch. When the external probe is connected correctly, the Probe led is on
- Turn the left knob to set the working temperature; the Set led indicates that the display is showing the set point temperature
- Push the left knob to start heating. The left display shows the temperature read by the probe for 3 seconds and the setpoint for 1 second

## 4. Setting mode

Press both knobs for 3 seconds when heating and stirring are OFF to access the setting mode.

Turn the left knob to browse through the available menus, listed below.

Press the right knob to see the set parameter. Then, turn the knob to modify the value of the parameter (when it is possible). Press the knob to confirm.

To exit the setting mode, do not press/turn the knobs for 10 seconds or press both knobs at the same time.

For REC 10 Digital, follow the procedure using the single knob.

| Parameter shown |           | Default value | Range                         | Description  |
|-----------------|-----------|---------------|-------------------------------|--|
| Display 1       | Display 2 |               |                               |  |
| tENP            | LINt      | 550           | OFF - 50 – 550<br>10 °C steps | <u>Heating plate temperature limit</u><br>It limits the maximum value of the temperature setpoint for the heating plate.<br>➤ If OFF is selected, the heating function is “disabled”.  |
| Strt            | NOdE      | StOP          | StOP - run                    | <u>Starting mode</u><br>➤ StOP: when the instrument is powered, the heating and the stirring functions are not active and both displays show OFF.<br>➤ run: when the instrument is powered, it starts working at the last temperature and speed setpoints.                           |
| Ctrl            | tyPE      | FlnE          | FlnE - FAST                   | <u>Thermoregulation mode with external probe</u><br>➤ FlnE: optimized temperature control for minimized overshoot and oscillations, slow rise in temperature.<br>➤ FAST: fast rise in temperature, increased overshoot and oscillations.   |
| PrOb            | ALAr      | On            | On - OFF                      | <u>External probe safety alarm</u><br>➤ It allows enabling/disabling AL 6 and AL 7 error messages (see chap. 5).   |
| PrOb            | CAL       | 0.0 °C        | -10 ÷ 10 °C<br>0,1 °C steps   | <u>Probe calibration</u><br>➤ It allows to align the external temperature probe reading to a reference thermometer.  |
| rOt             | dIr       | CC            | C - CC                        | <u>Stirring (not available for REC 10 Digital) direction:</u><br>➤ C: stirring direction is clockwise<br>➤ CC: stirring direction is counter clockwise   |
| ANbt            | CAL       | 0.0 °C        | -10 ÷ 10 °C<br>0,1 °C steps   | <u>Heating plate probe calibration</u><br>➤ It allows to align the heating plate probe reading to a reference thermometer.   |
| COun            | NOtO      | ---           |                               | <u>Motor (not available for REC 10 Digital) and heating element operating time</u><br>It shows motor and heating element operating times. Operating times are shown in hours up to 9999 hours (around 416 days). Then they are shown in days (a digits' decimal point is turned on). |
| COun            | HEAt      |               |                               |  |
| rSEt            |           | nO            | yES - nO                      | <u>Reset</u><br>➤ It allows to restore all system settings to their default values.  |

## 5. Error and warning messages

When the display shows an error message, the instrument functions automatically stop.

|            |   |                                 |
|------------|---|---------------------------------|
| <b>AL1</b> | Heating plate overtemperature (T > 570 °C)              |                                 |
| <b>AL2</b> | Excessive heating time                                  |                                 |
| <b>AL3</b> | The stirring system doesn't run correctly               | Not available on REC 10 Digital |
| <b>AL4</b> | External probe overtemperature                          |                                 |
|            |   |                                 |
| <b>AL6</b> | Slow increase in temperature read by the external probe |                                 |
| <b>AL7</b> | Fast decrease in temperature read by the external probe |                                 |

To remove the error message, restart the instrument from the main switch.

If the alarm persists on the display, please contact VELP Scientifica's technical assistance service.

## 6. Maintenance

No routine or extraordinary maintenance is required except periodic cleaning.

### Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

### Repair

Repairs must be carried out by authorized Velp personnel only.

The transport of the instrument by freight forwarders, couriers or others must be carried out using the original shockproof packaging. Follow any instructions on the original packaging (e.g., palletizing).

It is the responsibility of the user to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

## 7. Technical data

| AREC 4 DIGITAL          |                                     |  |               |
|-------------------------|-------------------------------------|--|---------------|
| <b>General features</b> | Model                               | F20500671  | F20510671     |
|                         | Voltage                             | 230 V – 50/60 Hz                                 | 115 V – 60 Hz |
|                         | Power input                         | 280 W  |               |
|                         | Current consumption                 | 1,2 A  | 2,4 A         |
|                         | Dimensions (WxHxD)                  | 160x100x270 mm (6,3x3,9x10,6 in)                 |               |
|                         | Weight                              | 1,9 kg (4,2 lb)                                  |               |
|                         | Construction material               | Aluminium body – Technopolymer enclosure         |               |
|                         | Working in continuous               | Admitted   |               |
|                         | Maximum load on the plate           | 25 kg (55 lb)                                    |               |
|                         | Noisiness                           | << 80 dBa  |               |
|                         | Environmental temperature admitted  | +5...+40 °C                                      |               |
|                         | Storage temperature admitted        | -10...+60 °C                                     |               |
|                         | Max humidity                        | 80%  |               |
|                         | Overvoltage category                | II   |               |
|                         | Pollution degree CEI EN61010-1      | 2  |               |
| Max altitude            | 4000 m                              |  |               |
| <b>Heating</b>          | Heating plate power output          | 250 W  |               |
|                         | Heating plate dimensions            | 100 x 100 mm (4"x4")                             |               |
|                         | Temperature range                   | 0 ÷ 550 °C (1 °C step)                           |               |
|                         | Temperature resolution              | 1 °C   |               |
|                         | Heating plate construction material | Ceramic  |               |
|                         | Hot plate alarm                     | Over 50°C  |               |
|                         | Overtemperature alarm               | Over 570°C                                       |               |
| <b>External probe</b>   | Type                                | Pt100 – Ø 3mm<br>or<br>Pt1000 – Ø 3mm (optional) |               |
|                         | Temperature range                   | 0 ÷ 350 °C (1 °C step)                           |               |
|                         | Temperature resolution              | 1 °C   |               |
|                         | Accuracy                            | ± 1 °C *   |               |
|                         | Electrical data                     | 3.3 VDC – 1 W (max)                              |               |
| <b>VTF</b>              | Temperature range                   | -10 ÷ 300 °C (1 °C step)                         |               |
|                         | Temperature resolution              | 0,2 °C   |               |
|                         | Accuracy                            | +/-0,5°C*  |               |
|                         | Electrical data                     | 12Vdc – 1.2W (max)                               |               |
| <b>Stir</b>             | Stirring capacity                   | 5 l H <sub>2</sub> O                             |               |
|                         | Speed range                         | 100 ÷ 1500 rpm (5 rpm step)                      |               |
|                         | Speed resolution                    | 1 rpm  |               |
|                         | Motor type                          | BLDC   |               |
|                         | Motor rating output                 | 10 W   |               |
| <b>Counters</b>         | Motor (not for REC 10)              | Operating hours                                  |               |
|                         | Heating plate                       |  |               |

\* in the following conditions: 800 ml of water in 1 liter glass beaker (diameter 105 mm), stirring bar 8 x 40 mm, 600 rpm, 50°C.

| AREC 10 DIGITAL         |                                     |  |               |                  |
|-------------------------|-------------------------------------|--|---------------|------------------|
| <b>General features</b> | Model                               | F20500691  | F20510691     | F20580691        |
|                         | Voltage                             | 230 V – 50/60 Hz                                 | 115 V – 60 Hz | 100 V – 50/60 Hz |
|                         | Power input                         | 1530 W   | 1520 W        | 1320 W           |
|                         | Current consumption                 | 6,7 A  | 13,2 A        | 13,2             |
|                         | Dimensions (WxHxD)                  | 280x100x420 mm (11x3,9x16,5 in)                  |               |                  |
|                         | Weight                              | 5,6 kg (12,4 lb)                                 |               |                  |
|                         | Construction material               | Aluminium body – Technopolymer enclosure         |               |                  |
|                         | Working in continuous               | Admitted   |               |                  |
|                         | Maximum load on the plate           | 30 kg (66 lb)                                    |               |                  |
|                         | Noisiness                           | << 80 dBa  |               |                  |
|                         | Environmental temperature admitted  | +5...+40 °C                                      |               |                  |
|                         | Storage temperature admitted        | -10...+60 °C                                     |               |                  |
|                         | Max humidity                        | 80%  |               |                  |
|                         | Overvoltage category                | II   |               |                  |
|                         | Pollution degree CEI EN61010-1      | 2  |               |                  |
|                         | Max altitude                        | 4000 m   |               |                  |
| <b>Heating</b>          | Heating plate power output          | 1500 W   | 1490 W        | 1290 W           |
|                         | Heating plate dimensions            | 260 x 260 mm (10"x10")                           |               |                  |
|                         | Temperature range                   | 0 ÷ 550 °C (1 °C step)                           |               |                  |
|                         | Temperature resolution              | 1 °C   |               |                  |
|                         | Heating plate construction material | Ceramic  |               |                  |
|                         | Hot plate alarm                     | Over 50°C  |               |                  |
|                         | Overtemperature alarm               | Over 570°C                                       |               |                  |
| <b>External probe</b>   | Type                                | Pt100 – Ø 3mm<br>or<br>Pt1000 – Ø 3mm (optional) |               |                  |
|                         | Temperature range                   | 0 ÷ 350 °C (1 °C step)                           |               |                  |
|                         | Temperature resolution              | 1 °C   |               |                  |
|                         | Accuracy                            | ± 1 °C *   |               |                  |
|                         | Electrical data                     | 3.3 VDC – 1 W (max)                              |               |                  |
| <b>VTF</b>              | Temperature range                   | -10 ÷ 300 °C (1 °C step)                         |               |                  |
|                         | Temperature resolution              | 0,2 °C   |               |                  |
|                         | Accuracy                            | +/-0,5°C*  |               |                  |
|                         | Electrical data                     | 12Vdc – 1.2W (max)                               |               |                  |
| <b>Stir</b>             | Stirring capacity                   | 25 l H <sub>2</sub> O                            |               |                  |
|                         | Speed range                         | 100 ÷ 1500 rpm (5 rpm step)                      |               |                  |
|                         | Speed resolution                    | 1 rpm  |               |                  |
|                         | Motor type                          | BLDC   |               |                  |
|                         | Motor rating output                 | 10 W   |               |                  |
| <b>Counters</b>         | Motor (not for REC 10)              | Operating hours                                  |               |                  |
|                         | Heating plate                       |  |               |                  |

\* in the following conditions: 800 ml of water in 1 liter glass beaker (diameter 105 mm), stirring bar 8 x 40 mm, 600 rpm, 50°C.

| REC 10 DIGITAL          |                                     |  |               |                  |
|-------------------------|-------------------------------------|--|---------------|------------------|
| <b>General features</b> | Model                               | F20700693  | F20710693     | F20780693        |
|                         | Voltage                             | 230 V – 50/60 Hz                                 | 115 V – 60 Hz | 100 V – 50/60 Hz |
|                         | Power input                         | 1510 W   | 1500 W        | 1300 W           |
|                         | Current consumption                 | 6,6 A  | 13 A          | 13               |
|                         | Dimensions (WxHxD)                  | 280x100x420 mm (11x3,9x16,5 in)                  |               |                  |
|                         | Weight                              | 5,3 kg (11,6 lb)                                 |               |                  |
|                         | Construction material               | Aluminium body – Technopolymer enclosure         |               |                  |
|                         | Working in continuous               | Admitted   |               |                  |
|                         | Maximum load on the plate           | 30 kg (66 lb)                                    |               |                  |
|                         | Noisiness                           | << 80 dBa  |               |                  |
|                         | Environmental temperature admitted  | +5...+40 °C                                      |               |                  |
|                         | Storage temperature admitted        | -10...+60 °C                                     |               |                  |
|                         | Max humidity                        | 80%  |               |                  |
|                         | Overvoltage category                | II   |               |                  |
|                         | Pollution degree CEI EN61010-1      | 2  |               |                  |
|                         | Max altitude                        | 4000 m   |               |                  |
| <b>Heating</b>          | Heating plate power output          | 1500 W   | 1490 W        | 1290 W           |
|                         | Heating plate dimensions            | 260 x 260 mm (10"x10")                           |               |                  |
|                         | Temperature range                   | 0 ÷ 550 °C (1 °C step)                           |               |                  |
|                         | Temperature resolution              | 1 °C   |               |                  |
|                         | Heating plate construction material | Ceramic  |               |                  |
|                         | Hot plate alarm                     | Over 50°C  |               |                  |
|                         | Overtemperature alarm               | Over 570°C                                       |               |                  |
| <b>External probe</b>   | Type                                | Pt100 – Ø 3mm<br>or<br>Pt1000 – Ø 3mm (optional) |               |                  |
|                         | Temperature range                   | 0 ÷ 350 °C (1 °C step)                           |               |                  |
|                         | Temperature resolution              | 1 °C   |               |                  |
|                         | Accuracy                            | ± 2 °C *   |               |                  |
|                         | Electrical data                     | 3.3 VDC – 1 W (max)                              |               |                  |
| <b>VTF</b>              | Temperature range                   | -10 ÷ 300 °C (1 °C step)                         |               |                  |
|                         | Temperature resolution              | 0,2 °C   |               |                  |
|                         | Accuracy                            | +/-1°C*  |               |                  |
|                         | Electrical data                     | 12Vdc – 1.2W (max)                               |               |                  |
|                         | Heating plate                       | Operating hours                                  |               |                  |

\* in the following conditions: 800 ml of water in 1 liter glass beaker (diameter 105 mm), 50°C.

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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|                |   |
|----------------|---|
| EN 61010-1     | Safety requirements for electrical equipment for measurement, control, and laboratory use<br>Part 1: General requirements |
| EN 61010-2-010 | Particular requirements for laboratory equipment for the heating of material  |
| EN 61010-2-051 | Particular requirements for laboratory equipment for mixing and stirring  |
| EN 61326-1     | Electrical equipment for measurement, control and laboratory use - EMC requirements -<br>Part 1: General requirements     |

and satisfies the essential requirements of the following directives:

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|                        |  |
|------------------------|--|
| 2006/42/EC             | Machinery directive                                    |
| 2014/30/EU             | Electromagnetic compatibility directive                |
| 2015/863/EU (RoHS III) | Restriction of the use of certain hazardous substances |
| 2012/19/EU (WEEE)      | Waste of electric and electronic equipment             |

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung /  
加修改.

#### 4. Declaration of conformity <sup>UK</sup>CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use  
Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -  
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic  
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

## Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO 14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

### Analytical instruments

Elemental Analyzers  
Digestion Units  
Distillation Units  
Solvent Extractors  
Fiber Analyzers  
Dietary Fiber Analyzers  
Oxidation Stability Reactor  
Consumables

### Laboratory Equipment

Magnetic Stirrers  
Heating Magnetic Stirrers  
Heating Plates  
Overhead stirrers  
Vortex mixers  
Dispersers  
COD Thermoreactors  
BOD and Respirometers  
Cooled Incubators  
Flocculators  
Overhead Shakers  
Turbidimeter  
Open Circulating Baths  
Pumps

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VELP Scientifica Srl