**EN - Instruction Manual**  
**IT - Manuale di istruzioni**  
**FR - Manuel d'instructions**  
**ES - Manual de instrucciones**  
**DE - Bedienungsanleitung**  
**CN - 操作指南**

## **MULTI-TX5 Multiposition Digital Vortex**

F202A0460

**General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 基本信息**

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2012/19/UE.

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**This unit must be used for laboratory applications indoor only.**

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a not specified way by the manufacturer or with not specified accessories, product's safety may be compromised.

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**This unit has been designed and manufactured in compliance with the following standards:**

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**IEC/EN 61010-1  
IEC/EN 61010-2-051**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

**CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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Use only the power supply cable provided with the instrument.

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The solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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Before placing the test tubes inside the unit, make sure that their caps are tightly closed to prevent sample leakage.

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Do not lift the instrument by the closing plate, shafts or bottom tray. The closing plate is removable.

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If the shaking movement of the device increases (resonance), reduce the speed or pass through the critical phase as quickly as possible. This is a physical phenomenon that depends on the load and speed selected.

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This unit is designed for intermittent use only. Do not operate with the unit for more than two hours continuously.

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This unit is designed to work automatically without the operator holding the test tube, in order to avoid hand/arm vibrations. Do not touch the instrument during vibration

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Do not use with explosive and dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres.

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The unit is fitted with two fuses (2xT1.6A L 250 V), found in the socket on the back. To replace one or both disconnect the mains cable and, using a screwdriver, lift up the small cover on the fuse box.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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

The MULTI-TX5 Multiposition Digital Vortex is the ideal solution for hands-free and high-throughput sample processing. The instrument is used in the laboratory for general use and applications requiring the precise regulation of the stirring speed. The powerful brushless motor stirs from 100 up to 2500 rpm.

Different interchangeable foam racks permit to work with different test tubes, from 10 to 29 mm in diameter. Powerful feet reduce motion and sliding on the laboratory bench leading to safe work conditions.

The instrument is equipped with advanced functions as auto-reverse of the stirring direction, intermittent mode programs, pulse mode and integrated timer to meet users' workflow. The safety lock feature avoids any possible alterations of the working conditions.

- 1 Closing plate
- 2 Sample rack for test tube
- 3 Display
- 4 Timer control knob
- 5 Speed control knob
- 6 Main switch

- A Intermittent mode led
- B Pulse mode led
- C Timer led
- D Menu and lock button
- E Autoreverse led

## 2. EN - Installation

- Unpacking
  - Check the integrity of the unit after unpacking
- The box includes
  - MULTI-TX5 Multiposition Digital Vortex
  - Instruction manual
  - Power supply cable (EU and US type)
  - Sample rack for test tube Ø12, 44 pos.
- Installation
  - Place the unit on a flat, stable, clean, non-slip, dry and non-flammable surface.
  - The feet of the device must be cleaned and undamaged
  - Make sure that the technical specifications of the unit and the power supply line are the same
  - Make sure that the electrical network is grounded
  - Make sure that the main switch is on position "0" (OFF)
  - Connect the instrument to an easily accessible socket (compliant with the current safety norms), using only the provided power supply cable

**Note:** keep the power supply cable far away from moving parts of the unit.

## 3. EN - Operation

### Power-up

- Turn the unit on using the main switch
- The display shows the software version and the last speed setpoint values

**Note:** If the set starting mode is *StOP* (see chap. 4), the stirring function is not active and the display shows OFF. Otherwise, the unit starts working at the last set speed setpoint.

### Preparing the foam rack with the test tubes

- Fill the test tubes. Best results will be obtained if all test tubes are filled to the same level and do not exceed half of the total capacity.
- Place test tubes in a foam rack of correct size and distribute them uniformly. Always be sure that a test tube is in each of the four corner holes on the foam rack, in order to support the closing plate. Also, empty test tubes may be used for this.

### Placing test tube rack in the unit

- Loosen both screw knobs on the closing plate (1)
- Lift the plate up, keeping it near the metal shafts. The closing plate can be removed or held on the longer shaft. To do this, lift until the closing plate comes off the short shaft. In this condition, the closing plate may be swung around the longer shaft. Tighten the screw knob.
- Place the foam rack inside the foam support.
- Lower the closing plate onto both shafts and down on top of the test tubes. Let the closing plate's weight rest on the tubes, without pressing down otherwise the motor could be damaged.
- Tighten both screw knobs securely.

**Note:** To enlarge the surface of the rack support, two fillers are delivered with the unit. They are to be inserted into the upper and lower rack support. Do not use adhesive if this necessity is only temporary.

### Stirring

- Rotate the right knob to set the stirring speed
- Push the knob to start stirring
- The display shows the current motor speed
- Push the knob to stop stirring

**Note:** If the shaking movement of the device increases (resonance), reduce the speed or pass through the critical phase as quickly as possible.

This unit is designed for intermittent use only. Do not operate with the unit for more than two hours continuously. This unit is designed to work automatically without the operator holding the test tube, in order to avoid hand/arm vibrations.

## Timer

- Rotate the left knob to set the timer
- Push the knob to confirm
- The display shows the timer countdown when stirring starts.
- Select hh.mm to remove timer function.

**Note:** If the timer is not selected, the display will show the count-up of time.

**Note:** when the timer is set, the dedicated led (C) is on. When the countdown starts, the led is blinking.

After the timer has elapsed, the stirring will be stopped, the stirring display will show End and the timer display "time". Three acoustic signals will be emitted.

## Autoreverse, intermittent mode and pulse functions

Press the central button repeatedly to access in sequence the following function menus:

- Autoreverse (display shows Ar Eu)
  - Rotate the right knob to set the autoreverse time (minutes/seconds)
  - Press it to confirm

**Note:** when the autoreverse is set, the dedicated led (E) is on.

**Note:** select nnSS on the display to disable the function.

The rotation direction will change after each set interval time.

- Intermittent mode (display shows It On)
  - Rotate the right knob to set the run time of the intermittent mode (minutes/seconds)
  - Press it to confirm
  - Rotate the knob to set the pause time of the intermittent mode (minutes/seconds) (display shows It off)
  - Press it to confirm

**Note:** when the intermittent mode is set, the dedicated led (A) is on during the run time and it is blinking during the pause time.

**Note:** select nnSS on the display to disable the function.

The stirring function will be interrupted cyclically.

- Pulse
  - Rotate the right knob to set Pulse ON or OFF
  - Press it to confirm

With this function the instrument will shake at full speed, with immediate ramp, the samples for 2 seconds and then stop for 1 second, repeating the operation until the unit is manually stopped or the timer reaches zero.

**Note:** when the pulse mode is set, the dedicated led (B) is on during the run time. Display will show Puls --- with the number of pulses performed.

**Note:** select OFF to disable the function.

To exit the function menu, push the central button again or do not press any buttons/knobs for 5 seconds.

**Note:** each function can be set individually or in combination with the others.

## Lock function

- While the unit is working, push the central button for 3 seconds to enable the lock function.

If the lock function is enabled, it is not possible to change the work conditions.

The led around the central button blinks if the knobs or the central button are pressed/rotated.

- Push the central button for 3 seconds to disable the function.

## 4. EN - Setting mode

Press both knobs for 3 seconds when stirring is OFF to access the setting mode.

Use the “menu” buttons to browse through the available menus, listed below.

Press the speed control knob to access the selected menu. Then, turn the knob to modify the value of the parameter (when it is possible). Press the knob to confirm.

To exit the setting mode, do not press any buttons/knobs for 10 seconds or press both knobs at the same time.

Text shown		Default value	Range	Menu's name and description
Display 1	Display 2			
<i>SPEd</i>	<i>LINt</i>	2500	100 ÷ 2500 rpm 5 rpm steps	<u>Stirring speed limit</u> It limits the maximum value of speed setpoint for the motor.
<i>Strt</i>	<i>NOdE</i>	<i>StOP</i>	<i>StOP – run</i>	<u>Starting mode</u> ➤ <i>StOP</i> : when the instrument is powered, the stirring function is not active and display shows OFF ➤ <i>run</i> : when the instrument is powered, it starts working at the last speed setpoints
<i>Strt</i>	<i>RAmP</i>	<i>MED</i>	<i>LOW-MED-HIGH</i>	<u>Increasing stirring speed ramp</u> From 0 to 2500, empty: LOW ~20 seconds MED ~10 seconds HIGH ~5 seconds
<i>beep</i>		Yes	Yes - No	<u>Acoustic signal</u> When timer reach 0, the unit emits three beeps.
<i>COun</i>	<i>NOtO</i>		---	<u>Motor operating time</u> It shows motor operating times. Operating times are shown in hours up to 9999 hours (around 416 days). Then they are shown in days (a digits' decimal point is turned on).
<i>rSEt</i>		<i>NO</i>	<i>YES – NO</i>	<u>Reset</u> It allows the reset of all the system settings to the default values.

## 5. EN - Error and warning messages

When the display shows an error message, the instrument's functions are stopped automatically.

<b>AL3</b>	The stirring system doesn't run correctly
------------	---

To remove the error message, restart the instrument.

If the alarm persists on the display, please contact VELP Scientifica's service department.

## 6. EN - Maintenance

No routine or extraordinary maintenance is necessary except periodically cleaning. The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes, and spills should be removed promptly in order to avoid damages.

Support rack and Sample rack for test tube are consumable parts that wear out during instrument operation. These parts are not covered by the warranty

### Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

### Repair

Repairs must be carried out by authorized VELP Scientifica's personnel only.

The instrument must be transported using the original shockproof packaging when freight forwarders, couriers or others are taking care of the unit transport. Follow the indications on the original packaging (e.g. place the instrument on a pallet).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

## 7. EN - Technical data

<b>General features</b>	Model	F202A0460
	Voltage	100-240 V 50/60 Hz
	Power input	40 W
	Current consumption	< 1A
	Dimensions (WxHxD)	245 x 380 x 490 mm (9.6 x 14.9 x 19.3 in)
	Weight	19 kg (42 lb)
	Construction material	Painted steel and stainless steel
	Working in continuous	Max 2 hours
	Support plate dimension	294 x 194 mm (11.6 x 7.6 in)
	Maximum load on the support plate	4,5 kg
	Support system	4 anti-sliding suction cup feet
	Noisiness	<< 80 dBa
	Environmental temperature admitted	+5...+40 °C
	Storage temperature admitted	-10...+60 °C
	Max humidity	80%
	Overvoltage category	II
	Pollution degree CEI EN61010-1	2
	Max altitude	2000 m
	<b>Stir</b>	Speed setting
Speed control		Electronic
Speed range		100 ÷ 2500 rpm (5 rpm steps)
Speed resolution		1 rpm
Motor type		BLDC
Autoreverse / Intermittent mode		5 s ÷ 99 min 59 s (1 s steps)
Type of movement		Orbital
<b>Timer</b>	Orbital diameter	3.6 mm
	Countdown range	1 min ÷ 99 h 59 min (1 min steps)
<b>Counters</b>	Motor	Operating hours

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use  
Part 1: General requirements

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -  
Part 1: General requirements

and satisfies the essential requirements of the following directives:

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2006/42/EC Machinery directive

2015/863/EU (RoHS III) Restriction of the use of certain hazardous substances

2012/19/EU (WEEE) Waste of electric and electronic equipment

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung / 加修改.

## 11. Declaration of conformity <sup>UK</sup> <sub>CA</sub>

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured according to the relevant designated standards:

EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements
EN 61010-2-051	Particular requirements for laboratory equipment for mixing and stirring
EN 62479	Assessment of the compliance of low-power electronic and electrical equipment with the basic restrictions related to human exposure to electromagnetic fields (10 MHz to 300 GHz)

and satisfies the essential requirements of regulations:

S.I. 2008/1597	Supply of Machinery (Safety) Regulations 2008
S.I. 2012/3032	Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012
S.I. 2013/3113	Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

## Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO 14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

### Analytical instruments

Elemental Analyzers  
Digestion Units  
Distillation Units  
Solvent Extractors  
Fiber Analyzers  
Dietary Fiber Analyzers  
Oxidation Stability Reactor  
Consumables

### Laboratory Equipment

Magnetic Stirrers  
Heating Magnetic Stirrers  
Heating Plates  
Overhead stirrers  
Vortex mixers  
Dispersers  
COD Thermoreactors  
BOD and Respirometers  
Cooled Incubators  
Flocculators  
Overhead Shakers  
Turbidimeter  
Open Circulating Baths  
Pumps

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VELP Scientifica Srl