IKA


designed for scientists

RET basic *safety control*

| | | |
|---------------------------------------|----|----|
| Betriebsanleitung Ursprungssprache | DE | 4 |
| Operating instructions | EN | 14 |
| Mode d'emploi | FR | 24 |
| 使用说明 | ZH | 34 |

| | | |
|----------------------------|----|----|
| Indicaciones de seguridad | ES | 44 |
| Veiligheidsaanwijzingen | NL | 46 |
| Awertenze per la sicurezza | IT | 48 |
| Säkerhetsanvisningar | SV | 50 |
| Sikkerhedshenvisninger | DA | 52 |
| Sikkerhedsanvisninger | NO | 54 |
| Turvallisuusohjeet | FI | 56 |
| Instruções de segurança | PT | 58 |
| Wskazówki bezpieczeństwa | PL | 60 |
| Bezpečnostní pokyny | CS | 62 |
| Biztonsági utasítások | HU | 64 |
| Varnostna navodila | SL | 66 |
| Bezpečnostné pokyny | SK | 68 |
| Ohutusjuhised | ET | 70 |
| Drošības norādes | LV | 72 |
| Saugos nurodymai | LT | 74 |

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| EN | |
|----------|---|
| A | Mains switch |
| B | Adjustable safety circuit |
| C | Rotating / pressing knob - Temperature setting |
| D | Rotating / pressing knob - Speed setting |
| E | LED display, heater |
| F | LED display, motor |
| G | LED, heating plate |
| H | LED, external temperature sensor (PT 1000) |
| I | LED (set=set point value) |
| J | LED, external temperature sensor (ETS-D) |
| K | Power socket |
| L | Connection for PT 1000 temperature sensor series, temperature probes or contact plugs |
| M | Threaded hole for stand |

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| | Page |
|--------------------------------------|------|
| EU Declaration of conformity | 14 |
| Safety instructions | 14 |
| Unpacking | 16 |
| Intended use | 16 |
| Assembling the stand | 16 |
| Operation | 17 |
| Setting the operating mode | 18 |
| Er 5 | 18 |
| Setting the safety temperature limit | 19 |
| Setting the temperature control mode | 19 |
| Maintenance and cleaning | 20 |
| Accessories | 20 |
| Error codes | 21 |
| Technical data | 22 |
| Warranty | 23 |

EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61010-2-010, EN 61010-2-051, EN 61326-1, EN 60529 and EN ISO 12100.

A copy of the complete Declaration of Conformity or further declarations of conformity can be requested at sales@ika.com.

Safety instructions

/// General information

• Read the operating instructions in its entirety before using the device and follow the safety instructions.

- Keep the operating instructions in a place where it can be accessed by everyone.
- Ensure that only trained staff work with the device.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Electrical outlet must be grounded (protective ground contact).
- **Attention – Magnetism!**
Effects of the magnetic field have to be taken into account (e.g. data storage media, cardiac pacemakers ...).
- **Risk of burns!**
Exercise caution when touching parts of the housing and the heating plate. The heating plate can reach dangerous temperatures. Pay attention to the residual heat on the heating plate after switching off the stirrer. The device may only be transported when the heating plate has cooled down.

/// Device design

- Do not use the device in explosive atmospheres, it is not EX-protected.
- With substances capable of forming an explosive mixture, appropriate safety measures must be applied, e.g. working under a fume hood.
- To avoid body injury and property damage, observe the relevant safety and accident prevention measures when processing hazardous materials.
- Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- The feet of the device must be clean and undamaged.
- Ensure that the power cord set / temperature sensor cable does not touch the heating plate.
- Check the device and accessories for damage before each use. Do not use damaged components.

/// Permissible medium / contaminants / side reactions

- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Beware of hazards due to:
 - flammable materials,
 - combustible media with a low boiling temperature,
 - glass breakage,

- incorrect container size,
- overfilling of media,
- unsafe condition of container.

- Process pathogenic materials only in closed vessels under a suitable fume hood.

• **Warning!**

The safety temperature must be set in accordance with EN 61010-2-010 Chapter "Requirements for devices containing or using flammable liquids".

- The surface temperature of the flammable medium that is exposed to air may not exceed its flash point.

A danger usually arises if a medium is heated in open vessels.

- The surface temperature of the heating device (e.g. the mounting plate) may not exceed the value of $(t - 25) \text{ }^\circ\text{C}$ (= set value of the safety circuit) on the surface of the flammable medium and in contact with air, whereby t is the fire point of the liquid.

A danger usually arises if a medium is heated in glass vessels (glass breakage).

- If a setting made by the user (medium temperature or safety temperature) could bring a flammable medium into a state in which the conditions mentioned above could be exceeded, additional measures must be introduced that will protect the user from this danger.

- The heating plate can heat up due to the action of the magnets at high motor speeds, even if the heater is switched off.

- Please consider any possible contaminations and unwanted chemical reactions.

- It may be possible for wear debris from rotating accessory parts to reach the material being processed.

- When using PTFE-coated magnetic bars, the following has to be noted: Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above $300 \text{ }^\circ\text{C} - 400 \text{ }^\circ\text{C}$. Only elementary fluorine, chlorotrifluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.

(Source: Römpps Chemie-Lexikon and "Ulmann", Volume 19)

/// Procedures during sample runs

- Wear your personal protective equipment in accordance with the hazard category of the media to be processed. There may be a risk from:

- splashing and evaporation of liquids,
- ejection of parts,
- release of toxic or combustible gases.

- Reduce speed if:

- medium splashes out of vessel because the speed is too high,

- device is not running smoothly,
- container moves on the base plate.

/// Accessories


- Safe operation is guaranteed only with the use of original IKA accessories.
- Ensure that the external temperature sensor is inserted into the medium to a depth of at least 20 mm when connected.
- Always disconnect the plug before attaching accessories.
- Accessories must be securely attached to the device and cannot come off by themselves. The centre of gravity of the assembly must lie within the surface on which it is set up.
- Observe the operating instructions of the accessories.

/// Power supply / Switching off the device

- The device will automatically restart in mode B following any interruption to the power supply.
- The voltage stated on the type plate must correspond to the mains voltage.
- The outlet for the mains plug must be easily accessible.
- The device can only be disconnected from the mains outlet by pulling out the power cord set or the connector plug.

/// For protection of the equipment

- The device may only be opened by qualified and IKA approved experts.
- Do not cover the device, even partially e.g. with metallic plates or film. This may result in overheating.
- Protect the device and accessories from bumps and impacts.
- Keep the base plate clean.
- Observe the minimum distances:
 - between devices min. 100 mm,
 - between device and wall min. 100 mm,
 - above the device min. 800 mm.

Unpacking

- Unpacking
 - Please unpack the device carefully
 - In the case of any damage a detailed report must be set immediately (post, rail or forwarder)
- Scope of delivery
 - Heating magnetic stirrer
 - Mains cable
 - Operating instructions
 - Screwdriver
 - Temperature sensor PT 1000
 - Protection cover

Intended use

- Use
 - The magnetic stirrer is suitable for mixing and / or heating substances.
- Area of use

Indoor environments similar to that a laboratory of research, teaching, trade or industry area.

The safety of the user cannot be guaranteed:

 - if the device is operated with accessories that are not supplied or recommended by the manufacturer,
 - if the device is operated improperly or contrary to the manufacture's specifications,
 - if the device or the printed circuit board are modified by third parties.

Assembling the stand


- ☞ Remove screw plugs (M)
- ☞ Remove the protective cap from the support rod
- ☞ Put the washer between housing and nut
- ☞ Screw the support rod onto the device by hand until the end stop is reached
- ☞ Use an A/f 17 spanner to tighten the M10 nut
- ☞ Accessories should be attached using cross sleeves

Note: When using bath attachments with a diameter over 180 mm, please use support rod in conjunction with an extension. (See "Accessories")

Risk of tipping!

Make sure that the center of mass of the attached device does not protrude over the safe area indicated with a dashed rectangle.

Operation

Ensure that the protective film is removed from the base plate before use!

| | |
|---|--|
| Commissioning | ☞ Move device switch (A) to the OFF position |
| | ☞ Insert the mains power cable into the power socket (K) |
| Stirring | ☞ Move device switch (A) to the ON position |
| | ➤ Operating mode A is preset (see "Setting the operating mode") |
| | ☞ Adjust the motor speed using the rotary knob (D) |
| Heating | ➤ The value selected will be shown on the display (F) |
| | ☞ Set the start point for the agitation function by pressing the rotary knob (D) |
| | ⓘ ➤ The displayed value will flash until the desired motor speed is reached |
| | ☞ Set the safe temperature limit (see "Setting the safe temperature limit") |
| | ☞ Adjust the set-point temperature using the rotary knob (C) |
| | ➤ The value selected will be shown on the display (E) |
| Connecting external thermometers (direct temperature control in the media) | ☞ Set the temperature control mode (see "Setting the temperature control mode") |
| | ☞ Set the start point for the heating function by pressing the rotary knob (C) |
| | ⓘ The set-point and actual temperatures will be shown alternately on the display (E) |
| Operation | ➤ When the heating is switched on, the LED "Heating plate" (G) and the LED "Set-point value" (I) will be lit. |
| | ⓘ ➤ During agitation and standby operation, the display (E) will show $\overline{H.O.}$ if the base plate temperature is above 50 °C. |
| | ☞ Move device switch (A) to the OFF position |
| | ☞ Detach contact plug (L) |
| | ☞ Attach a safety contact thermometer or a PT 1000 temperature sensor to the socket (L) |
| Temperature sensor PT 1000 | ☞ Move device switch (A) to the ON position |
| | ⓘ Temperature sensor PT 1000 ➤ The actual temperature for the temperature sensor shown on display (E) will correspond to the temperature of the media. The LED "external temperature sensor" (H) will be lit |
| Contact thermometer e.g. ETS-D5 | ⓘ Contact thermometer e.g. ETS-D5 ➤ Follow the operating instructions for the contact thermometer |
| | The LED "decimal point for the temperature display" (J) will be lit |
| When a contact thermometer is connected, the display (E) will only show the set-point temperature that has been set | |

Setting the operating mode

Operating the device in mode A, B or D

Mode A

All settings will be stored if the device is switched off or disconnected from the power supply. The agitation and heating functions will be set to OFF when the device is powered on.

Mode B

All settings will be stored if the device is switched off or disconnected from the power supply. The agitation and heating functions will be set to ON or OFF when the device is powered on, depending on the previous status of the device.

Mode D

In Mode D the device behaves the same as in Mode A with the exception that:

- The temperature setting must be confirmed by pressing the temperature knob. To change the temperature setting, turn the temperature knob until the desired value appears.
- The new value flashes for 5 seconds on the display. Confirm the new temperature setting by pressing the temperature knob, otherwise the temperature setting value will skip back to its previous value.

Factory setting: mode A

The mode selected will be shown on the display when the device is started up.

Change the mode

- ☞ Move device switch (A) to the OFF position
 - ☞ Press and hold rotary knobs (C and D)
 - ☞ Move device switch (A) to the ON position
 - ☞ Release rotating knobs (C and D)
- ⇒ The set value is indicated on the display (F)

xx A

Note: Mode D is available from software version 40

Er 5

Providing safety protection, Error 5 is taken to recognise that the sensor has not been placed in the medium when the heating is switched on.

Note: The user can set a value from 1 to 30 min for this time limit depending on the application.


If time limit is set to 0 min, the Error 5 detection is disabled.

This function will only be active if:

- Sensor temperature is < 50 °C
- Difference target temperature / sensor temperature > 5 K

Factory setting: 5 min

The time selected will be shown on the display when the device is started up if the sensor is connected.

Changing the Error 5 settings

- ☞ Set target temperature to 5 °C
- ☞ Switch off device
- ☞ Hold down temp knob
- ☞ Switch on device
- ☞ Select time value from 0...30 min (in 1 min steps) by turning the speed knob
- ☞ Release the temp knob

Note: Error 5 is available from software version 40.

Setting the safe temperature limit

The maximum achievable heating plate temperature is restricted by an adjustable safety temperature limit. Once this limit has been reached, the device stops heating.

The safe temperature limit must always be set at least 25 °C lower than the fire point of the media to be processed!

The adjustable maximum heating plate temperature is at least 10 °C under the set safety temperature limit.

Setting range: [50 °C] to [max. SET temperature + 50 °C]

Factory setting: [max. SET temperature + 50 °C]

Setting the safety temperature limit

After switching on the device, the safety temperature limit can be adjusted using a screwdriver.

Do not turn the setting screw beyond the clockwise or anticlockwise stop. This will cause irreparable damage to the potentiometer.

- Set the main switch to the ON position
- Using the screwdriver supplied, turn the "Safe Temp" setting screw to the clockwise stop.
- Use the temperature rotary knob to set the target temperature to the desired "Safe Temperature" and wait until the temperature is achieved.

- Turn the "Safe Temperature" setting screw slowly anticlockwise until the heating function switches off and the display shows [Er 25].
- The "Safe Temp" value is displayed on the screen.

Setting the safety temperature limit as to be done as described.

The value "Safe Temperature" shown on the display only serves for visualization.

Functional test: safety circuit switch-off

- Heat the device to a temperature of over 50 °C.
- Set the safety temperature limit as far to the left as possible (50 °C) and set the mains switch to the ON position.
- The display will show: [Er 25]

Setting the temperature control mode

When using an external PT1000 temperature sensor, the user can choose between two types of control:

PI mode

Good control results, minimised overshooting, slow rise in temperature.

2P mode (two-point controller)

Maximum heating rate, increased overshooting

Factory setting: PI mode

Changing the temperature control mode

1. Set target temperature to 2 °C
2. Switch off device
3. Hold down temp knob
4. Switch on device

This activates or deactivates the two-point controller (2P) for the external PT 1000 temperature sensor (2P appears on the display).

When switching on the device, 2P is also displayed next to the operating mode A/B when the two-point controller is activated.

When the heating function is active in 2P mode, (with external PT 1000 temperature sensor) the actual temperature flashes - that notifies the user that the process temperature has been greatly overshot!

Maintenance and cleaning

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

Cleaning

- For cleaning disconnect the mains plug!
- Use only cleaning agents which have been approved by IKA to clean the devices:
 - Water containing surfactant / isopropyl alcohol.
- Wear protective gloves during cleaning the devices.
- Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- Do not allow moisture to get into the device when cleaning.
- Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

Ordering spare parts

When ordering spare parts, please give:

- Device type.
- Serial number, see type plate.
- Item number and designation of the spare part, see www.ika.com.

Repair

Please send instrument in for repair only after it has been cleaned and is free from any materials which may constitute a health hazard.

For this you should request the "Safety Declaration (Decontamination Certificate)" from IKA, or use the download printout of it from the IKA website www.ika.com.

Return the instrument in its original packaging. Storage packaging is not sufficient. Also, please use suitable shipping package materials.

Accessories

- Stirring bars:
 - ø 6 mm, length upto 15 mm
 - ø 7 mm, length upto 60 mm
 - ø 10 mm, length upto 80 mm
- RS 1 Set of magnetic stirring bars
- RSE Stirring bar remover
- H 16 V Support rod
- H 16.3 Extension
- H 38 Holding rod
- H 44 Cross sleeve
- ETS-D5 Contact thermometer
- ETS-D6 Contact thermometer

for further accessories see www.ika.com

Error codes

Any malfunctions during operation will be identified by an error message on the display (E and F).

Proceed as follows in such cases:

- ☞ Switch off device using the main switch (A)
- ☞ Carry out corrective measures
- ☞ Restart device

| Error code | Cause | Effect | Solution |
|------------|--|--------------------------|--|
| E3 | Temperature inside device too high | Heating off | - Switch off device and allow to cool down |
| E4 | Motor blockage | Heating off Motor off | - Switch off device - <i>Warning! Only to be carried out by authorized service personnel: Carry out an internal test on the device to check the plug-in connector for the motor</i> |
| E5 | No temperature increase in sensor at continuous temperature difference | Heating off | - Place the sensor in the medium - Reduce the volume of the media - Use a carrier fluid with better heat conductivity properties - Replace the glass vessel with a metal pot - Increase the "Time-out" period |
| E6 | Break in safety circuit | Heating off | - Plug in contact plug (L) - Plug in PT 1000 contact thermometer / temperature sensor - Replace faulty connecting cable, plug, or contact thermometer |
| E24 | Surface temperature (temperature of control sensor): of the base plate is higher than the set safe temperature limit | Heating off | - Switch off device until the surface temperature of the base plate is lower than the selected safe temperature limit - Set a higher safe temperature limit |
| E44 | Surface temperature (temperature of safety sensor): of the base plate is higher than the set safe temperature limit | Heating off | - Switch off device until the surface temperature of the base plate is lower than the selected safe temperature limit - Set a higher safe temperature limit |
| E25 | Heating and switching element monitoring | Heating off | - Switch off device - Set the safe temperature limit > 55 °C see also "Functional check of inactivating the safety circuit" - <i>Warning! Only to be carried out by authorized service personnel: Carry out an internal test on the device to check the plug-in connector for the heating element</i> |
| E26 | Difference between temperature of safety sensor and temperature of control sensor control temperature > (safety temperature + 40 K) | Heating off | - Switch off device - <i>Warning! Only to be carried out by authorized service personnel: Carry out an internal test on the device to check the plug-in connector for the temperature sensor</i> |
| E46 | Difference between temperature of safety sensor and temperature of control sensor safety temperature > (control temperature + 40 K) | Heating off | - Switch off device - <i>Warning! Only to be carried out by authorized service personnel: Carry out an internal test on the device to check the plug-in connector for the temperature sensor</i> |

If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:

- Contact the service department,
- Send the device for repair, including a short description of the fault.

Technical Data

| | | | | | |
|---|--------------|--|--|-----------|---------------------------|
| Number of stirring positions | | 1 | Connection for ext. temperature sensor | | PT 1000 |
| Stirring quantity max. per stirring position (H ₂ O) | | 20 | Temperature setting resolution of medium | K | 1 |
| Motor rating output | W | 9 | Adjustable safety circuit | °C | 50 - 360 |
| Direction of rotation | | right | Set-up plate material | | stainless steel 1.4301 |
| Speed display set-value | | LED | Set-up plate dimensions | mm | Ø 135 |
| Speed display actual-value | | LED | Automatic reverse rotation | | no |
| Speed control | | Turning knob | Intermittent mode | | no |
| Speed range | rpm | 50 - 1700 | Viscosity trend measurement | | no |
| Setting accuracy speed | rpm | 10 | Break detection stirring bar | | no |
| Stirring bar length | mm | 20 - 80 | Timer | | no |
| Self-heating of the hotplate by max. stirring (RT:22°C/duration:1h) | K | 15 | pH measurement | | no |
| Heat output | W | 600 | Graph function | | no |
| Temperature display set-value | | LED | Programs | | no |
| Temperature display actual-value | | LED | Sensor in medium detection (Error 5) | | yes |
| Temperature unit | | °C | Weighing function | | no |
| Heating temperature | °C | (Room temp. + device self heating) - 340 | Dimensions (W x H x D) | mm | 160 x 95 x 270 |
| Heat control | | Turning knob | Weight | kg | 2.5 |
| Temperature setting range | °C | 0 - 340 | Permissible ambient temperature | °C | 5 - 40 |
| Heating rate heating plate | K/min | 7 | Permissible relative humidity | % | 80 |
| Temperature setting resolution of heating plate | K | 1 | Protection class according to DIN EN 60529 | | IP 42 |
| | | | RS 232 interface | | no |

| | | |
|---------------|-----------|-----------|
| USB interface | | no |
| Analog output | | no |
| Voltage | V | 220 - 230 |
| Frequency | Hz | 50/60 |
| Power input | W | 630 |

Subject to technical changes!

Warranty

In accordance with IKA Terms and Conditions of Sale, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

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IKA India Private Limited
IKA England LTD.
IKA Vietnam Company Limited
IKA Works (Thailand) Co. Ltd.