

• BelloCell-500 Technical Report VIII



Cultivation of HEK 293 cells

**Description** BelloCell-500 provides a powerful tool to achieve high cell density and high productivity of target bioproducts in a cell culture because it has a unique feature of offering high oxygen transfer and low shear stress culture environment. Users can easily collect highly concentrated cells, virus or secreted products from one 500 ml BelloCell-500 bottle. In this study, the application of BelloCell-500 for growth of HEK-293 cells is illustrated.  $1.27 \times 10^8$  cells/bottle was seeded and obtain a total of  $3.26 \times 10^9$  cells counted by crystal violet dye nuclei count method at 358 hours, with a total 26 folds increase of cell population. It took 7 days to grow from  $1.27 \times 10^8$  cells to  $2 \times 10^9$  cells. However, it took another 7 days to grow from  $2 \times 10^9$  cells to  $3 \times 10^9$  cells. Glucose concentration in the culture medium was monitored and kept above 1.0 g/L. This technical sheet provides a general protocol for users to start up their culture. However, the optimum condition of each cell culture for each case may require the users to determine.

**Material**

Device	Cell Line	Medium	Seed
BelloCell-500	HEK-293	$\alpha$ MEM/10%FBS + 2.5 g/L glucose + 2.5 mM glutamine + 2.2 g $\text{NaHCO}_3$	$1.27 \times 10^8$ cells/bottle

Protocol \*Detail protocol is in General Instruction Manual

**Inoculum preparation** Prepare one roller bottle. Seed with  $2.5 \times 10^7$  cells total. Culture at 37°C, 5%  $\text{CO}_2$  for total 5 days. Replenish medium at day 3<sup>rd</sup>. Harvest cells by standard trypsinization protocol. Prepare  $1.27 \times 10^8$  suspend cells with viability of 97.42%, and concentrate cells in 50 ml culture medium. **Inoculation** Pre-warm  $\alpha$ MEM/10%FBS medium in 37°C water bath. Take out one BelloCell-500 bottle aseptically and place in a biosafety cabinet. Open the cap and add 450 ml culture medium in each bottle. Dispense 50 ml media containing  $1.27 \times 10^8$  suspend cells on top of the matrix of BelloCell-500.



Bring the bottle and lock up on the BelloStage console in incubator at 37°C, 5% CO<sub>2</sub> and start the run immediately. Avoid swirling or shaking the bottle before start compression.

**Immobilization** Set up operation parameters on the BelloStage control box and start the controller by pressing “START” button. The inoculation parameters are set as below:

Rising rate	Top holding time	Down rate	Bottom holding time
2.0 mm/s	20 sec	2.0 mm/s	0 sec

**Culture** After 3.5 hours, switch the parameters to culture parameters. The culture control parameters are set as below:


Up rate	Top holding time	Down rate	Bottom holding time
1.5 mm/s	0 sec	1.5 mm/s	1 min 30 sec

Monitor the residual glucose concentration and the color of medium in order to predict the time to change culture medium. *The setup parameters are only for reference. It does not necessary to be optimum parameters.*

**Cell Harvest** The cell harvest was followed according to the protocol on the CD manual.


## Results

### Glucose and Lactate profile






Glutamine and Ammonia profile




Glucose uptake rate and Lactate production rate profile






Cell grow curve by crystal violet dye nuclear count method




PH/CO2





### Cell Harvest



The culture of HEK 293 cells in BelloCell-500 is very smooth and the grow rate is slower than the other commonly used cell line, it require 7 days to have nearly 20 folds increase of cell population. The maximum cell density in BelloCell-500 system for HEK293 cells is around  $3.5 \times 10^9$  cells/bottle, and will require 12 days culture to reach this value. For adenovious production, we suggest to seed cells with  $2 \times 10^8$  cells/bottle, and culture for 6~7 days until cell density reach above  $2 \times 10^9$  cells before start infection.

#### Note

HEK293 cell is very sensitive to trypsin and easy to detach. Over trypsinizing the cells will cause the difficulty for cells to be immobilized in the bottle and cause a result of slow growth or even fail to growth. To avoid this, try to minimize cell dissociation process by shortening the trypsinization incubation time (within 3 mins) and terminate the enzymatic reaction by adding serum or trypsin inhibitor. CESCO also develops another bottle to enhance cell immobilization efficiency, i.e. BelloCell-500AP. If users are interested with the product, please contact CESCO Bioengineering ([info@cescobio.com.tw](mailto:info@cescobio.com.tw)) directly or your local distributors.

- BelloCell-500 Technical Report VIII •  
Cultivation of HEK 293 cells



### Summary

Seed	Inoculum volume	Medium volume	Medium
$1.27 \times 10^8$ cells/bottle	50 ml/bottle	500 ml/bottle	$\alpha$ MEM/10%FCS
Total culture age	Total medium consumed	Total medium replenish frequency	Final cell density (nuclear count)
356 hrs	3500 ml	6	$3.25 \times 10^9$ cells/bottle

Please contact Cescobio Bioengineering Technical support for any questions or comments.

http:// [www.cescobio.com.tw](http://www.cescobio.com.tw)

e-mail: [info@cescobio.com.tw](mailto:info@cescobio.com.tw)

