Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
操作指南

OHS 70RS Digital

F201A0723

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 基本信息

Before using the unit, please read the following instruction manual carefully.

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Attention, surface chaude!

Prudencia, superficie caliente!

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2012/19/UE.

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The product can be used with flammable liquids.

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This unit must be used for indoor laboratory applications only.
The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a way not specified by the manufacturer or with unspecified accessories, the product's safety may be compromised.

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IEC/EN 61010-1
IEC/EN 61010-2-051

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improve their quality.
Nell'impegno di migliorare costantemente la qualità dei prodotti, VELP si riserva la facoltà di variarne le caratteristiche.
Dans le but d'améliorer constamment la qualité de ses produits, VELP se réserve le droit d'apporter des modifications aux caractéristiques de ceux-ci.
VELP se reserva el derecho de modificar las características de productos con el fin de mejorar constantemente su calidad.
VELP behält sich zum Zwecke der ständigen Verbesserung der Produktqualität das Recht auf Änderung der Geräteeigenschaften vor.
VELP 保留修改其产品特性的权利，以不断提高其质量。

Safety Regulations / Norme di Sicurezza / Consignes de Sécurité / Advertencias de Seguridad / Sicherheitshinweise / 安全法规

The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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Use only the power cable provided with the instrument.

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The values indicated on the rating plate of the instrument must correspond to those of the a.c. mains supply.

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Fasten the control unit to the support rod (accessory) using the fixing block or placing it into specific drive unit holder (accessory). Position the instrument with a distance from the wall of 30 cm (at least).

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Safe working conditions are ensured only when the accessories described in the dedicated chapter are used.

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The working speed set on the instrument must be such as to avoid wobbling and/or splashes.

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Do not use with explosive or dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie or to stir harmful liquids prior using protective measures according to the safety standards of the processed products and/or in force in the laboratories including personal protective equipment and the presence of an extraction hood which ensures at least 10-fold air change in accordance with the standards EN 14175 and DIN 12924.

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It is dangerous to run the unit with the stirring blade turning in free air. Always place the stirring shaft in the receptacle before turning the unit on.

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It is responsibility of the user to appropriately decontaminate the instrument in the case that dangerous substances fall on or in it, according to the safety datasheet of substances used and to the current laboratories' safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts about the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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The solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the placement in a hood and personal protective equipment (masks, gloves, goggles, etc.).

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Switch off the instrument before removing stirring shaft.

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The OHS 70 RS Digital is designed to deliver powerful and precise mixing for fluid samples across low to high viscosities.

This instrument is used in a wide range of laboratory applications, operated either with stirring shafts or integrated with mini-jacketed reactors, particularly where constant speed and accurate torque control are required.

To maximize flexibility, safety, and durability, the motor is decoupled from the electronic control unit. Thanks to a 2-meter connection cable, the drive unit can be positioned conveniently on the Stand support, or even outside a fume hood, protecting the electronics and ensuring maximum instrument robustness and long-term reliability.


In addition, a broad range of interchangeable stirring shafts allows the OHS 70 RS Digital to adapt seamlessly to various container shapes, volumes, and sample characteristics.

Through the high-contrast LCD display, the instrument allows simple and easy programming of time and stirring speed.

Through the digital display, it is possible to monitor the test at any time, with the analysis time constantly shown.

Using default parameters, to ensure maximum safety, the unit automatically stops at the end of the analysis time.

Integrated safety features, including overload and over-temperature protection, ensure maximum operational safety by automatically stopping the unit if operating limits are exceeded

1 OHS 70RS Digital – Control Unit

2 OHS 70RS Digital – Drive Unit

A Display

B Timer button

C Control knob

D Fixing block

E Chuck

F USB-C port interface

G RS-232 port interface

H Motor connector

I Air vent

J Power connector

K Main switch

2. Installation

Unpacking

- Check the integrity of the unit after unpacking.

The box includes

- OHS 70RS Digital – Control Unit (1);
 - OHS 70RS Digital – Drive Unit (2);
 - Allen key and Allen head screw
 - Support shaft
 - Instruction manual;
 - EU (Schuko) power cable;
 - US (NEMA 5-15P) power cable.
- Note:** stirring shafts are not included and must be purchased separately

Installation

- Place the Control Unit (1) and the Drive Unit (2) on a non-flammable surface;
- Insert the Support shaft into the Drive Unit (2) and secure it tightly using the supplied Allen screw and key;
- Firmly secure the Drive Unit (2) to a suitable support (e.g. A00000484 Double Rod Mini H-Stand Set) using the Support shaft;
- Insert the Control Unit (1) into the dedicated optional accessory (A00000602 - Drive unit holder for OHS 70RS Digital) and place it on a stable, non-flammable surface, or secure it firmly to a suitable support rod using the Fixing block (D);
- Connect the Drive Unit (2) cable to the Motor Connector (H) of the Control Unit (1);
Note: Lock the connection in place by tightening the male connector's threaded ring onto the female connector.
- Position a suitable receptacle under the Drive Unit (2) and secure it;
- Slide the stirring shaft into the chuck (E) and tight it
- Adjust the Drive Unit (2) height to ensure the stirring blade is correctly positioned inside the sample
- Make sure that the technical specification of the unit and of the power supply line are the same;
- Make sure that the electrical network is grounded;
- Make sure that the Main switch (K) is OFF;
- Connect the instrument to an easily accessible socket (compliant with the current safety norms), using only the provided power cable;
Note: Keep the power cable far away from the moving elements.

3. Operation

Power-up

- Turn the unit on using the Main switch (K);
- The display shows the software version and the last settings.
Note: if the set starting mode is STOP ([see chap. 5](#)), the stirring and timer functions are not active. Otherwise, the instrument starts working at the last set speed setpoint with the last saved timer value.

Stirring

- Rotate the Control knob (C) to set the stirring speed;
- Push the Control knob (C) to start stirring;
- The real speed increases until set point achievement;
- A microprocessor ensures constant speed even if the viscosity changes (SpeedServo);
- Rotate the Control knob (C) to adjust the stirring speed;
- Push the Control knob (C) to stop stirring.

Timer

Count-Up:

By default, at first power-up and after resetting the unit ([see chap. 5](#)), the timer is set to Count-Up mode with an initial value of 00:00:00 (hours : minutes : seconds of operation).

When the timer operates in Count-Up mode, the Count-Up icon  is displayed.

Count-Up operation:

- Press the Timer button (B) to start counting, regardless of whether stirring is active or inactive;
- Once Count-Up has started, press the Timer button (B) once to pause it;
- To reset the Count-Up, double-press the Timer button (B) quickly;
- Upon reaching 99:59:59, the Count-Up timer automatically restarts from 00:00:00;
- Count-Up also automatically starts and stops when stirring is activated or stopped by pressing the Control knob (C).

Count-Down:

When the timer is stopped, hold the Timer button (B) for 2 seconds to switch from Count-Up to Count-Down mode. The

Count-Down icon  will appear on the display, the display shows the previously saved Count-Down value: The default Count-Down value after resetting the unit ([see chap. 5](#)) is 1 minute (00:01:00).

Edit and start/stop Count-Down:

- Press and hold the Timer button (H) for 2 seconds. The display digits (hh:mm:ss) will start flashing, indicating that the time can now be adjusted.
- Rotate the Control knob (C) to adjust the duration of the Count-Down up to 99:59:59.
- Confirm and save the new timer value by either pressing the Control knob (C) or the Timer button (B)
- The digits will stop flashing to confirm that the setpoint has been saved.
Note: to change the setpoint again, simply repeat the procedure described above.
- Press the Timer button (B) to start counting, regardless of whether stirring is active or inactive;
- Once Count-Down has started, press the Timer button (B) once to pause it;
- To reset the Count-Down, double-press the Timer button (B) quickly;
- Count-Down also automatically starts and stops when stirring is activated or stopped by pressing the Control knob (C).

When Count-Down expires:

- If the “TIME END” parameter is set to “OFF” ([see chap. 5](#)):
 - Stirring stops and display shows “OFF” and “END”;
 - The instrument beeps 3 times.
- If the “TIME END” parameter is set to “ON” ([see chap. 5](#)):
 - Stirring continues;
 - Display shows “END” without any audible beep.

Clearing the “END” message:

- Press the Timer button (B) to mute the beeps (if active), clear the “END” message, and restore the Count-Down setpoint.
- If stirring was stopped (TIME END = OFF), press the Control knob (C) to mute the beeps and immediately start a new stirring and Count-Down cycle.

Switching from Count-Down to Count-Up:

- While the Count-Down timer is stopped (paused or not started), press and hold the Timer button (B) for 2 seconds until the digits begin flashing.
- Press and hold the Timer button (B) for another 2 seconds.
- Display **B** will show 00:00:00.
- the Count-Up icon  is displayed
- The previously saved Count-Down setpoint and remaining time are preserved in memory for future Count-Down cycles.

4. Connection to PC

The instrument can be connected to a PC using a suitable cable for efficient data logging and remote control via dedicated software.

Hardware Interfaces:

- USB-C port interface (F): standard USB-C interface for PC connection;
- RS-232 port interface (G): 9 pin DE-9 female connector.

Communication Protocols ([see chap. 5](#)):

- “VELP”: select this protocol when using VELP ControllerSoft™;
- “OTHE” (Other): select this protocol when communicating with external third-party software or custom devices (using the commands described below).

“Other” protocol:

Command	Value (range)	Type*	Function
SS	0 (OFF)	Q/N	Sets or reads motor speed (RPM). • Sending SS starts stirring at the specified speed and locks speed adjustment using the Control knob (C) (start/stop via knob remains active); • Send RM to restore manual speed adjustment via the Control knob (C).
	30 - 800		
RM	None	N	Releases PC speed control lock and restores manual adjustment via the Control knob (C). Sending a new SS command locks knob adjustment again.
PI	None	Q	Returns instrument info (software version and serial number). <i>Example: OHS 15 Digital v1.0.0 SN_00609379</i>
MS	Nessuno	Q	Returns current motor status: • MS1: Motor stopped • MS4: Motor running at target speed • MS6: Alarm status
SN	None	Q	Returns instrument serial number

* Command Types: **Q** = Query / Read, **N** = Action / Write

Protocol Syntax

Commands are case sensitive and must be in upper case.

The end-of-line character is a carriage return <CR>, or Hex 0x0D. Do not add a linefeed to the end of line.

Query Format: <COMMAND><CR>	Command Format: < COMMAND ><Value><CR>
This example queries the instrument serial number. Example: • Command = SN<CR> • Reply = SN00123456<CR>	This example sets the motor rotational speed. Example: • Command = SS350<CR> • Reply = SS350<CR>

Communication parameters:

Baudrate [bps]	9600
Data bits	8
Stop bit	1
Parity	None

5. Setting mode

Press and hold the Control knob (C) for 3 seconds when stirring is not active to access the instrument menu.

Rotate the Control knob (C) to browse through the available parameters, listed below.

Press and turn the knob to modify the value of the parameter (except for motor counter). Press the knob to confirm.

To exit the setting mode, press and hold the Control knob (C) for 3 seconds or do not press/turn the knob for 10 seconds.

Instrument menu:

Parameter shown		Default value	Range	Description
LIMIT		800	OFF / 30 – 800	<u>Speed limit</u> It limits the maximum value of speed set point for the motor. If OFF is selected, stirring is “disabled”.
STRT	MODE	STOP	STOP - RUN	<u>Starting mode</u> - STOP: when the instrument is powered, stirring is not active. - RUN: when the instrument is powered, it starts working at the last working conditions.
TIME	END	OFF	ON - OFF	<u>Continue stirring after timer end</u> ➤ If ON, stirring continues after Count-Down ends.
COMM	PROTOD	VELP	VELP - OTHE	<u>Communication Protocol</u> - VELP: for Velp ControllerSoft™; - OTHE: for third-party software or devices.
COUN		-	-	<u>Motor operating time</u> It shows motor operating time. Operating time is shown in hours up to 9999 hours (around 416 days). Then they is shown in days.
CAL	MOTOR	NO	YES - NO	<u>Reset</u> It allows to perform a calibration of the motor (necessary operation if the motor has been replaced).
RSET		NO	YES - NO	<u>Reset</u> It allows to restore all system settings to their default values.

6. Error and warning messages

When the display shows an error message, the instrument functions automatically stop.

AL1	The motor doesn't start or not reach the set speed
AL3	Motor overload (torque out of range) / overcurrent
AL4	High internal temperature
AL5	Electronic board malfunction
AL8	Speed variation / oscillations too high
AL9... - AL15	Hardware malfunction

To remove the error message, restart the instrument from the main switch.

If the alarm persists on the display, please contact VELP Scientifica's technical assistance service.

7. Maintenance

No routine or extraordinary maintenance is required except periodic cleaning.

Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

Repair

Repairs must be carried out by authorized Velp personnel only.

The transport of the instrument by freight forwarders, couriers or others must be carried out using the original shockproof packaging. Follow any instructions on the original packaging (e.g., palletizing).

It is the responsibility of the user to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

8. Technical data

		OHS 70RS Digital
General features	Model	F201A0723
	Voltage and Frequency	115/230 V - 50/60 Hz
	Power input	90 W
	Control Unit - Dimensions (WxHxD)	93 x 95 x 130 mm (3,66 x 3,74 x 5,12 in)
	Control Unit - Weight	0,85 kg (1,87 lb)
	Drive Unit - Dimensions (WxHxD)	67 x 187 x 140 mm (2,64 x 7,37 x 5,5 in)
	Drive Unit - Weight	1,9 kg (4,18 lb)
	Construction material	Stainless steel enclosure – Tempered glass control panel
	User interface	LCD display and touch button
	Working in continuous	Admitted
	Noisiness	<< 80 dBa
	Environmental temperature admitted	+5...+40 °C
	Storage temperature admitted	-10...+60 °C
	Max humidity	80%
	Overvoltage category	II
	Pollution degree CEI EN61010-1	2
	Max altitude	2000 m
Stirring	Stirring capacity	Up to 40 l H ₂ O – Up to 50.000 mPa*s
	Max torque admitted	70 Ncm
	Programmable speed range	0 / 30 – 800 rpm
	Speed setting resolution	1 rpm
	Motor type	BLDC
	Motor rating output	60 W
	Chuck range diameter	2 – 13 mm
	Shaft diameter Pass-Through	N/A
Timer	Count-Up and Count-Down range	1 s ÷ 99 h 59 min 59 s (1 s steps)
Connectivity	USB interface	USB-C
	RS 232 interface	Yes
Counters	Motor counter	Operating hours/days

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We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements
EN 61010-2-051	Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-051: Particular requirements for laboratory equipment for mixing and stirring
EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements

and satisfies the essential requirements of the following directives:

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2006/42/EC	Machinery directive
2014/30/EU	Electromagnetic compatibility directive
2015/863/EU (RoHS III)	Restriction of the use of certain hazardous substances
2012/19/EU (WEEE)	Waste of electric and electronic equipment

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung /
加修改.

11. Declaration of conformity

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2017/1206 Radio Equipment Regulations 2017

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen Velp!

Established in 1983, Velp is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. Velp works according to **ISO 9001**, **ISO 14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortexers and Shakers
Dispersers
Dry Block Heaters and COD
BOD and Respirometers
Cooled Incubators
Flocculators
Turbidimeter
Open Circulating Baths
Pumps

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Velp Scientifica Srl