

## PROCESS FILTER REACTOR, CYLINDRICAL, JACKETED, ELECTRIC OR AIR MOTOR



Our jacketed filter reactors can be used for single or multi-step reactions and filtrations in the same vessel, production and process scale peptide synthesis. Filter base eliminates material transfers and significantly reduces losses. Complete reactor system on our **NEW** open air, auto-center stand. Supplied with the components listed and numbered on the schematic drawing inside (*see inside panels A and B*).

### 10L COMPLETE JACKETED FILTER REACTORS:

Part Number	Description	Motor Type	Overall Clearance		
			Height With Motor(in)	Below Valve(in)	Price
CG-1959-F610	10L Process Filter Reactor, Complete	Electric	74"	19.4"	13,150.00*
CG-1958-F611	10L Process Filter Reactor, Complete	Air	70.2"	19.4"	11,990.00*

### 20L COMPLETE JACKETED FILTER REACTORS:

Part Number	Description	Motor Type	Overall Clearance		
			Height With Motor(in)	Below Valve(in)	Price
CG-1959-F620	20L Process Filter Reactor, Complete	Electric	73.8"	19.4"	13,380.00*
CG-1958-F621	20L Process Filter Reactor, Complete	Air	70"	19.4"	12,435.00*

\* Not Subject To Discount



- Contact us for more information -

3800 North Mill Road, Vineland, NJ 08360

US & Canada Phone: 1-800-843-1794

Fax: 1-800-922-4361

Customer Service: customer-service@cglifesciences.com



# huber

- Sales and Installation
- Service and Support

WWW.CGLIFESCIENCES.COM

# FILTER REACTOR SYSTEMS

Jacketed and Unjacketed

Standard and custom filter reactor systems

Air or electric motors

Filter support system can be customized for your filter media



## PROCESS FILTER REACTOR, CYLINDRICAL, UNJACKETED, ELECTRIC OR AIR MOTOR

Our unjacketed filter reactors can be used for single or multi-step reactions and filtrations in the same vessel, production and process scale peptide synthesis. Filter base eliminates material transfers and significantly reduces losses. Complete reactor system on our **NEW** open air, auto-center stand. Supplied with the components listed and numbered on the schematic drawing inside (*see inside panels C and D*).

### 10L COMPLETE UNJACKETED FILTER REACTORS:

Part Number	Description	Motor Type	Overall Clearance		
			Height With Motor(in)	Below Valve(in)	Price
CG-1959-F810	10L Process Filter Reactor, Complete	Electric	74"	19.4"	11,635.00*
CG-1958-F811	10L Process Filter Reactor, Complete	Air	70.2"	19.4"	10,520.00*

### 20L COMPLETE UNJACKETED FILTER REACTORS:

Part Number	Description	Motor Type	Overall Clearance		
			Height With Motor(in)	Below Valve(in)	Price
CG-1959-F820	20L Process Filter Reactor, Complete	Electric	73.8"	19.4"	12,240.00*
CG-1958-F821	20L Process Filter Reactor, Complete	Air	70"	19.4"	11,100.00*

\* Not Subject To Discount



