

When Precision and Uniformity are Key

Optional Digital Thermoregulator with Timer

- Accurate to ± 0.5 °C
- Timer adjustable from 0 to 24:59 (h:min)



HEATER TIMER

Precise temperature regulation and uniform heat distribution in a safe, reliable stirrer.

Stirrer Features

- 135mm diameter aluminum alloy plate for even heat distribution
- Protective coating for easy cleaning and chemical resistance
- Elevated front panel and integral spill channel
- Two-year warranty included

Best Value

Save when you buy the kit!

Part Number	Description	Price
CG-1995-V-50	Magnetic Hot Plate Stirrer Kit with Digital Contact Thermoregulator and Vertical Support Rod	738.00*

Kit Components

CG-1995-V-10	Magnetic Hotplate, Analog	406.85*
CG-1995-V-20	Digital Contact Thermoregulator	314.15*
CG-1995-V-22	Vertical Support Rod	36.05*

If purchased separately **757.05***

Specifications

Magnetic Stirrer

H ₂ O Stirring Quantity (liters)20
Speed Display Analog
Speed Range (RPM) Up to 1200

Heating Function

Timer Function with CG-1995-V-20 0 to 24:59 (h:min)
Heating Power (watts) 630
Temperature Range (°C) RT to 370
Heat Control Accuracy with CG-1995-V-20 ± 0.5
Overtemp Protection Yes
Use With Temperature Sensor CG-1995-V-20
CG-1995-V-20 Display Resolution 0.2

General Data

Dimensions W x H x D (mm) 165 x 115 x 280
Weight (kg) 2.6
Power Requirements 115V, 50/60 Hz

Contact Us About



Where chemistry meets biology.

Order by phone, fax or web:

T 1-800-843-1794 • F 1-800-922-4361
 customer-service@cglifesciences.com
 www.cglifesciences.com

*Not subject to discount.

FN5M0314CG-13065