

CHEM-SPIN SCRAPER ROTARY TOOLS

New Sizes!
Seven different sizes now available.


Laser Etched
Part number on each tool.

CG-1981-C2000

304 Stainless Steel Chain
For rapid sample removal.



Chem-spin in Action - see the video on line!



Before ▶▶▶▶▶ During ▶▶▶▶▶ After



Glass Bushing Adapter included with Complete Kits

See **Complete Kit** details on the back ▶

Chem-spin Scraper Rotary Tools for Round Bottom Flasks

These tools provide an easy method of removing the majority of desired, dried-on solids from the interior walls of round bottom flasks. Motor driven, the spinning chain rapidly and effectively brushes the surface as the loose solids collect on the bottom. Minimal physical effort required. Tedious, time-consuming, manual scraping with a spatula is eliminated.

Chem-spin Rotary Scraper Tools

Description	Part No.	Price
For 100mL Round Bottom Flasks ◀ NEW	CG-1981-C100	48.00*
For 250mL Round Bottom Flasks ◀ NEW	CG-1981-C250	48.00*
For 500mL Round Bottom Flasks	CG-1981-C500	48.00*
For 1000mL Round Bottom Flasks	CG-1981-C1000	49.00*
For 2000mL Round Bottom Flasks	CG-1981-C2000	51.00*
For 3000mL Round Bottom Flasks ◀ NEW	CG-1981-C3000	56.00*
For 5000mL Round Bottom Flasks ◀ NEW	CG-1981-C5000	56.00*

Rotary Tool Motors with 12v Lithium-ion Batteries

Motor with 2 Batteries	CG-1981-C-20	135.00*
Motor with 1 Battery	CG-1981-C-21	89.15*
Spare Battery Only	CG-1981-C-22	50.20*

Borosilicate Glass ⌘ Bushing Adapters with ⌘ 14/35 Inner Bores

⌘ 24/40 Inner Joint	CG-1016-05	32.90
⌘ 29/42 Inner Joint	CG-1016-08	38.95



Patent pending design.



1-800-843-1794

www.cglifesciences.com

customer-service@cglifesciences.com

* Not subject to discount.

Scan to watch
Chem-spin video ▶



CHEM-SPIN SCRAPER ROTARY TOOLS

Kits provide everything you need ▶

Chem-spin **MOTOR KITS** with Rotary Tools, Bushing Adapters, Motors and Two 12V Batteries

Rotary Tools Included	Bushing Adapter Included	Complete Kit Part No.	Complete Kit Price
100, 250 and 500mL	₹ 24/40 with ₹ 14/35 Inner Bore	CG-1981-C-28	311.10*
	₹ 29/42 with ₹ 14/35 Inner Bore	CG-1981-C-29	317.00*
500, 1000 and 2000mL	₹ 24/40 with ₹ 14/35 Inner Bore	CG-1981-C-30	315.10*
	₹ 29/42 with ₹ 14/35 Inner Bore	CG-1981-C-32	321.00*

Chem-spin **TOOL KITS** with Rotary Tools Only

Rotary Tools Included	Description	Complete Kit Part No.	Complete Kit Price
100, 250 and 500mL	Three-piece Tool Kit	CG-1981-C-24	144.00*
500, 1000 and 2000mL	Three-piece Tool Kit	CG-1981-C-25	148.00*
100, 250, 500, 1000 and 2000mL	Five-piece Tool Kit	CG-1981-C-26	244.00*

**New Kits
Make Ordering
Easy!**



Chem-spin **MOTOR KITS**

Include Rotary Tools, Adapters, Motors and Batteries



Bushing Adapters

Available for ₹ 24/40 and ₹ 29/42 Round Bottom Flasks



Chem-spin **TOOL KITS** and **MOTOR KITS**

Multiple Chem-spin rotary tools to fit popular flask sizes.

All components available separately. See other side ▶

Chem-spin **TOOL KITS** include a variety of Chem-spin Rotary Tools.

Chem-spin **Motor Kits** add Motors, Batteries and Bushing Adapters.

Close-up view of PTFE tip that aids in solids removal.



1-800-843-1794
www.cglifesciences.com
customer-service@cglifesciences.com

Scan to watch Chem-spin video ▶



* Not subject to discount. BC1M0916CG-16055r3