




CHEMRxnHUB
COMPACT



ChemRxnHub Compact System, 50mL, Jacketed Lab Reactor

ChemRxnHub Compact process systems are complete with benchtop stand, overhead stirrer and glassware. Available with the optional manifold connection and heat transfer fluid draining system. The unique design of the ChemRxnHub Compact support stand and clamping system allows for quick and easy installation/removal of glassware. Compact design fits into most benchtop hoods.

Part Number	Description	Price
CG-1949-M-0050	50mL Mini Reaction Vessel, Jacketed, Complete with CAT motor	*
CG-1949-M-0050M	50mL Mini Reaction Vessel, Jacketed, Complete with CAT motor, Manifold system	*
CG-1949-M-0050V	50mL Mini Reaction Vessel, Jacketed, Complete with Velp Mini Overhead motor	*
CG-1949-M-0050VM	50mL Mini Reaction Vessel, Jacketed, Complete with Velp Mini Overhead motor, Manifold system	*

Item No.	Qty	Part No.	Description
1	1	CG-2077-U-006	Chem-Stir 6mm, Bearing/Guide 45mm, 19/22
2	1	CG-2079-M-050	Stirrer Shaft, Ground, 6mm, Glass Agitator, Four-Blade
3	1	CG-2044-C-10	Compact Tru-Stir™ Stirrer Shaft Coupling 6mm
4a	1	CG-2037-M-060	Overhead motor CAT
4b	1	CG-2039-V-10	Velp Overhead motor with separate controller
5	1	CG-1929-M-050	50mL Mini Reaction Vessel, Jacketed
6	1	CG-1943-M-60	Lid, Reaction Vessel, 5-Neck 19/22 CN, 14/20 SN's, 60mm
7	1	CG-141-T-11	60mm Quick Release Clamp, PTFE Coated
8	1	CG-147-80	O-ring, Perfluoro (FFKM), 60mm Flange WHITE
9	1	CG-1213-M-19	Condenser, Reflux, 205mm, 19/22 Joint, Removable Hose Connections, 130mm Jacket Length
10	1	CG-1014-19	Adapter, Inlet, 19/22 Inner Joint, 90°
11	1	CG-3002-L-14	14/20 Stopper, Solid
12	1	CG-3002-L-19	19/22 Stopper, Solid
13	1	CG-1044-M-19	19/22 offset thermocouple holder 1/8" Thermocouple Adapter
14	1	CG-1949-M-10	Reactor Support Frame, ChemRxnHub Compact
15	1	CG-1947-B061	60mm Stainless Steel Reactor Clamp, ChemRxnHub Compact
16	2	CG-1969-V-10	90° Adapter, Stainless Steel, 30mm Tool-Free™ O-ring Flange to M16 Male Thread
17	1	CG-2040-C-20	Heavy Duty Clamp Holder, Fits 3/4" Support Rod, Silver Epoxy Coated (Supplied with Velp motor)
18	1	CG-3498-03	1/2" OD Stainless Steel Support Rod, Bent 90° (Supplied with Velp motor)
Manifold Option			
	1	CG-1969-X-10	ChemRxnHub Manifold System to Include: 1ea CG-1969-X-12 (Outlet Manifold Assembly Complete) and 1ea CG-1969-X-14 (Inlet Manifold Assembly Complete)

*Not Subject To Discount

To order, call 1.800.843.1794 - email: cs@cglifesciences.com