


# REACTOR DISTILLATION KITS

## AVAILABLE CUSTOM OPTIONS FOR 30/50L REACTORS

- Jacketed RBF Receivers
- Cylindrical Style Receivers
- Dual Receiver Shelf Kits
- Different Style Condensers
- Detachable Drains for Receiving Vessels



# REACTOR DISTILLATION KITS

## CG-1968-D DISTILLATION KIT FOR 30/50L PROCESS REACTOR SYSTEM, SIDE SHELF


Distillation kit for Chemglass 30 and 50L reactor systems. Supplied complete with support shelf to fit CG-1968-X style frame and the components highlighted in the photo below. This add-on kit makes it easy to convert any of our standard reactors over to a complete distillation system. **PLEASE NOTE: REACTOR SYSTEM MUST BE PURCHASED SEPARATELY.**

### SIDE SHELF DISTILLATION KIT ADD-ON FOR EXISTING REACTOR SYSTEM:

Part Number	Description	Receiving Flask Capacity
CG-1968-D-410	Distillation Shelf and Glassware Kit	10L
CG-1968-D-420	Distillation Shelf and Glassware Kit	20L

### Standard Components Supplied With Complete Distillation Kits:

No.	Qty	Part No.	Description
1	1	CG-1235-A-45	Adapter, Distillation 45/50
2	1	CG-1215-C-10	Condenser, High Efficiency, 45/50
3	1	CG-1014-C-02	Adapter, inlet, 29/42
4	1	CG-9253-10	Kwik Clamp, Small
5	1	CG-9220-E-01	Chain Clamp, 12" Arm
6	1	CG-9253-20	Kwik Clamp, Large
7	1	CG-1968-D-128	Support Rod, Distillation for CG-1968-X frame
8	1	CG-1975-22	3/4 to 1/2 compression fitting
9	2	CG-1975-21	1/2 to 1/2 compression fitting
10	1	CG-9205-E-10	Clamp, 3-finger, 12"
11a	1	CG-1965-D-225	10L receiving flask, 2 Neck, 45/50 (for 10L Distillation system)
11b	1	CG-1965-D-230	20L receiving flask, 2 Neck, 45/50 (for 20L Distillation system)
12	1	CG-1965-D-210	Adapter, Distillation, 45/50 to 0-14 Valve
13	1	CG-1965-D215	Adapter, Distillation, Vent, 45/50 to 0-14 Valve
14	1	CG-1963-D-50	Shelf, Distillation
15	1	CG-1975-03	FEP Tubing; 1/2" OD x 3/8" ID x 10ft Long
--	1	CG-1975-30	Tubing Cutter



Highlighted Items are Included in Kit