

CLS-3788-R42

PREMIER BUILT-IN UNDERCOUNTER REFRIGERATOR



The Premier Built In Undercounter Refrigerators features a microprocessor temperature controller and display, which allows for precise temperature control, verification and recovery versus mechanical style thermostats. An interior fan driven, forced-draft air circulation system provides excellent temperature uniformity and recovery after door openings.

Critical samples and supplies are security backed with an extensive package of audible and visual alarms.



Model	Cu. Ft.	Defrost	Doors	Shelves	W"	D"	H"	RL
CLS-3788-R42	4.6	Cycle	1 Solid	2 Adj./1 Fixed	23 3/4	24 1/2	33 3/8	R600a

- Microprocessor temperature controller
- Digital temperature display
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- 2 probes; 1 in air and 1 in glycol bottle
- Adjustable operating temperature range:
Refrigerator: 1°C – 10°C
- Keyed door lock
- Recessed door handle
- Magnetic door gasket for positive seal

- Leveling legs
- Probe access port (3/8")
- Powder coat exterior
- High density urethane foam cabinet and door insulation
- HFC-free refrigerant (R600a), foam insulation and packaging
- Two year parts and labor warranty
- Defrost
Refrigerator: Cycle defrost
- UL/C-UL listed