

Chemglass Adapter System for use with Fritted, Disposable Filter Funnels

Typically used for solid-liquid separations in chemistry and biological laboratories, set-ups consist of a funnel, vacuum adapter and optional threaded connector for vials.


The vacuum adapters (CG-1049-J-series) are designed so that vacuum can be controlled with a PTFE stopcock to prevent filtrates from being pulled into the vacuum line. Stems on the filter funnels are 80mm long and prevent contamination of the vacuum adapter during the filtration process.

Chemglass Adapter System for use with Fritted, Disposable Filter Funnels

Purchase components separately.

Part No.	Funnel Adapters with Vacuum Stopcocks	Price
CG-1049-J-14	Outer to 14/20 Inner	\$115.95 *
CG-1049-J-24	Outer to 24/25 Inner	\$119.50 *
CG-1049-J-29	Outer to 29/26 Inner	\$119.50 *
CG-1049-J-30	Vial, Outer to 20-400 GPI Thread	\$116.15 *

Part No.	Accessory Threaded Connectors	Price
CG-1318-23	20-400 to 20-400 GPI Thread	\$24.65 *
CG-1318-24	20-400 to 22-400 GPI Thread	\$25.25 *




Disposable Filter Funnels with Polyethylene Frits

Available in Bulk Bags or Small-Quantity Trays

Bulk Packed in ZipLock Bags

Part No.	Funnel Cap. (mL)	Frit Porosity (µm)	Top Style	Case Qty.	Price
OP-6602-10	18	10	Flared	100	\$298.75 *
OP-6602-12	40-60	10	Flared	100	\$328.75 *
OP-6602-13	40-60	40	Flared	100	\$331.25 *
OP-6602-14	110	10	Straight	50	\$251.25 *
OP-6602-16	300	20	Straight	30	\$282.50 *
OP-6602-18	500	20	Straight	20	\$263.75 *



Shelf Pack Convenience Trays

Packaged in a smaller, convenient, shrink-wrapped, SBS-coated, white tray complete with SBS partitions inserted between the top and bottom layers.

Part No.	Funnel Cap. (mL)	Frit Porosity (µm)	Top Style	Case Qty.	Price
OP-6602-C10	18	10	Flared	24	\$88.00 *
OP-6602-C12	40	10	Flared	24	\$97.75 *
OP-6602-C14	110	10	Straight	12	\$76.95 *



US Patent Pending

*Not Subject to Discount

Disposable Filter Funnels & Adapter System

Disposable Filter Funnels with Polyethylene and Diatomaceous Earth Frits

Available in Bulk Bags or Small-Quantity Trays

Bulk Packed in ZipLock Bags

Part No.	Funnel Cap. (mL)	DE Content	Top Style	Case Qty.	Price
OP-6603-10	18	500mg	Flared	50	\$252.50 *
OP-6603-12	40-60	1,500mg	Flared	50	\$268.75 *
OP-6603-13	150	5g	Straight	30	\$237.50 *

Shelf Pack Convenience Trays

Packaged in a smaller, convenient, shrink-wrapped, SBS-coated, white tray complete with SBS partitions inserted between the top and bottom layers.

Part No.	Funnel Cap. (mL)	DE Content	Top Style	Case Qty.	Price
OP-6603-C10	18	500mg	Flared	24	\$151.95 *
OP-6603-C12	40	500mg	Flared	24	\$162.75 *
OP-6603-C13	150	500mg	Straight	12	\$119.30 *

Cone-Shaped, Disposable Filter Funnels, Molded Polypropylene with Polyethylene Frits



Disposable Filter Funnels, Cone-Shaped, PP with PE Frits

Part No.	Funnel Cap. (mL)	Frit Porosity (µm)	Case Qty.	Price
OP-6605-08	8	10	100	\$250.00 *
OP-6605-10	20	10	100	\$303.75 *
OP-6605-12	60	10	100	\$332.50 *

Barrel Shaped, Disposable Filter Funnel, Polyethylene Frit, 2L Funnel

Disposable Filter Funnels, Barrel Shaped. Adapter, Vacuum

Part No.	Description	Price
OP-6602-J-20	Filter Funnel, Disposable, 2L, Polyethylene Frit, Barrel Shaped, Funnel Only	\$295.00 *
OP-6602-J-21	Filter Funnel, Disposable, 2L, Polyethylene Frit, Barrel Shaped, Supplied with Lid	\$310.00 *
CG-1049-M-24	Adapter, Vacuum, 24/25 Outer to 24/25 Inner	\$65.75 *
CG-1049-M-45	Adapter, Vacuum, 24/25 Outer to 45/50 Inner	\$87.55 *



US Patent Pending

*Not Subject to Discount