



Replacement GLASSWARE

Compatible with Mettler EasyMax™
Automated Reactor Systems
2025



One-Piece Reaction Vessels

One-Piece Reaction Vessel, 5-Neck, 50mL, Mettler-Style

One-piece vessel with 19/22 outer CN, (1) 10/18 outer SN, (1) 14/20 outer SN, and (2) 19/22 outer SN's. Vessel is graduated from 10 to 50mL.

Equivalent to Mettler EasyMax® 102 Part # 51161792.

Part No.	Mettler Equivalent	Cap. (mL)	No. of Necks	Body OD (mm)	2024 Price
CG-MR-50V1P	51161792	50	5	43	\$747.45*

One-Piece Reaction Vessel, 5-Neck, 100mL, Mettler-Style

One-piece vessel with 19/22 outer CN, (1) 10/18 outer SN, (1) 14/20 outer SN, and (2) 19/22 outer SN's. Vessel is graduated from 20 to 100mL. Body outside diameter is 57mm.

Equivalent to Mettler EasyMax 102 / 402 Part # 51161786.

Part No.	Mettler Equivalent	Cap. (mL)	No. of Necks	Body OD (mm)	2024 Price
CG-MR-100V1P	51161786	100	5	57	\$624.95*



CG-MR-50V1P ►



◄ CG-MR-100V1P

*Not subject to discount.

CG-MR-50V1 ▶



◀ CG-MR-50V2

CG-MR-100V1 ▶
CG-MR-400V1 ▶

◀ CG-MR-250V5



◀ CG-MR-250L5

Reaction Vessel, Flat Bottom, Flat Flange, 50mL, Mettler-Style

Part No.	Mettler Equivalent	Cap. (mL)	Flange OD (mm)	OAH (mm)	2024 Price
CG-MR-50V1	51161252	50	45	74.5	\$275.50*

Reaction Vessel, Round Bottom, Flat Flange, 50mL, Mettler-Style

Part No.	Mettler Equivalent	Cap. (mL)	Flange OD (mm)	OAH (mm)	2024 Price
CG-MR-50V2	51161255	50	45	74.5	\$273.85*

Reaction Vessel, Shallow Dish Bottom, EasyMax-Style Flange, 100mL and 400mL, Mettler-Style

Mettler-style, shallow dish bottom, vessels with EasyMax-style flanges.

Part No.	Mettler Equivalent	Cap. (mL)	Flange OD (mm)	OAH (mm)	2024 Price
CG-MR-100V1	51161620	100	57	90	\$327.20
CG-MR-400V1	30072099	400	98	140	\$620.15

Reaction Vessel, Dish Bottom, 250mL, Mettler-Style

Part No.	Mettler Equivalent	Cap. (mL)	Flange OD (mm)	OAH (mm)	2024 Price
CG-MR-250V5	51162635	250	67	127	\$1,762.80*

Reaction Vessel Lid, 7-Neck, 250mL, Mettler-Style

Mettler-style vessel Lid with 24/25 Outer CN, (2) 14/20 Outer SN, and (4) 19/22 Outer SN's, 88mm OD flange.

Part No.	Mettler Equivalent	Cap. (mL)	No. of Necks	Flange OD (mm)	2024 Price
CG-MR-250L5	51162699	250	7	88	\$1,530.45*

Reaction Vessel Lid, 4-Neck, 50mL, Mettler-Style

Mettler-style vessel lid with 14/20 outer CN, (1) GL-14 SN, (1) 10/18 outer SN, (1) 10/10 outer SN, 45mm OD flange.

Part No.	Mettler Equivalent	Cap. (mL)	No. of Necks	Flange OD (mm)	2024 Price
CG-MR-50L1	51161032	50	4	45	\$500.15*

**Reaction Vessel Lid, 6-Neck, 50mL, Mettler-Style**

Mettler-style vessel lids with 14/20 outer CN with 10/18 outer, (2) GL-14 SN's, (1) 10/18 outer SN, (1) 7/12 outer or (1) 7/15 outer SN and 45mm OD flange.

Part No.	Mettler Equivalent	Cap. (mL)	No. of Necks	Flange OD (mm)	OAH (mm)	2024 Price
CG-MR-50L2	51161100	50	6	45	73	\$508.35*
CG-MR-50L3	51160817	50	6	45	63	\$516.55*

**Reaction Vessel Lid, 6-Neck, for 57mm OD Reactor Body, Mettler-Style**

Mettler-style vessel lid with 14/20 outer CN, (3) 14/20 outer SN's, and (2) 10/18 outer SN's. Fits 57mm OD reactor body.

Part No.	Mettler Equivalent	No. of Necks	Fits Reactor Body OD (mm)	2024 Price
CG-MR-100L1	51161619	6	57	\$559.20*

**Reaction Vessel Lid, Short, 6-Neck, 250mL, Mettler-Style**


Mettler-style vessel lid with 14/20 outer CN, (1) 14/20 outer SN, (2) 19/22 outer SN, and (2) GL-14 SN's, 98mm OD flange with notch.

Part No.	Mettler Equivalent	Cap. (mL)	No. of Necks	Flange OD (mm)	2024 Price
CG-MR-250L1	51161188	250	6	98	\$828.65*

**Reaction Vessel Lid, 6-Neck, 250mL, Mettler-Style**

Mettler-style vessel lid with 14/20 outer CN, (3) 14/20 outer SN's, and (2) 19/22 outer SN's, 98mm OD flange without notch.

Part No.	Mettler Equivalent	Cap. (mL)	No. of Necks	Flange OD (mm)	2024 Price
CG-MR-250L2	51161072	250	6	98	\$809.00*



*Not subject to discount.

Micro Coil Condenser, Mettler-Style

140mm OAL, 14/20 top outer joint, 14/20 lower inner joint, GL-14 thread inlet / outlet. Supplied complete with caps with apertures and hose connections.

Mettler-Style Micro Coil Condenser

Part No.	Mettler Equivalent	Joint Size	Thread Size	OAL (mm)	2024 Price
CG-MR-250C1	51161029	14/20	GL-14	140	\$205.00*

Micro Coil Condenser, Jacketed, Mettler-Style

160mm OAL, 14/20 top outer joint, 14/20 lower inner joint, GL-14 thread inlet / outlet. Supplied complete with caps with apertures and hose connections.

Part No.	Mettler Equivalent	Joint Size	Thread Size	OAL (mm)	2024 Price
CG-MR-250C2	—	14/20	GL-14	160	\$382.10*

Micro Allihn Condenser, Mettler-Style


101mm OAL, GL-14 top thread, 10/18 lower inner joint, GL-14 thread inlet/outlet. Supplied complete with caps with apertures and hose connections.

Part No.	Mettler Equivalent	Joint Size	Thread Size	OAL (mm)	2024 Price
CG-MR-50C1	—	10/18	GL-14	101	\$206.15*

Micro Graham Condenser, Mettler-Style

101mm OAL, GL-14 top thread, 10/18 lower inner joint, GL-14 thread inlet/outlet. Supplied complete with caps with apertures and hose connections.

Part No.	Mettler Equivalent	Joint Size	Thread Size	OAL (mm)	2024 Price
CG-MR-50C2	51160813	10/18	GL-14	101	\$288.50*



3800 North Mill Road
Vineland NJ 08360

*Not subject to discount.

†Mettler, Easymax and Optimax are trademarks of Mettler Toledo.

© 2025 Chemglass Life Sciences.
Prices subject to change without notice.
Printed in US.