



3800 N. Mill Road  
Vineland NJ 08360  
Phone: 800.843.1794  
Fax: 800.922.4361  
[www.cglifesciences.com](http://www.cglifesciences.com)

## PRODUCT CHANGE NOTIFICATION

DOC-PCNWEB-002, Rev. A

**PCN Number** 220303

**PCN Title** Groove Profile in the Glass O-Ring Joints, CG-124 Series

**Notification Date** 3/3/2022

**Implementation Date** 3/3/2022

**Part Description** Glass O-Ring Joints

**Affected Part Number(s)** CG-124-Series

**Change Description** Remove the "indent" within the bottom of the O-Ring Groove

**Reason for Change** To ensure the O-Ring is made correctly to O-Ring Groove Standards

**Performance Impact** Better contact with O-Ring

**PCN Correspondence** [QC@cglifesciences.com](mailto:QC@cglifesciences.com)

# PRODUCT CHANGE NOTIFICATION

## FORM - FIT - FUNCTION DEFINITIONS

**FORM:** Refers to such characteristics as external dimensions, weight and size of a part or assembly.

**FIT:** Refers to the ability of the part or feature to connect to, mate with, or join to another feature or part within an assembly.

**FUNCTION:** Is a criteria that is met when the part performs its stated purpose effectively and reliably.