



3800 N. Mill Road
Vineland NJ 08360
Phone: 800.843.1794
Fax: 800.922.4361
www.cglifesciences.com

PRODUCT CHANGE NOTIFICATION

DOC-PCNWEB-002, Rev. A

PCN Number 230328

PCN Title 20 Liter Jacketed Reactor

Notification Date 3/28/2023

Implementation Date 3/28/2023

Part Description 20 Liter Jacketed Reactor with Fixed Hard Lower Drain Seal

Affected Part Number(s) CG-1930-34 to new part number CG-1930-X-34

Change Description Chemglass is changing to a fixed hard seal on the lower outlet instead of an open silicone seal.

Reason for Change Supply Chain issues with getting the high temperature silicone.

Performance Impact None

PCN Correspondence QC@cglifesciences.com

PRODUCT CHANGE NOTIFICATION

FORM - FIT - FUNCTION DEFINITIONS

FORM: Refers to such characteristics as external dimensions, weight and size of a part or assembly.

FIT: Refers to the ability of the part or feature to connect to, mate with, or join to another feature or part within an assembly.

FUNCTION: Is a criteria that is met when the part performs its stated purpose effectively and reliably.