



3800 N. Mill Road
Vineland NJ 08360
Phone: 800.843.1794
Fax: 800.922.4361
www.cglifesciences.com

PRODUCT CHANGE NOTIFICATION

DOC-PCNWEB-002, Rev. A

PCN Number 260211

PCN Title Cleanroom Packaging Update

Notification Date 2/11/2026

Implementation Date 2/13/2026 (Phased in as old bags get used up starting 3/1/26.)

Part Description NP500 Material, 3 Seal Pouch 3/8" Seals, Tear Notches

Affected Part Number(s) All cleanroom assembled Feed Bottles, Lines and Bags.

Change Description Switching to Tear Notched Tabs and Heavier Bag Material for durability.

Reason for Change Tear Notched Top for Customer convenience and also a more durable bag.

Performance Impact None

PCN Correspondence QC@cglifesciences.com

PRODUCT CHANGE NOTIFICATION

FORM - FIT - FUNCTION DEFINITIONS

FORM: Refers to such characteristics as external dimensions, weight and size of a part or assembly.

FIT: Refers to the ability of the part or feature to connect to, mate with, or join to another feature or part within an assembly.

FUNCTION: Is a criteria that is met when the part performs its stated purpose effectively and reliably.