

Fall Savings

# Single Channel, Adjustable Volume, Micropipettes

20%  
OFF\*



## BRANDTECH Transferpette S

- Comfortable, one-handed, all-function operation
- Slim shafts for working with narrow tubes
- Slow force — short pipetting stroke, easy mounting and tip ejection
- Volume lock is clearly visible
- Universal tip cones fit most quality tip brands
- Completely autoclavable at 121 °C

Part No.	Volume Range	Price	Sale Price
CLS-3490-B01	0.1-1µL	<del>\$357.00</del>	\$285.60 ◇
CLS-3490-B02	0.1-2.5µL	<del>\$357.00</del>	\$285.60 ◇
CLS-3490-B03	0.5-10µL	<del>\$357.00</del>	\$285.60 ◇
CLS-3490-B04	2-20µL	<del>\$357.00</del>	\$285.60 ◇
CLS-3490-B05	5-50µL	<del>\$357.00</del>	\$285.60 ◇
CLS-3490-B06	10-100µL	<del>\$357.00</del>	\$285.60 ◇
CLS-3490-B07	20-200µL	<del>\$357.00</del>	\$285.60 ◇
CLS-3490-B08	100-1,000µL	<del>\$357.00</del>	\$285.60 ◇
CLS-3490-B09	500-5,000µL / 0.5-5mL	<del>\$357.00</del>	\$285.60 ◇
CLS-3490-B10	1,000-10,000µL / 1-10mL	<del>\$357.00</del>	\$285.60 ◇

◇ LS Consumable Product Line

### \*RULES - Fall Savings - 20% OFF with Promo Code LS09224

Participation in this offer is restricted to purchases of **CLS-3490-B-series** BRANDTECH Transferpette S, Single Channel, Adjustable Volume Micropipettes. Offer valid from 9/23/2024 through 11/30/2024. Purchaser must specify promo code **LS09224** when placing the order.

**Valid on US and Canadian purchases only.** May not be combined with any other offer or discount. This offer is non-transferable and void where prohibited by law or company policy. Offer is subject to Chemglass Life Sciences' standard terms, conditions and contracts. **Offer not valid through distributors.**

See complete details  
on our website.



www.cglifesciences.com  
tech@cglifesciences.com  
T 1.800.843.1794  
F 1.800.922.4361