

# Magnetic Hot Plate Stirrer

**BUY 2  
GET 1  
FREE**

- **Integral timer** - 1min to 99h:59min
- Ceramic-coated, aluminum top is **scratch-proof & chemical resistant**.
- **Industry leading** magnetic strength with extremely accurate stirring.
- **Bright white LED** and illuminated icons.
- **Fast**, precise heating and regulation.
- **Torque compensation and auto reverse** of stirring direction.
- **PT100 Probe** included.
- **Adjustable safety temperature** circuit and overheat protection
- **Temp Range:** 50-370°C +/- 1°C
- **Speed Range:** 30-1700rpm
- **C Nemko US** (CE, CSA, UL)



Item #	Description	Price
CG-1997-V-50	COMPLETE Digital Hotplate Kit, AREX-6 with probe, with Vertical Support Rod, Boss Head, and Support/Holding Rod, 115V	\$1,144.80

**\*RULES - 3rd Hot Plate Stirrer FREE - Promo Code CM01226**

Participation in this offer is restricted to purchases at list price of CG-1997-V-50, Digital Magnetic Hot Plate Stirrer. Purchase two hot plate stirrers and get a third hot plate stirrer for free. **Offer valid while supplies last.** Purchaser must specify promo code CM01226 when placing the order.

**Valid on US and Canadian purchases only.** May not be combined with any other offer or discount. This offer is non-transferable and void where prohibited by law or company policy. Offer is subject to Chemglass Life Sciences' standard terms, conditions and contracts. **Offer not valid through distributors.**



T 1.800.843.1794  
F 1.800.922.4361  
www.cglifesciences.com  
cs@cglifesciences.com