Convenience Comes to Non-Sterile Syringes!



Bulk-Packed Dispenser Cartons are More Economical and Easier on the **Environment**



prevent spills.

Sterile, Strip-Packed Syringes also available.

See CG-3080-series and CG-3081-series on the web. Convenient dispenser cartons

- Reduced waste no individual wraps to dispose of
- Require less storage space in stockrooms and labs

Polypropylene Syringes, Non-Sterile

These all polypropylene syringes are latex-free and contain no rubber or silicone. The plastic plunger has a smooth draw and a positive safety stop to prevent accidental spills. Bulk-packed in dispenser cartons.

with Luer Locks

Part No.	Vol	ume (mL)	Quantity p	er Carton	Price
CG-3083-01		3	10	00	32.00*
CG-3083-02		5	10	00	33.00*
CG-3083-03		10	10	00	34.00*
CG-3083-04	NFW	20	50	0	27.00*
CG-3083-05	Maraa -	30	30	0	28.00*
CG-3083-06	Sizes!	50	20	0	28.50*

with Luer Slips

Luer tips for sizes 1, 3 and 5mL are centric to the barrel. The 10 through 50mL sizes have luer tips which are offset from the barrel. The 1mL syringe is a long, thin, tuberculin style. Bulk-packed in dispenser cartons.

Part No.	Volu	ume (mL)	Quantity	per Carton	Price
CG-3082-01		1		100	22.00*
CG-3082-02		3		100	27.00*
CG-3082-03		5		100	28.00*
CG-3082-04		10		100	29.00*
CG-3082-05	NEW	20		50	25.80*
CG-3082-06		30		30	28.50*
CG-3082-07	Sizes!	50		20	23.20*



Place your order today!

www.cglifesciences.com customer-service@cglifesciences.com technical-service@cglifesciences.com

> T 1.800.843.1794 F 1.800.922.4361