

96-Well Vial Trays


NEW!
Pre-loaded with 8x30mm Glass Vials


Reusable
Polypropylene
Holder with
Borosilicate
Glass Vials



- high throughput screening
- solubility screening
- polymorph screening
- powder dispensing
- reaction screening
- biological testing

Ideal for Automated Work Stations
Compatible with:

Unchained Labs (Freeslate)
CM2 and CM3 Chemspeed
Waters Acquity UPLC
Shimadzu Prominence
Agilent BenchCel



Sealing mats
available slit
or unslit.



Sealing Mat
CG-1914-05

Chemglass 96-well vial tray / reaction block with the same footprint as a standard 96-well plate. Compatible with 96-well automated work stations. The reusable polypropylene holder accommodates 8x30mm OD borosilicate vials in standard 8x12, alpha-numeric format: A-H, 1-12. Accepts standard SBS (7mm round diameter plug) cap mat; made with 100% premium silicone with PTFE spray coating, available in packs of 5, non- or pre-slit.

Part No.	Description	Qty.	Price
CG-1914-01	Reaction Block Vial Holder, 96-Well with 8x30mm Vials, Complete	1	\$25.20 ▲
CG-1914-02	Reaction Block Vial Holder ONLY, 96-Well	1	\$12.65 ▲
CG-1914-03	Vial Loader, with (96) 8x30mm Vials	1	\$19.65 ▲
CG-1914-04	Sealing Mat, Silicone / PTFE, 6mm Plugs	5	\$48.95 ▲
CG-1914-05	Sealing Mat, Silicone / PTFE, 6mm Plugs, Pre-Slit	5	\$65.30 ▲
CV-2100-0830	Vial Only, 0.75mL, Clear, Glass, Shell, 8x30mm	100	\$11.30 ▲

▲ Subject to Vials / Containers Discount

Freeslate is a registered trademark for Unchained Labs. CM2 and CM3 are registered trademarks for Chemspeed Technologies. Acquity UPLC is a registered trademark for Waters Corp. Prominence is a registered trademark for Shimadzu Scientific Instruments. BenchCel is a registered trademark for Agilent Technologies.

Printed in USA. Prices subject to change. CG-18025IHpdf0120



Vial Loader
CG-1914-03



T 800.843.1794 • F 800.922.4361
info@cglifesciences.com
www.cglifesciences.com