


Sterile, Serological Pipets - Individually Wrapped in Paper / Plastic Wrappers



NEW!
Buy 3 Cases
Get a 4th FREE!
 (Buy at list price. Free case of equal or lesser value.*)
 Use Promo Code
CC2018



barrier plugged
 color-coded bands for quick identification
 easy-to-read, black graduations
 high-clarity polystyrene



NEW CELLcompanion™ Individually Wrapped Serological Pipets are manufactured to ASTM Standard Specification for Serological Pipets E-934. Polystyrene pipets with color bands for quick size identification. Easy-peel, paper / plastic wraps are also color-coded. Pipets are plugged and available in sizes 1, 2, 5, 10, 25 and 50mL. Supplied sterile (gamma irradiated), non-pyrogenic, Rnase- and Dnase-free. Bulk packed in plastic bags.

Individually Wrapped Serological Pipets

Part No.	Capacity (mL)	Color Code	Bag Qty.	Case Qty.	Price
CLS-3481-X001	1	yellow	100	800	\$105.00 ♣
CLS-3481-X002	2	green	100	600	\$99.00 ♣
CLS-3481-X005	5	blue	50	200	\$35.00 ♣
CLS-3481-X010	10	orange	50	200	\$39.50 ♣
CLS-3481-X025	25	red	50	200	\$74.00 ♣
CLS-3481-X050	50	purple	30	90	\$85.00 ♣

Ideal for Cell Culture and other Lab Procedures

- high-clarity pipets
- bi-directional, positive and negative, graduations
- negative graduations allow for additional capacity
- color-coded for quick identification
- premium material barrier plug
- fit standard pipet controllers
- printed indicator arrows on 25 and 50mL sizes as visual aids when using pipet controllers
- anti-drip tips on 25 and 50mL



Contact us to order

T 800.843.1794
 F 800.922.4361

www.cglifesciences.com
 info@cglifesciences.com