

Trade in your old plate today.

# It's worth \$100\*

And you'll get a NEW Digital Stirring Hotplate with **Advanced Features** Only Found in Higher Priced Plates!

Hurry!  
ENDS 8/31/18

## Laser-Verified High Performance Control

- Fast, precise temperature control from 0 to 370°C
- Extremely accurate stirring from 30 to 1700rpm
- Auto-reverse option for enhanced mixing
- Torque compensation maintains set RPM as viscosity changes
- Chemical- and scratch-resistant, CerAlTop™, ceramic-coated top

## Advanced Features

- Integral timer - settable from 1 minute to 99 hours, 59 minutes
- Service life meter - track service interval or service lifetime
- Bright white LED display for high-contrast easy reading
- 135mm diameter, ceramic-coated, aluminum alloy plate

## Engineered for Safety


- Overheat protection & adjustable safety temperature circuit
- Timer automatically turns heat off after elapsed time

Use Promo Code  
**TRADE0818**

See web for required redemption form.



See Details and Case Study on Back



Illuminated Icons for Heat, Probe, Setpoint, Timer and Autoreverse

Easy to read display



Maintains set RPM even with viscosity changes

Scan QR code to watch the video




1-800-843-1794

[www.cglifesciences.com](http://www.cglifesciences.com)  
[tech@cglifesciences.com](mailto:tech@cglifesciences.com)

\* See web for redemption form: [https://chemglass.com/pages/promo\\_offers](https://chemglass.com/pages/promo_offers)  
For each used hotplate stirrer received, take a \$100 discount on a new CG-1997-V-series, Magnetic Hotplate Stirrer. Participate as often as you like. Not to be combined with any other discount. Include promo code **TRADE0818**. Offer ends 8/31/18.

# Fast, Precise Heating and Temperature Control

Set to 100°C, Pie Wedge CG-1991-P-10, 20ml Vials



Part No.	Description	List Price	Trade-In Price
CG-1997-V-50TRA	COMPLETE Digital Hotplate Kit, ARES-6, Temp. Probe, Vertical Support Rod, Boss Head and Support / Holding Rod, 115V	\$799.00*	\$699.00 *
CG-1997-V-10TRA	Digital Hotplate and Temperature Probe, ARES-6, Temperature Probe, 115V	\$745.00*	\$645.00 *

\*Not subject to discount. CG-180681Hpdf1518ou



Offer Ends  
8/31/18



**Digital Stirring Hot Plate Trade In Procedure**

To take advantage of our offer, please send your old plate along with a copy of your proof-of-purchase for a Chemglass, CG-1997-V-series, Digital Stirring Hot Plate (Promo Code **TRADE0818**) and this completed form to:

**Chemglass Trade In / Trade Up**  
Chemglass Life Sciences  
3800 N Mill Road  
Vineland, NJ 08360

Name \_\_\_\_\_

Job Title \_\_\_\_\_

Organization Name \_\_\_\_\_

Division \_\_\_\_\_ Dept. \_\_\_\_\_

Bldg. / Lab / Room \_\_\_\_\_ P.O. Box \_\_\_\_\_

Street Address \_\_\_\_\_

City \_\_\_\_\_ State / Province \_\_\_\_\_ Zip / Postal Code \_\_\_\_\_

Phone \_\_\_\_\_ E-mail \_\_\_\_\_

Model / Brand Returned	Serial No.
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____

**Equipment Decontamination Verification**

Returned unit(s) must be accompanied by a completed and signed **EQUIPMENT DECONTAMINATION VERIFICATION** to ensure each item is free of RADIOACTIVE, TOXIC or BIOHAZARD contamination.

**Use separate forms for different types of contamination. Please choose only one of the following:**

- This equipment has not been used in a laboratory / location exposing it to biohazards or radioactive materials.
- This equipment was decontaminated\* from biohazardous materials with: \_\_\_\_\_  
\*Example: autoclave, 10% bleach, ethylene oxide
- This equipment was decontaminated and tested\* for radioactivity with: \_\_\_\_\_  
\*Example: wipe test with results (3H, 14C, (-emitters), Geiger counts, etc.

**SIGNATURE** \_\_\_\_\_ **DATE** \_\_\_\_\_

**NAME** (print) \_\_\_\_\_ **PHONE** \_\_\_\_\_

Federal regulations require that all instruments sent for service shall be free of radioactivity and biohazards. Chemglass reserves the right to refuse to accept this equipment if these procedures are not performed.