

PRODUCT CHANGE NOTIFICATION

PCN Number 112119

PCN Title Fernbach Flask Neck Wall Thickness Change

Notification Date 11/21/2019

Implementation Date 11/25/2019

Part Description Fernbach Flasks

Affected Part Number(s) CG-1565-02, 03, 04 ,05 ,06, 07 ,08, CLS-2020-11 & 16,
CG-2026-22, CLS-2036-04, 05, 06, 07, 08, 09,10 & 11,
CLS-2044-03, 05, 06, & 08, CLS-2046-03, 04, 05, 06, 08, 09, & 10 ,
CLS-2046-10, CLS-2048-04, 05, 06, 08, 09, & 10, CLS-2050-02, 04, 06

Change Description Wall Thickness of neck for Fernbach Flasks changed from
1 1/2" medium wall (38mm x 2.4mm) tubing to 1 1/2" heavy
wall (38mm x 4.0mm) tubing.

Reason for Change Increase wall thickness for better durability.

Performance Impact NONE

Distinguishable Characteristics Increased wall thickness with tubing.

PCN Correspondence QC@cglifesciences.com