



3800 N. Mill Road
Vineland NJ 08360
Phone: 800.843.1794
Fax: 800.922.4361
www.cglifesciences.com

PRODUCT CHANGE NOTIFICATION

DOC-PCNWEB-002, Rev. A

PCN Number	<u>200722</u>
PCN Title	<u>CLS-1481-21</u>
Notification Date	<u>7/22/2020</u>
Implementation Date	<u>7/23/2020</u>
Part Description	<u>Cap, GL-45, Blue, Open Top, (1-5/16" Dia. Hole) Polypropylene</u>
Affected Part Number(s)	<u>ALL CLS-1437 Series Port Assemblies, CLS-1427 Series, CLS-1437-F Series, and CLS-1437-H Series</u>
Change Description	<u>Chemglass to manufacture cap. Color to be slightly different and Chemglass Bio Sciences Logo on cap.</u>
Reason for Change	<u>Improved lead time and part tolerance due to one step molding process and eliminates the machining process the previous part needed.</u>
Performance Impact	<u>N/A</u>
PCN Correspondence	<u>QC@cglifesciences.com</u>

PRODUCT CHANGE NOTIFICATION

FORM - FIT - FUNCTION DEFINITIONS

FORM: Refers to such characteristics as external dimensions, weight and size of a part or assembly.

FIT: Refers to the ability of the part or feature to connect to, mate with, or join to another feature or part within an assembly.

FUNCTION: Is a criteria that is met when the part performs its stated purpose effectively and reliably.