



Disposable Filter Funnels & Adapter System

**Q U I C K
E A S Y
C L E A N**

Chemglass Adapter System for use with Fritted, Disposable Filter Funnels

Typically used for solid-liquid separations in chemistry and biological laboratories, set-ups consist of a funnel, vacuum adapter and optional threaded connector for vials.

The vacuum adapters (CG-1049-J-series) are designed so that vacuum can be controlled with a PTFE stopcock to prevent filtrates from being pulled into the vacuum line. Stems on the filter funnels are 80mm long and prevent contamination of the vacuum adapter during the filtration process.

Part No.	Funnel Adapters with Vacuum Stopcocks	Price
CG-1049-J-14	Outer to 14/20 Inner	86.00 *
CG-1049-J-24	Outer to 24/25 Inner	88.60 *
CG-1049-J-29	Outer to 29/26 Inner	88.60 *
CG-1049-J-30	Vial, Outer to 20-400 GPI Thread	86.15 *

Part No.	Accessory Threaded Connectors	Price
CG-1318-23	20-400 to 20-400 GPI Thread	18.30 *
CG-1318-24	20-400 to 22-400 GPI Thread	18.70 *
CG-1318-25	20-400 to 24-400 GPI Thread	24.40 *

Disposable Filter Funnels with Polyethylene Frits

Available in Bulk Bags or Small-Quantity Trays


Bulk Packed in ZipLock Bags

Part No.	Funnel Cap. (mL)	Frit Porosity (µm)	Top Style	Case Qty.	Price
OP-6602-10	18	10	Flared	100	246.40 *
OP-6602-12	40-60	10	Flared	100	271.40 *
OP-6602-13	40-60	40	Flared	100	272.90 *
OP-6602-14	110	10	Straight	50	207.60 *
OP-6602-16	300	20	Straight	30	240.85 *
OP-6602-18	500	20	Straight	20	224.00 *

Shelf Pack Convenience Trays

Packaged in a smaller, convenient, shrink-wrapped, SBS-coated, white tray complete with SBS partitions inserted between the top and bottom layers.

Part No.	Funnel Cap. (mL)	Frit Porosity (µm)	Top Style	Case Qty.	Price
OP-6602-C10	18	10	Flared	24	65.25 *
OP-6602-C12	40	10	Flared	24	72.50 *
OP-6602-C14	110	10	Straight	12	57.05 *



US Patent Pending

*Not Subject to Discount

Disposable Filter Funnels & Adapter System

Disposable Filter Funnels with Polyethylene and Diatomaceous Earth Frits

Available in Bulk Bags or Small-Quantity Trays

Bulk Packed in ZipLock Bags

Part No.	Funnel Cap. (mL)	DE Content	Top Style	Case Qty.	Price
OP-6603-10	18	500mg	Flared	50	215.80 *
OP-6603-12	40-60	1,500mg	Flared	50	231.00 *
OP-6603-13	150	5g	Straight	30	200.80 *

Shelf Pack Convenience Trays

Packaged in a smaller, convenient, shrink-wrapped, SBS-coated, white tray complete with SBS partitions inserted between the top and bottom layers.

Part No.	Funnel Cap. (mL)	DE Content	Top Style	Case Qty.	Price
OP-6603-C10	18	500mg	Flared	24	112.70 *
OP-6603-C12	40	500mg	Flared	24	120.70 *
OP-6603-C13	150	500mg	Straight	12	88.50 *

Cone-Shaped, Disposable Filter Funnels, Molded Polypropylene with Polyethylene Frits

Disposable Filter Funnels, Cone-Shaped, PP with PE Frits

Part No.	Funnel Cap. (mL)	Frit Porosity (µm)	Case Qty.	Price
OP-6605-08	8	10	100	195.00 *
OP-6605-10	20	10	100	244.95 *
OP-6605-12	60	10	100	241.10 *

Chemglass Adapter System for use with Fritted, Disposable Filter Funnels

Purchase components separately.

Part No.	Funnel Adapters with Vacuum Stopcocks	Price
CG-1049-J-14	Outer to 14/20 Inner	86.00 *
CG-1049-J-24	Outer to 24/25 Inner	88.60 *
CG-1049-J-29	Outer to 29/26 Inner	88.60 *
CG-1049-J-30	Vial, Outer to 20-400 GPI Thread	86.15 *

Part No.	Accessory Threaded Connectors	Price
CG-1318-23	20-400 to 20-400 GPI Thread	18.30 *
CG-1318-24	20-400 to 22-400 GPI Thread	18.70 *
CG-1318-25	20-400 to 24-400 GPI Thread	24.40 *

