



pH METERS



OHAUS Bench pH Meter with Two Point Calibration - OHAUS ST2100-F

- minimal training
- quick setup
- easy operation



- easy operation with large LCD and user-friendly software
- built-in electrode holder offers convenient, all-in-one operation
- accurate results ensured with auto buffer recognition, electrode condition icon, automatic temperature compensation and reference

The Starter 2100 is an economical choice for simple pH measurement. A large LCD display, clearly marked keys and user-friendly software make the ST2100 perfect for even novices. Accurate results are driven by auto buffer recognition, electrode condition icon, automatic temperature compensation and reference input. One year warranty.

Complete unit includes: 2100 meter, ST210 2-in-1 plastic refillable pH electrode, STTEMP30 temperature electrode and pH buffer powder set: 4.01, 7.00, 10.01.

CG-1876-2100 pH Meter, Digital, Bench Meter, 2 Point Calibration **\$476.00 ***

OHAUS Bench pH Meter with Three Point Calibration - OHAUS ST3100-F

- reliability
- accuracy
- connectivity



- intuitive software includes three point calibration, three buffer groups, automatic and manual endpoint functions as well as a 99 item library
- adjustable, standalone electrode holder
- documenting measurement data is made simple through the RS232 which can be used to connect peripheral devices such as printers or computers to the meter

The Starter 3100 pH meter includes a variety of user-friendly features providing advanced operational efficiency. Intuitive software and a well-organized, large LCD display make it easy to use. Durable and intelligently designed. One-year warranty.

Complete unit includes: 3100 meter, ST310 3-in-1 plastic refillable pH electrode, pH buffer powder set: 4.01, 7.00, 10.01 and in-use cover.

CG-1877-3100 pH Meter, Digital, Bench Meter, 3 Point Calibration **\$656.00 ***

OHAUS Portable pH Meter - OHAUS ST300-B

- convenience
- durability
- simplicity



- automatic and manual temperature compensation for accurate readings
- ORP or temperature electrode connectivity for additional functionality
- IP54 protection shields water / dust damage
- 30 measurement library stores data for future reference
- auto buffer recognition stores calibration data to help avoid calibration process errors

The ST300 portable pH meter offers convenience, reliability and durable construction in an ergonomic compact design. Intuitive software and a well-organized display make it user friendly. Convenient, built-in, retractable stand. Five, clearly marked keys offer simplicity for novices and professionals alike. IP54 protection shields the portable meter from damage by dust and water. Quickly recall calibration data and information with a single touch. Connection of an ORP or temperature electrode adds additional functionality. One-year warranty.

Complete unit includes: 300 Meter, IP54 sets, electrode clip, wrist strap and 4 AAA batteries (electrode not included).

CG-1875-300 pH Meter, Digital, Portable **\$395.00 ***

OHAUS is a reg. trademark.
*Not subject to discount.

See Probes and Buffers on back ►



pH METERS



OHAUS pH Meter Electrodes



CG-1875-T310



CG-1875-T320



CG-1875-T350

OHAUS 3-in-1 pH Electrodes with Integrated Temperature Sensors

- sturdy and durable electrodes built to withstand daily use
- temperature sensors power automatic temperature compensation and ensure accurate measurements
- electrodes fit perfectly in the electrode holders on OHAUS bench meters and electrode clips on all OHAUS portable meters.

CG-1875-T310	pH Probe, 3-in-1 Refillable pH Electrode, Integrated Temperature Probe, Plastic Shaft	\$123.00 *
CG-1875-T320	pH Probe, 3-in-1 Non-Refillable pH Electrode, Gel Electrolyte, Integrated Temperature Probe, Plastic Shaft	\$123.00 *
CG-1875-T350	pH Probe, 3-in-1 Refillable pH Electrode, Integrated Temperature Probe, Glass Shaft	\$153.00 *



OHAUS pH Buffer Solution, 250mL Bottles

CG-1875-B004	4.01 pH Buffer Solution, 250mL Bottle	\$14.00 †
CG-1875-B007	7.00 pH Buffer Solution, 250mL Bottle	\$14.00 †
CG-1875-B010	10.01 pH Buffer Solution, 250mL Bottle	\$14.00 †

*Not subject to discount. †Reusable Glass / Plastic.
Prices subject to change without notice.

CG-18031BC1.5M0218