

HandyStep® touch

Limited Time
Introductory Price!

The next generation repetitive pipette



The perfect team

Choose the HandyStep® touch for your repetitive pipetting tasks—dispensing, auto-dispensing and the ability to store favorites. The HandyStep® touch *S* adds such advanced features like multi-aspiration, sequential dispensing and titration. Paired with the new PD-Tip™ // Dispenser Tips with tip size coding, they make the perfect team for easy pipetting.

Mode	Function	HandyStep® touch	HandyStep® touch <i>S</i>
Multi Dispensing (MULTI-DISP)	Successive dispensing of a set volume	✓	✓
Auto Dispensing (AUTO-DISP)	Automatic dispensing by the instrument	✓	✓
Pipetting (PIP)	Pipetting a defined volume	✓	✓
Sequential Dispensing (SEQ-DISP)	Successive dispensing of different volumes	—	✓
Multi Aspirating (MULTI-ASP)	Aspirating different volumes	—	✓
Titration	Manual or successive delivery of liquids for titration	—	✓
Favorites	Saving favorites for quick access	✓	✓
Settings	Individual customization of the device (e.g. brightness or device name)	✓	✓



Functions of the HandyStep® touch



Expanded functionality on the HandyStep® touch *S*



(800)-843-1794 | chemglass.com



www.brandtech.com

To Place Orders: Call -800-843-1794 • Email: cs@cglifesciences.com • Visit: www.chemglass.com

Technical data

HandyStep® touch

HandyStep® touch S

Repetitive pipette, includes universal AC adapter (for USA, Europe, UK, Japan, Australia and New Zealand), shelf/rack mount, sample PD-Tip™ II, operating manual and performance certificate.

Inductive charging with the optional charging stand



Support stand with HandyStep® touch

Description	BrandTech® Cat. No.	Chemglass Cat. No.	2020 List Price	Intro Price
HandyStep® touch	705200IP	CLS-3468-200	\$1,055.00	\$949.00
HandyStep® touch S	705210IP	CLS-3468-210	1,299.00	1,099.00

Inductive charging stand bundles

HandyStep® touch & Charging stand	705201IP	CLS-3468-200P	\$1,148.00	1,098.00
HandyStep® touch S & Charging stand	705211IP	CLS-3468-210P	1,525.00	1,320.00

Accessories

Description	BrandTech® Cat. No.	Chemglass Cat. No.	2020 List Price	Intro Price
Charging stand	705220IP	CLS-3468-300	\$345.00	\$275.00
Support stand	705230	CLS-3468-302	144.00	—
Shelf/rack mount	705235	CLS-3468-306	35.00	—
Universal AC adapter	705250	CLS-3468-304	85.00	—
Spare battery	705225	CLS-3468-308	80.00	—

PD-Tip™ II

Precision dispenser tips coded for automatic size recognition.

Capacity mL	Pack of	BrandTech® non-sterile	Chemglass Non-Sterile	2020 List Price	Intro Price	BrandTech® BIO-CERT®	Chemglass BIO-CERT®	2020 List Price	Intro Price
0.1	100	705700IP	CLS-3468-N-000	\$124.00	\$97.40	705730IP	CLS-3468-A-030	\$170.00	\$137.60
0.5	100	705702IP	CLS-3468-N-002	124.00	97.40	705732IP	CLS-3468-A-032	170.00	137.60
1.0	100	705704IP	CLS-3468-N-004	124.00	97.40	705734IP	CLS-3468-A-034	170.00	137.60
1.25	100	705706IP	CLS-3468-N-006	124.00	97.40	705736IP	CLS-3468-A-036	170.00	137.60
2.5	100	705708IP	CLS-3468-N-008	124.00	97.40	705738IP	CLS-3468-A-038	170.00	137.60
5	100	705710IP	CLS-3468-N-010	124.00	97.40	705740IP	CLS-3468-A-040	170.00	137.60
10	100	705712IP	CLS-3468-N-012	124.00	97.40	705742IP	CLS-3468-A-042	170.00	137.60
12.5	100	705714IP	CLS-3468-N-014	124.00	97.40	705744IP	CLS-3468-A-044	170.00	137.60
25	50 + 1 adapter	705716IP	CLS-3468-N-016	84.00	75.20	705746IP*	CLS-3468-A-046	52.00	48.40
50	25 + 1 adapter	705718IP	CLS-3468-N-018	42.00	39.80	705748IP	CLS-3468-A-048	52.00	49.40

Set of PD-Tip™ II, non-sterile, 20 pieces per size: 0.5/1/1.25/2.5/5/10/12.5 mL

705720IP CLS-3468-N-050 ~~168.00~~ 122.00

Adapter for 25 and 50 mL PD-Tip™ & PD-Tip™ II, PP, autoclavable

Pack of 10

702398 CLS-3468-K-020 96.80 —

Pack of 5

702399 CLS-3468-K-010 76.60 —

*BIO-CERT® : Pack of 25 + 1 adapter



For more information visit www.brandtech.com or scan the QR Code.

HandyStep® touch, HandyStep® touch S and PD-Tips II have patents & patent applied for features.

California Residents: For more information concerning California Proposition 65, please refer to www.brand.de/calprop65.

Subject to technical modification without notice. Prices subject to change. Introductory Pricing ends 12/31/2020.



(800)-843-1794 | chemglass.com

888-522-2726 | www.brandtech.com

To Place Orders: Call -800-843-1794 • Email: cs@cglifesciences.com • Visit: www.chemglass.com