



3800 N. Mill Road  
Vineland NJ 08360  
Phone: 800.843.1794  
Fax: 800.922.4361  
[www.cglifesciences.com](http://www.cglifesciences.com)

## PRODUCT CHANGE NOTIFICATION

DOC-PCNWEB-002, Rev. A

**PCN Number** 210708

**PCN Title** CG-1964-55

**Notification Date** 7/8/2021

**Implementation Date** 8/19/2021

**Part Description** Adapter, 1 1/2" Beaded Pipe to 3/4" Compression Fitting, PTFE

**Affected Part Number(s)** CG-1964-55, Adapter, 1 1/2" Beaded Pipe to 3/4" Compression Fitting, PTFE

**Change Description** Change material from 25% glass filled PTFE to Virgin PTFE

**Reason for Change** Glass filled material is cracking and for additional material durability

**Performance Impact** None

**PCN Correspondence** [QC@cglifesciences.com](mailto:QC@cglifesciences.com)

# PRODUCT CHANGE NOTIFICATION

## FORM - FIT - FUNCTION DEFINITIONS

**FORM:** Refers to such characteristics as external dimensions, weight and size of a part or assembly.

**FIT:** Refers to the ability of the part or feature to connect to, mate with, or join to another feature or part within an assembly.

**FUNCTION:** Is a criteria that is met when the part performs its stated purpose effectively and reliably.