



661 Route 23 South, Wayne, New Jersey 07470
Tel: 973-694-0500 * Fax: 973-694-7199
Internet Address: <http://www.bel-art.com>

CERTIFICATE OF ANALYSIS

This is to verify that all items noted below are in conformance with the Purchase Order Contract, drawings, specifications, and other applicable documentation; that all process certifications, material certifications, are on file at this or another Bel-Art facility and are available for review by the customer and or Government Representatives upon request.

Product #: 10794-0025 **Mfg. Date:** NA **Lot #:** General

Description: 2.5 Gallon Polypropylene Carboy

Material: Polypropylene, Thermo Plastic Elastomer Washer in cap

Dimensions: OD 9.75", HT w/cap 15.25", HT w/o cap 15", HT to base of neck 10.5",
Nick ID 2.75", 83mm Buttress Cap. Capacity 10 liters / 2.5 Gallons

Other: Ribbed cap for easy removal, Autoclavable at 121°C @ 15 p.s.i.g. for 20 minutes
(WITH CAP REMOVED). Handles for easy carrying and pouring, Based on the
process and materials involved in the manufacturing of this product, this product is
considered to be BSE/TSE protein free. Based on the information provided to us this
product does not contain or come in contact with lead, latex, Mercury, or other
To the best of our knowledge based on information provided to us this product does
or come in contact with any animal origin. Materials meet FDA regulation 21 CFR
177.1520 and USP Class VI, Manufactured: India

Reviewed and Approved By:

Signature: 

Printed Name: Carl Sisco

Title: Technical Support Manager

Date: 12/6/2018

SUBSIDIARIES: MADDAK, INC * APPLIED COATING, INC
ADDITIONAL PLANTS IN POCOMOKE, MD * KANSAS CITY, MISSOURI