

Industrial Evaporation



MADE IN
GERMANY

The average operational **lifespan of 10 years** is backed by a **3 year warranty** and makes your purchase a worthwhile investment.

Illuminated bath and evaporation flask allows for superior **visual control**

A guard hood consisting of a metal frame and safety glass covers the bath completely while providing **clear visible monitoring**

The brand-new **EASY LOCK** flask attachment allows for single-user safe flask exchange

Excellent user safety: condenser and receiving flasks are protected by high-impact transparent PMMA doors

Receiver cassettes **guard against glassware breakage** while providing a self-standing support

Increased throughput at reduced energy consumption due to an industry leading vacuum leakproof

Features **level regulation system** to prevent from spilling over and potential floods in lab

Unique self-aligning flask support system allows for immediate **"one-person operation"**

Non-visible heating coils guarantee the most convenient and **fastest cleaning of the bath**

Intuitive touch screen control for ease of operation

Unit **powers off completely** in case of a bath overtemperature situation, overpressure and glassware breakage

No hassle or safety concerns due to low bath fluid levels: the self-filling water bath takes the worry out of running dry

Save time with up to 250 pre-programmed profiles of **up to 999 process gradients**

Automatic module Distimatic Industrial - For batch and continuous evaporation without limits - an affordable and more flexible alternative to common 50 or 100-liter evaporators