

CG-1968-D-530 COMPONENTS

Qty	Part Number	Description
1	CG-1014-C-03	90° Inlet Adapter, 45/50 Inner, 1/2" Swagelock
1	CG-141-T-10	80mm Quick Release Clamp, PTFE Coated
1	CG-145-09	Keck Clamp, Brown, fits 45/50 Joint size, 10/Pkg
1	CG-147-81	80mm Perfluoro O-Ring
2	CG-150-09	#65 Stainless Steel Spherical Clamp
1	CG-1947-A100	100mm Clamp, Reactor Support
2	CG-1965-D-210	Adapter, 0-14mm Valve to 45/50 Inner joint, 1/2" Swagelock on Sidearm
2	CG-1965-D-215	Adapter, 0-14mm Valve to 45/50 Inner joint, 1/2" Swagelock on Sidearm, Air Inlet Valve on Side
2	CG-1965-D-230	20 Liter Receiving Flask, 2 Neck, 45/50 Outer CN, Angled 45/50 Outer SN, 0-20 (CG-511-02) Lower Drain with 3/4" Swagelock on Sidearm
1	CG-1968-D-102	Mobile Distillation Cart, Dual, Heavy Duty, PTFE Coated Aluminum and Stainless Steel
1	CG-1968-D-125	Side Swivel Extension Assembly, for CG-1968-X Support Stand
1	CG-1968-D-300	160° Adapter, 65/40 Polished Socket joint to 80mm Flat Flange
1	CG-1968-D-302	Distilling Column, Vacuum Jacketed, 45/50 Outer on Top, 65/40 O-Ring Ball joint on Bottom
1	CG-1968-D-304	Bent Connecting Adapter, 45/50 Inner Sealed to 65/40 O-Ring Ball
1	CG-1968-D-306	Double Cooling Coil Condenser, 65/40 Socket on Top, Angled 45/50 Inner on Bottom, CG-300-03 Hose Connection on Inlet/Outlet & Vacuum Take-Off
1	CG-1968-D-308	Distillation Splitter, 45/50 Outer on Top, 6mm Teflon Stopcock Sealed to 1/2" Swagelock on Side of Splitter, (2) 0-14 Valve with 3/4" Swagelock Side Arms
1	CG-1968-E-420	400mm Reaction Kettle Lid, 45/50 Center Neck, 4-45/50 Side Necks 2-80mm Duran Flange
1	CG-1975-04	PTFE Tubing, 1/2" O.D. x 3/8" I.D., 20ft Length
6	CG-1975-21	1/2" to 1/2" PFA Compression Fitting
2	CG-1975-22	1/2" to 3/4" PFA Compression Fitting
1	CG-1975-30	Tubing Cutter
2	CG-2040-C-01	Heavy Duty Clamp Holder, Fit 1/4" to 5/8" Support Rod
2	CG-9205-E-11	Clamp, Three Prong Extension, Large, Dual Adjust, 1/2" O.D. x 12" Arm
2	CG-9220-E-01	Small Chain Clamp, 12" Extended Arm
4	CG-9253-10	Small Kwik Mount Klamp