

CLS-3788-F51



STANDARD FREESTANDING UNDERCOUNTER FREEZER

The Standard Free-Standing Undercounter Freezers are intended for general purpose applications only where precise temperature control and maintenance are not essential and there are no requirements for temperature alarms.

Freezers feature a manual defrost operating system to limit product-damaging temperature spikes found in auto defrost units. Periodic manual defrosting is required.



Model	Cu. Ft.	Defrost	Doors	Shelves	W"	D"	H"	RL
CLS-3788-F51	4	Manual	1 Solid	4 Drawers	20	22 5/8	33 3/4	R600a

- Mechanical thermostat
- Adjustable operating temperature range:
Freezer: -15°C – -25°C
- Keyed door lock
- Magnetic door gasket for positive seal
- Recessed door handle
- Leveling legs
- Probe access port (3/8")

- White cabinet and door
- High density urethane foam cabinet and door insulation
- HFC-free refrigerant (R600a), foam insulation and packaging
- One year parts and labor warranty
- Note: This unit must have 4" clearance on sides and back for adequate ventilation
- Defrost:
Freezer: Manual defrost
- UL/C-UL listed