

## CG-305 O-RING, VITON

Viton o-ring is stocked in brown 75 durometer. They are suitable for acids, aliphatic hydrocarbons, oxidizing agents and all petroleum products. Viton has a low compression set and high temperature resistance. Viton is especially good for vacuum service and low gas permeability and has a useful temperature range of -30 to 200° C. For complete chemical resistance and fluid compatibility, see chart in the front of this catalog. **FOR VITON O-RINGS USED WITH SCHOTT FLANGES, SEE CG-147.**



CG-305

Part Number	O-Ring Size	ID x Cross Section (mm)	Pkg Qty
CG-305-5-101	5-101	2.54 x 0.097	12
CG-305-006	006	2.90 x 1.78	12
CG-305-007	007	3.68 x 1.78	12
CG-305-008	008	4.47 x 1.78	12
CG-305-009	009	5.28 x 1.78	12
CG-305-010	010	6.07 x 1.78	12
CG-305-011	011	7.65 x 1.78	12
CG-305-012	012	9.25 x 1.78	12
CG-305-013	013	10.82 x 1.78	12
CG-305-014	014	12.42 x 1.78	12
CG-305-017	017	17.17 x 1.78	12
CG-305-018	018	18.77 x 1.78	12
CG-305-019	019	20.35 x 1.78	12
CG-305-020	020	21.95 x 1.78	12
CG-305-022	022	25.12 x 1.78	12
CG-305-024	024	28.30 x 1.78	12
CG-305-027	027	33.05 x 1.78	12
CG-305-104	104	2.84 x 2.62	12
CG-305-105	105	3.63 x 2.62	12
CG-305-107	107	5.23 x 2.62	12
CG-305-108	108	6.02 x 2.62	12
CG-305-109	109	7.59 x 2.62	12
CG-305-110	110	9.19 x 2.62	12
CG-305-111	111	10.77 x 2.62	12
CG-305-112	112	12.37 x 2.62	12
CG-305-113	113	13.94 x 2.62	12
CG-305-114	114	15.54 x 2.62	12
CG-305-115	115	17.12 x 2.62	12
CG-305-116	116	18.72 x 2.62	12
CG-305-117	117	20.29 x 2.62	12
CG-305-118	118	21.89 x 2.62	12
CG-305-120	120	25.07 x 2.62	12
CG-305-121	121	26.64 x 2.62	6
CG-305-122	122	28.24 x 2.62	6
CG-305-126	126	34.59 x 2.62	6
CG-305-128	128	37.77 x 2.62	6

LISTING CONTINUED ON NEXT PAGE...

LISTING CONTINUED FROM PREVIOUS PAGE...

## CG-305 O-RING, VITON

Part Number	O-Ring Size	ID x Cross Section (mm)	Pkg Qty
CG-305-129	129	39.34 x 2.62	6
CG-305-132	132	44.12 x 2.62	6
CG-305-136	136	50.47 x 2.62	6
CG-305-137	137	52.07 x 2.62	6
CG-305-146	146	66.34 x 2.62	6
CG-305-201	201***	4.00 x 3.00	6
CG-305-204	204	9.12 x 3.53	6
CG-305-205	205	10.69 x 3.53	6
CG-305-206	206	12.29 x 3.53	6
CG-305-208	208	15.47 x 3.53	6
CG-305-210	210	18.66 x 3.53	6
CG-305-212	212	21.82 x 3.53	6
CG-305-214	214	25.00 x 3.53	6
CG-305-215	215	26.58 x 3.53	6
CG-305-217	217	29.75 x 3.53	6
CG-305-218	218	31.35 x 3.53	6
CG-305-219	219	32.92 x 3.53	6
CG-305-220	220	34.52 x 3.53	6
CG-305-222	222	37.70 x 3.53	6
CG-305-223	223	40.87 x 3.53	3
CG-305-225	225	47.22 x 3.53	3
CG-305-226	226	50.39 x 3.53	3
CG-305-229	229	59.92 x 3.53	3
CG-305-232	232	69.44 x 3.53	3
CG-305-240	240	94.84 x 3.53	3
CG-305-242	242	101.19 x 3.53	3
CG-305-243	243	104.37 x 3.53	3
CG-305-258	258	151.99 x 3.53	3
CG-305-267	267	209.14 x 3.53	1
CG-305-317	317	23.16 x 5.33	1
CG-305-331	331	56.52 x 5.33	1
CG-305-338	338	78.74 x 5.33	1
CG-305-339	339	81.92 x 5.33	1
CG-305-341	341	88.27 x 5.33	1
CG-305-348	348	110.49 x 5.33	1
CG-305-349	349	113.67 x 5.33	1
CG-305-359	359	145.62 x 5.33	1

\*\*\* NOTE: NOT A STANDARD SIZE #201 O-RING



**CG-305**